Liquick Cor - CREATININE

Kit name	(EN) Cat. No
Liquick Cor-CREATININE mini	2-220
Liquick Cor-CREATININE 30	2-232
Liquick Cor-CREATININE 60	2-233
Liquick Cor-CREATININE 500	2-299

INTENDED USE

Diagnostic kit for determination of creatinine concentration used both for manual assay (Sample Start and Reagent Start method) and in several automatic analysers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

METHOD PRINCIPLE

Modified Jaffe's method, without deproteinization. In alkaline solution picrate reacts with creatinine to form a yellow-red 2,4,6-trinitrocylohexadienate. The colour intensity is proportional to the creatinine concentration.

REAGENTS

Package	Liquick Cor-CREATININE mini	Liquick Cor-CREATININE 30
1-CREATININE	2 x 24 ml	5 x 24 ml
2-CREATININE	1 x 12 ml	1 x 30 ml
3-STANDARD	1 x 1 ml	1 x 2 ml
	Liquick Cor-CREATININE 60	Liquick Cor-CREATININE 500
1-CREATININE	5 x 48 ml	3 x 400 ml
2-CREATININE	1 x 60 ml	1 x 300 ml

3-STANDARD is creatinine standard solution: 177 µmol/l (2.0 mg/dl).

The reagents when stored at 15-25°C are stable up to expiry date printed on the package. The reagents are stable for 2 weeks on board the analyser at 2-10°C.

Working reagent preparation and stability

Assay can be performed with use of separate 1-CREATININE and 2-CREATININE reagents or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-CREATININE with 1 part of 2-CREATININE. Avoid foaming.

Stability of working reagent: 4 weeks at 2-8°C
7 days at 15-25°C

Working reagent should be stored in tightly closed vials. If stored in open vial, it is stable for 1 day at 15-25°C!

Concentrations in the test

sodium hydroxide ≤ 450 mmol/l
carbonate buffer ≤ 150 mmol/l
picric acid ≤ 38.8 mmol/l

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents are usable when the absorbance of the working reagent is less than 0.225 (read against distilled water, wavelength λ = 500 nm, cuvette l = 1 cm, at temp. 25°C).
- 1-CREATININE and 2-CREATININE meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-CREATININE contains sodium hydroxide.

2-CREATININE contains picric acid.

Danger

H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330+P331 IF SWALLOWED: rinse mouth.
Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 500 nm (492 nm);
- thermostat at 25°C;
- general laboratory equipment;

SPECIMEN

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine, without preservatives⁹.

Urine preparation: before analysis urine sample should be diluted 100-fold with 0.9% NaCl, and the results multiplied by 100. Mix well probes before measurement. Specimen can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

Applications for analyzers are available on request.

Manual procedure

wavelength 500 nm (492 nm)
temperature 25°C
cuvette 1 cm
reaction type Fixed time

Sample Start method

Pipette into the cuvettes:

	test (T)	standard (S)
working reagent	1000 µl	1000 µl
standard	-	100 µl
sample	100 µl	-

Bring up to the temperature of determination. Then add:

standard	-	100 µl
sample	100 µl	-

Mix well and after exactly 30 sec read absorbance A1 of the test (T) and standard (S) against air. After next 60 sec repeat absorbance reading (A2) and calculate ΔA (A2 - A1) for the test and standard.

Reagent Start method

The determination can be also performed with use of separate 1-CREATININE and 2-CREATININE reagents.

Pipette into the cuvettes:

	test (T)	standard (S)
1-CREATININE	1000 µl	1000 µl
standard	-	100 µl
sample	100 µl	-

Mix well, then add:

2-CREATININE	250 µl	250 µl
--------------	--------	--------

Mix well and after exactly 30 sec read absorbance A1 of the test (T) and standard (S) against air. After next 60 sec repeat absorbance reading (A2) and calculate ΔA (A2 - A1) for the test and standard.

Calculation

$$\text{creatinine concentration} = \frac{\Delta A(T)}{\Delta A(S)} \times \text{standard concentration}$$

REFERENCE VALUES⁷

serum / plasma	mg/dl	µmol/l
female	0.6 - 1.1	53 - 97
male	0.7 - 1.3	62 - 115
24-hours urine	mg/kg/24h	µmol/kg/24h
female	11 - 20	97 - 177
male	14 - 26	124 - 230

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine with each batch of samples.

For calibration when using the manual methods CREATININE STANDARD 2 (Cat. No 5-123) is recommended.

For calibration of the automatic analysers systems CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174 and 5-176) and LEVEL 2 (Cat. No 5-175 and 5-177) are recommended.

The calibration curve should be prepared every week, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 0.2 mg/dl (17.68 µmol/l).
- Linearity:** up to 25 mg/dl (2210 µmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, triglycerides up to 500 mg/dl, ascorbate up to 62 mg/l and bilirubin up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.40	0.05	3.33
level 2	4.76	0.05	0.96
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.35	0.07	5.07
level 2	4.48	0.22	4.82

Method comparison

A comparison between creatinine values determined at **Biolis 24i Premium** (y) and at **ADVIA 1650** (x) using 28 samples gave following results:

$$y = 0.908 x + 0.061 \text{ mg/dl;}$$

$$R = 0.990 \quad (R - \text{correlation coefficient})$$

A comparison between creatinine values determined at **Biolis 24i Premium** (y) and at **COBAS INTEGRA** (x) using 43 samples gave following results:

$$y = 0.925 x + 0.006 \text{ mg/dl;}$$

$$R = 0.991 \quad (R - \text{correlation coefficient})$$

TRACEABILITY

CREATININE STANDARD 2 is traceable to the SRM 909B reference material.

WASTE MANAGEMENT

Please refer to local legal requirements.

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