LIPASE

	(PL)
Nazwa zestawu	Nr kat.
CORMAY LIPASE	1-309
HC-LIPASE	4-556
OS-LIPASE	9-426
B50-LIPASE	5-525

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LITERATURA

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

Data wydania: 07. 2023.

LIPASE

Kit name	(EN)
CORMAY LIPASE	Cat. No
HC-LIPASE	1-309
OS-LIPASE	4-556
B50-LIPASE	9-426
	5-525

INTENDED USE

Diagnostic kit for determination of lipase activity, used both for manual assay and in several automatic analysers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lipase is a digestive enzyme released into the intestine from the pancreas where it breaks down triglycerides into fatty acids and glycerol prior to absorption. Lipase measurements are used in the diagnosis and treatment of diseases of the pancreas such as acute pancreatitis, obstruction of the pancreatic duct and pancreatic tumours.

METHOD PRINCIPLE

The colorimetric method is based on a lipase specific degradation of a chromogenic substrate. The specific lipase substrate-DGGM [1,2-o-dilauryl-racglycero-3-glutaric acid-(6'-methylresorufin) ester] is cleaved by the catalytic action of lipase to form 1,2-o-dilauryl-racglycerol and an unstable intermediate, glutaric acid-(6-methyl resorufin) ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The lipase activity in the specimen is proportional to the production of methylresorufin in the reaction and can be determined photometrically.

REAGENTS

Package	CORMAY	HC-	OS-	B50-
	LIPASE	LIPASE	LIPASE	LIPASE
1-REAGENT	4 x 25 ml	4 x 88 ml	2 x 26 ml	2 x 35 ml
2-REAGENT	2 x 25 ml	4 x 49 ml	2 x 14 ml	2 x 18.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

Concentrations in the test

1-REAGENT	
TAPS [N-Tris(hydroxymethyl)methyl-3-aminopropanesulfonic acid]	100 mM
sodium hydroxide	40 mM
sodium deoxycholate	34 mM
2-REAGENT	
tartaric acid	9.5 mM
sodium hydroxide	19 mM
colipase	460 IU/ml
2-propanol	0.65 M
DGGM [1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)-ester]	0.4 mM

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- EUH210 Safety data sheet available on request.

SPECIMEN

Serum, heparinized plasma free from hemolysis. Sample may be stored for up to 5 days at 2-8°C or 24 hours at 20-25°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

ADDITIONAL EQUIPMENT

- automatic analyser or photometer able to read wavelength at 570 nm;
- thermostat at 37°C;
- general laboratory equipment;

PROCEDURE

Application guides for analysers are available on request. The application guide should be verified and validated by the user prior to testing patient samples.

Actions required:

When performing assays in analysers, there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT - LIPASE, LDL DIRECT - LIPASE, TG - LIPASE, TGmono - LIPASE, CHOL - LIPASE, CALCIUM ARSENAZO - LIPASE. To minimize this effect, the preventive actions should be taken, such as ordering an extra washing program or performing tests in a separate order.

Manual procedure

wavelength	570 nm
temperature	37°C
cuvette	1 cm

Pipette into the cuvettes:

	test (T)	calibrator (C)
1-REAGENT	800 µl	800 µl

Bring up to the temperature of determination. Then add:

sample	20 µl	-
calibrator	-	20 µl

Mix well and incubate for 5 minutes. Then add:

2-REAGENT	400 µl	400 µl
-----------	--------	--------

Mix well and after 1 min. of incubation at 37°C read the absorbance of calibrator (C) and test (T) against deionized water. After next 1 and 2 minutes repeat absorbance reading and calculate the mean absorbance change per minute for calibrator $\Delta A/\text{min. (C)}$ and for test sample $\Delta A/\text{min. (T)}$.

Calculation

- Calculate the mean absorbance change per minute $\Delta A/\text{min.}$ for the calibrator (C) i.e. for CORMAY MULTICALIBRATOR LEVEL 2. Then plot the calibration curve.

- Calculate the mean absorbance change per minute $\Delta A/\text{min.}$ for each of the test sample $\Delta A/\text{min. (T)}$. Then read the lipase activity by interpolation from the calibration curve.

REFERENCE VALUES ⁴

Normal range	13 – 60 U/l
--------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of manual assay and automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using Multi+ for manual assay and automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.2 U/l (0.003 µkat/l)
- LoD (Limit of Detection):**
0.4 U/l (0.007 µkat/l)
- LoQ (Limit of Quantitation):**
15 U/l (0.25 µkat/l) – Multi+
6 U/l (0.10 µkat/l) – Biolis 30i
- Linearity:**
up to 250 U/l (4.17 µkat/l) - Multi +
up to 330 U/l (5.50 µkat/l) - Biolis 30i

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

Precision (Multi+)

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	64.1	0.83	1.3
level 2	121.5	1.15	0.95

Precision (Biolis 30i)

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	62.6	0.69	1.09
level 2	112.6	0.97	0.86

Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	57.5	2.18	3.8
level 2	104.1	4.38	4.2

Method comparison

A comparison between lipase values determined at **Multi+** (y) and at **ATTELICA SIEMENS 930** (x), using 16 serum samples gave following results:

$$y = 1.2937 x - 10.903 \text{ U/l;}$$

$$R = 0.975 \quad (R - \text{correlation coefficient})$$

A comparison between lipase values determined at **Biolis 30i** (y) and **ADVIA SIEMENS 1800** (x), using 30 serum samples gave following results:

$$y = 1.1131 x - 4.6415 \text{ U/l;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

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