

MAGLUMI Intact PTH (CLIA)



130211001M



100



**Shenzhen New Industries
Biomedical Engineering Co., Ltd**
4F, Wearnes Tech Bldg,
Science & Industry Park,
Nanshan, Shenzhen, 518057 CHINA
Tel. + 86-755-86028224
Fax. + 86-755-26654850



Lotus Global Co., Ltd
15 Alexandra Road
London
NW8 0DP
UK
Tel. + 44-20-75868010
Fax. + 44-20-79006187




FOR PROFESSIONAL USE ONLY
Store at 2-8°C



**COMPLETELY READ THE INSTRUCTIONS BEFORE
PROCEEDING**



SYMBOLS EXPLANATIONS



Authorized Representative in Europe



Manufacturer



Attention. See Instructions For Use



Contents of kit



In vitro diagnostic medical device
(In vitro diagnostic use)




Lot number



Catalogue Code



Expiry date (Use by...)



Temperature limitation
(store at 2-8°C)



Number of tests



Keep away from sunlight



Keep upright

INTENDED USE

The kit has been designed for the quantitative determination of Intact Parathyroid hormone (Intact PTH) in human serum.

The method can be used for samples over the range of 0-5000 pg/ml.

The test has to be performed on the MAGLUMI chemiluminescence immunoassay (CLIA) fully auto analyzer (Including MAGLUMI 1000, MAGLUMI 2000, MAGLUMI 2000 Plus and new developed models).

SUMMARY AND EXPLANATION OF THE TEST

PTH (parathyrin), a single-chain polypeptide (with a molecular mass of approximately 9,500 daltons) containing 84 amino acids, exerts significant influence in the maintenance of optimal calcium ion concentrations. PTH raises serum ionized calcium levels through direct action on bone and the kidneys: it increases the rate of calcium ion flow from bone to the extracellular fluid, and increases both the renal tubular reabsorption of ionized calcium and the renal excretion of phosphate. Long-term regulation of total body calcium by PTH occurs through its stimulation of vitamin D metabolism, which results in enhanced intestinal absorption of ionized calcium.

In healthy individuals, PTH is secreted in response to circulating calcium ion levels. Any dip below an individual's normal level triggers a pronounced increase in PTH secretion. Calcium levels returning to normal exert a negative feedback effect, thus inhibiting PTH secretion by the parathyroid glands.

PTH undergoes proteolysis to a lesser extent in the parathyroid glands but mostly peripherally – especially in the liver but also in the kidneys and bone – to yield N-terminal fragments and longer lived C-terminal and midregion fragments. The N-terminal fragment contains the region that confers bioactivity. C-terminal and N-terminal fragments are initially generated in equivalent amounts, but the N-terminal fragments disappear rapidly. The C-terminal fragment has a half-life of several hours. In renal failure, C-terminal fragment clearance by glomerular filtration is impaired, so that high levels are found. C-terminal assays (as well as midregion assays) are consequently likely to be especially unreliable in chronic renal failure, where increased PTH is typically just a reflection of impaired renal clearance.

For the intact hormone, the in vivo half-life is 2 to 5 minutes. Intact PTH clearance is accomplished by both peritubular uptake and glomerular filtration followed by reabsorption. In normal renal function, intact PTH is the greatest part of circulating PTH-like bioactivity⁴ and is present in the circulation at concentrations of 10-11 to 10-12 mol/l.

In hypercalcemia due to primary hyperparathyroidism or to ectopic PTH production (pseudohyperparathyroidism), the majority of patients have elevated PTH levels. By contrast, in hypercalcemia due to malignancy or other causes, the concentration of PTH in circulation is typically low or within normal reference range limits. PTH levels are also characteristically high in secondary hyperparathyroidism – usually associated with renal failure – as a result of constant stimulation of the parathyroid gland by low calcium levels. Hypocalcemia accompanied by a low PTH level, on the other hand, is to be expected in hypoparathyroidism, either postsurgical or idiopathic.

Immunoassays specific for various PTH fragments have been developed. Most rely on antisera specific for a discrete region: the C-terminal, N-terminal, or midmolecule. The antisera employed in such assays recognize not only the specific region, but similar fragments as well.

Recent assays for intact PTH have the necessary sensitivity for detecting circulating intact PTH in normals and for discriminating between normals and those with primary hyperparathyroidism.

These assays also appear to discriminate better between primary hyperparathyroidism and hypercalcemia of malignancy compared with previous assays, and do so virtually without any significant overlap between these groups.

Much improved clinical sensitivity is reported for PTH assays when dynamic reference intervals (based on a range of serum PTH values obtained by acute modification of serum calcium concentrations in healthy subjects) are used, rather than a gaussian reference range (based on PTH values seen in normocalcemic individuals). Using an immunoradiometric assay for intact PTH and applying a dynamic reference range, Lepage, et al. obtained average clinical sensitivity of up to 100 percent with primary hyper- and hypoparathyroid samples. Moreover, only the intact PTH assay allowed complete separation between primary hyperparathyroid and nonparathyroidal hypercalcemic patients.

PRINCIPLE OF THE TEST

Sandwich immunoluminometric assay;

Use an anti-PTH monoclonal antibody to label ABEI, and use another monoclonal antibody to coat magnetic microbeads. Sample, Calibrator, or Control, with ABEI Label and magnetic microbeads are mixed thoroughly and incubated at 37°C, forming a sandwich; After sediment in a magnetic field, decant the supernatant, then cycle washing for 1 time. Subsequently, the starter reagents are added and a flash chemiluminescence reaction is initiated. The light signal is measured by a photomultiplier as RLU within 3 seconds and is proportional to the concentration of PTH present in controls or samples.



KIT COMPONENTS

Material Supplies

Reagent Integral for 100 determinations	
Nano magnetic microbeads: contains BSA, 0.2%NaN ₃ , coated with anti-PTH monoclonal antibody.	2.5ml
Calibrator low: bovine serum, 0.2%NaN ₃ .	3.0ml
Calibrator high: bovine serum, 0.2%NaN ₃ .	3.0ml
ABEI Label: anti-PTH monoclonal antibody labeled ABEI contains BSA, 0.2%NaN ₃ .	12.5ml
All reagents are provided ready-to-use.	

Reagent Vials in kit box	
Internal Quality Control: containing BSA, 0.2%NaN ₃ . (target value refer to Quality Control Information date sheet)	2.0ml

Accessories Required But Not Provided

MAGLUMI Reaction Module	REF: 630003
MAGLUMI Starter 1+2	REF: 130299004M
MAGLUMI Wash Concentrate	REF: 130299005M
MAGLUMI Light Check	REF: 130299006M



Preparation of the Reagent Integral

Before the sealing is removed, gentle and careful horizontal shaking of the Reagent Integral is essential (avoid foam formation!) Remove the sealing and turn the small wheel of the magnetic microbeads compartment to and fro, until the colour of the suspension has changed into brown. Place the Integral into the reagent area and let it stand there for 30 min. During this time, the magnetic microbeads are automatically agitated and completely resuspended.

Do not interchange integral component from different reagents or lots!

Storage and Stability

- Sealed: Stored at 2-8 °C until the expiry date.
- Opened: Stable for 4 weeks. To ensure the best kit performance, it is recommended to place opened kits in the refrigerator if it's not going to be used on board during the next 12 hours.



- Keep upright for storage.



- Keep away from direct sunlight.

CALIBRATION AND TRACEABILITY

1) Traceability

To perform an accurate calibration, we have provided the test calibrators standardized against the Non WHO Reference Material Parathyroid Hormone, Human-Type, Synthetic (1-34) fragment NIBSC code: 82/508.

2) 2-Point Recalibration

Via the measurement of calibrators, the predefined master curve is adjusted (recalibrated) to a new, instrument-specific measurement level with each calibration.

3) Frequency of Recalibration

- After each exchange of lot (Reagent Integral or Starter Reagents).
- Every 2 weeks and/or each time a new Integral is used (recommendation).
- After each servicing of the MAGLUMI Fully Auto analyzer.
- If controls are beyond the expected range.

SPECIMEN COLLECTION AND PREPARATION

- Sample material: serum
- Collect samples using standard procedures.
- Store at 2-8 °C: 24 hours, for long term storage below -20 °C
- Avoid repeated freezing and thawing cycles, stored samples should be thoroughly mixed prior to use (Vortex mixer).
- Please ask local representative of SNIBE for more details if you have any doubt.

Vacuum Tubes

- Blank tubes are recommended type for collecting samples.
- Please ask SNIBE for advice if special additive must be used in sample collecting.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated specimens;
 - Cadaver specimens or body fluids other than human serum;
 - Obvious microbial contamination.
- Use caution when handling patient specimens to prevent cross contamination. Use of disposable pipettes or pipette tips is recommended.
- Inspect all samples for bubbles. Remove bubbles with an applicator stick prior to analysis. Use a new applicator stick for each sample to prevent cross contamination.
- Serum specimens should be free of fibrin, red blood cells or other particulate matter.
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.

Preparation for Analysis

- Patient specimens with a cloudy or turbid appearance must be centrifuged prior to testing. Following centrifugation, avoid the lipid layer (if present) when pipetting the specimen into a sample cup or secondary tube.
- Specimens must be mixed **thoroughly** after thawing by **low** speed vortexing or by gently inverting, and centrifuged prior to use to remove red blood cells or particulate matter to ensure consistency in the results. Multiple freeze-thaw cycles of specimens should be avoided.
- All samples (patient specimens or controls) should be tested within 3 hours of being placed on board the MAGLUMI System. Refer to the SNIBE service for a more detailed discussion of onboard sample storage constraints.

Storage

- If testing will be delayed for more than 8 hours, remove serum or plasma from the serum or plasma separator, red blood cells or clot. Specimens removed from the separator gel, cells or clot may be stored up to 24 hours at 2-8°C.
- Specimens can be stored up to 30 days frozen at -20°C or colder.

Shipping

Before shipping specimens, it is recommended that specimens be removed from the serum or plasma separator, red blood cells or clot. When shipped, specimens must be packaged and labeled in compliance with applicable state, federal and international regulations covering the transport of clinical specimens and infectious substances. Specimens must be shipped frozen (dry ice). Do not exceed the storage time limitations identified in this section of the package insert.

WARNING AND PRECAUTIONS FOR USERS



- For use in *IN-VITRO* diagnostic procedures only.
- Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Safety Precautions

CAUTION: This product requires the handling of human specimens.

- The calibrators in this kit are prepared from bovine serum products. However, because no test method can offer complete assurance that HIV, Hepatitis B Virus or other infectious agents are absent; these reagents should be considered a potential biohazard and handled with the same precautions as applied to any serum or plasma specimen.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an overkill approach (USP 24, 2000, p.2143). A minimum of one hour at 121 °C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.
- It is recommended that all human sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Blood borne Pathogens 13. Biosafety Level 214 or other appropriate biosafety practices

should be used for materials that contain or are suspected of containing infectious agents.

- This product contains Sodium Azide; this material and its container must be disposed of in a safe way.
- Safety data sheets are available on request.

Handling Precautions

- Do not use reagent kits beyond the expiration date.
- Do not mix reagents from different reagent kits.
- Prior to loading the Reagent Kit on the system for the first time, the microbeads requires mixing to re-suspend microbeads that have settled during shipment.
- For microbeads mixing instructions, refer to the KIT COMPONENTS, Preparation of the Reagent Integral section of this package insert.
- To avoid contamination, wear clean gloves when operating with a reagent kit and sample.
- Over time, residual liquids may dry on the kit surface, please pay attention the silicon film still exists on the surface of the kit.
- For a detailed discussion of handling precautions during system operation, refer to the SNIBE service information.

TEST PROCEDURE

To ensure proper test performance, strictly adhere to the operating instructions of the MAGLUMI Fully Auto analyzer. Each test parameter is identified via a RFID tag on the Reagent Integral. For further information please refer to the MAGLUMI Chemiluminescence Analyzer Operating Instructions.

100µl +100µl +20µl	Sample, calibrator or controls ABEI Label Nano magnetic microbeads
30 min	Incubation
400µl	Cycle washing
3 s	Measurement

DILUTION

Sample dilution by analyzer is not available in this reagent kit.

Samples with concentrations above the measuring range can be diluted manually. After manual dilution, multiply the result by the dilution factor.

Please choose applicable diluents or ask SNIBE for advice before manual dilution must be processed.

QUALITY CONTROL

Observe quality control guidelines for medical laboratories.

Use suitable controls for in-house quality control. Controls should be run at least once every 24 hours when the test is in use, once per reagent kit and after every calibration. The control intervals should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined ranges. Each laboratory should establish guidelines for corrective measures to be taken if values fall outside the range.

LIMITATIONS OF THE PROCEDURE

1) Limitations

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. Bacterial contamination of samples or repeated freeze-thaw cycles may affect the test results. Assay results should be utilized in conjunction with other clinical and laboratory data to assist the clinician in making individual patient management decisions.

2) Interfering Substances

No interference with test results is seen by concentrations of bilirubin <0.06mg/ml, haemoglobin <16mg/dl or triglycerides <12.5mg/ml.

3) HAMA

Patient samples containing human anti-mouse antibodies (HAMA)

may give falsely elevated or decreased values. Although HAMA-neutralizing agents are added, extremely high HAMA serum concentrations may occasionally influence results.

4) High-Dose Hook

No high-dose hook effect was seen for intact PTH concentrations up to 10 pg/ml. If RLU value of the samples is higher than Calibrator high on the curve, the concentration calculated by the instrument is not necessarily accurate. For those samples, it is recommended to dilute them until an RLU value ranging between Standard A and calibrator high is exhibited. After that, send them for a measurement. The output value multiplied by dilution ratio equals the final RLU value of samples.

RESULTS

1) Calculation of Results

The analyzer automatically calculates the PTH concentration in each sample by means of a calibration curve which is generated by a 2-point calibration master curve procedure. The results are expressed in pg/ml. For further information please refer to the MAGLUMI Chemiluminescence Analyzer Operating Instructions.

2) Interpretation of Results

- Results of study in clinical centers with group of individuals, 95% of the results were < 80 pg/ml.
- Results may differ between laboratories due to variations in population and test method. If necessary, each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

1) Precision

Intra-assay coefficient of variation was evaluated on 3 different levels of control serum repeatedly measured 20 times in the same run, calculating the coefficient of variation.

Intra-assay precision			
Control	Mean(pg/ml)	SD(pg/ml)	CV%
Level 1	150.32	8.62	5.74
Level 2	496.15	26.94	5.43
Level 3	2739.72	141.36	5.16

Inter-assay coefficient of variation was evaluated on three batches of kits. Repeatedly measured 3 different levels of control serum 21 times, calculating the coefficient of variation.

Inter-assay precision			
Control	Mean(pg/ml)	SD(pg/ml)	CV%
Level 1	148.36	13.77	9.28
Level 2	502.55	47.69	9.49
Level 3	2716.49	253.18	9.32

2) Analytical Sensitivity

The sensitivity is defined as the concentration of PTH equivalent to the mean RLU of 20 replicates of the zero standard plus two standard deviations corresponding to the concentration from the standard curve. The sensitivity is typically less than 12.5 pg/ml.

3) Specificity

The specificity of the PTH assay system was assessed by measuring the apparent response of the assay to various potentially cross reactive analytes.

Compound	Concentration	Cross reactivity
ACTH	100pg/ml	0.3%

4) Recovery

Consider calibrator high of known concentration as a sample, dilute it by 1:2 ratio with diluents, and measure its diluted concentration for 10 times. Then calculate the recovery of measured concentration and expected concentration. The recovery should be within 90% -110%.

Expected	Mean Measuring	Recovery
1405.85pg/ml	1386.31pg/ml	98%

5) Linearity

Use PTH calibrator to prepare the six-point standard curve, measuring all points' RLU except point A, and then do four-parameter linear fitting in double logarithm coordinate, the absolute linear correlation coefficient(r) should be bigger than 0.9800.

Calibrator Point	Concentration pg/ml	Absolute linear correlation coefficient (r)
A	0	
B	50	r=0.9804
C	200	
D	5000	
E	2000	
F	5000	

6) Method comparison

A comparison of MAGLUMI PTH (y) with a commercially available PTH test (x) using clinical samples gave the following correlations (pg/ml):

Linear regression

$$y = 0.98x + 114.0$$

$$r = 0.969$$

$$S_{y.x} = 28.9$$

Number of samples measured: 200

The sample concentrations were between 50 and 4430 pg/ml.

REFERENCES

1. Armitage EK. Parathyrin (parathyroid hormone): metabolism and methods for assay. Clin Chem 1986;32: 418-24.
2. Health H. Tests of parathyroid function: utility and limitations. American Association for Clinical Chemistry: Continuing Education in Endocrinology and Metabolism. Endo(Feb) 1984;2(8):10 pages
3. Kao PC. Parathyroid hormone assay. Mayo Clin Proc 1982; 57:596-7.
4. Kao PC, Grant CS, Klee GG, Khosla S. Clinical performance of Parathyroid immunometric assays. Mayo Clin Proc 1992; 67: 637-45.
5. Parathyroid hormone-related protein (PTH-rP)-producing lung cancer with psaos abscess-like metastasis and humoral hypercalcemia of malignancy. International Journal of Clinical Oncology, 2000, Volume 5, Number 6, 410-413.
6. Retention prediction and computer-assisted optimization for the separation of PTH-amino acids in isocratic reversed-phase liquid chromatography. Chromatographia, 1988, Volume 25, Number 11, 974-982.
7. Sequence analysis of phosphotyrosine-containing peptides. Determination of PTH-phosphotyrosine by capillary electrophoresis. Chromatographia, 1990, Volume 30, Numbers 11-12, 691-695.
8. 5,000 Parathyroid Operations Without Frozen Section or PTH Assays: Measuring Individual Parathyroid Gland Hormone Production in Real Time. Annals of Surgical Oncology, 2009, Volume 16, Number 3, 656-666