

Dia-WASH

WASHING SOLUTION



Cat. No.: 13180 **12 x 15 ml**

PRODUCT NAME

Dia-WASH washing solution.

INTENDED USE

(For In Vitro Diagnostic Use Only)

Dia-WASH solution is used for daily maintenance and cleaning of automated coagulation analysers for diagnostic measurements.

SUMMARY AND EXPLANATIONS

Dia-WASH solution is an acidic washing solution, which effectively cleans compartments of coagulation analysers.

ACTIVE INGREDIENTS

Dia-WASH solution contains hydrochloric acid (<1,0%).

PRECAUTIONS

- Person installing the Dia-WASH solution must be a trained laboratory professional!
- All reagents, waste and utilized disposable laboratory equipment should be considered as hazardous waste! Their handling and disposal should be done according to the valid hazardous material processing regulation.
- Do not use the reagent beyond the expiration date printed on the label!

PREPARATION

Dia-WASH solution is ready to use. Wait until the solution reaches the working temperature!

STORAGE AND STABILITY

Dia-WASH solution in intact vial is stable until the expiration date given on the vial, when stored at 2-8°C. Stability after opening in the original vial is shown in below table:

T (°C)	20-25	15-19
Day	7	7

EXPECTED RESULTS

Dia-WASH solution will effectively clean the needles in the instrument used on daily basis and keeping carryover values within specification.

MATERIALS REQUIRED BUT NOT PROVIDED

- For the usage of the following tests on automated Coag Line analyser:
 - Dia-PROTEIN C (chromogenic); Cat. No.: 07108
 - Dia-PROTEIN S (clotting); Cat. No.: 07204
- Optical or mechanical automated coagulation analyser for measuring, Diagon analysers (Coag Line) are recommended.

DIAGON LTD.
 Baross u. 48-52, 1047 Budapest, Hungary
 Tel.: +36 1 3696500
 Fax.: +36 1 3696301
 Web: www.diagon.com
 E-mail: diagon@diagon.com

SYMBOLS			
	Manufacturer		Use-by date
	Batch code		Catalogue number
	Do not use if package is damaged		Fragile, handle with care
	Keep dry		Temperature limit
	Caution		Consult instruction for use
	Contains sufficient for <n> tests		In vitro diagnostic medical device
	CE mark		This side up