UA

	(PL)
Nazwa zestawu	Nr kat.
Liquick Cor-UA 30	2-235
Liquick Cor-UA 60	2-208
Liquick Cor-UA 120	2-209
HC-UA	4-508
OS-UA	9-409
B50-UA	5-515

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INTENDED USE

Diagnostic kit with ascorbate oxidase for determination of uric acid concentration used both for manual assay (Sample Start and Reagent Start method) and in several automatic analysers.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

REAGENTS

Package	Liquick Cor - UA 30	Liquick Cor - UA 60	Liquick Cor - UA 120
1-REAGENT	5 x 24 ml	5 x 48 ml	5 x 96 ml
2-REAGENT	1 x 30 ml	1 x 60 ml	1 x 120 ml
3-STANDARD	1 x 2 ml	-	-
	HC-UA	OS-UA	B50-UA
1-REAGENT	6 x 78.5 ml	2 x 56 ml	3 x 58.5 ml
2-REAGENT	6 x 20 ml	2 x 18.5 ml	3 x 17.5 ml

3-STANDARD is a standard solution of uric acid with a concentration within the range 4.545 – 5.555 mg/dl (270.34 – 330.41 µmol/l). The exact concentration is printed on the label of the each vial.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

Working reagent preparation and stability

Assay can be performed with use of separate 1-REAGENT and 2-REAGENT reagents or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-REAGENT with 1 part of 2-REAGENT. Avoid foaming.

Stability of working reagent: 3 months at 2-8°C
2 weeks at 15-25°C

Concentrations in the reagent

1-REAGENT
ascorbate oxidase ≤ 104 µkat/l
peroxidase (POD) ≤ 22.4 µkat/l

Reagent Start method

The determination can be also performed with use of separate 1-REAGENT and 2-REAGENT reagents.

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
1-REAGENT	1000 µl	1000 µl	1000 µl

Bring up to the temperature of determination. Then add:

	-	-	20 µl
standard	-	-	20 µl
sample	-	20 µl	-

Mix well, incubate for 5 min. Then add:

2-REAGENT	250 µl	250 µl	250 µl
-----------	--------	--------	--------

Mix well; perform measurement as described for Sample Start method.

Calculation

$$\text{uric acid concentration} = \frac{A(T)}{A(S)} \times \text{standard concentration}$$

REFERENCE VALUES⁵

serum / plasma	mg/dl	µmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458
24-hours urine	mg/24h	mmol/24h
	250 – 750	1.49 – 4.46

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For calibration when using the manual methods URIC ACID STANDARD 5 (Cat. No 5-125) is recommended.

For calibration of the automatic analysers systems CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using a Multi + analyser for manual assay (Sample Start method) and/or an automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.04 mg/dl (2.38 µmol/l) – Biolis 30i
- LoD (Limit of Detection):**
0.09 mg/dl (5.35 µmol/l) – Biolis 30i
- LoQ (Limit of Quantitation):**
0.3 mg/dl (17.84 µmol/l) – serum/plasma – Biolis30i
0.14 mg/dl (8.33 µmol/l) – urine – Biolis30i
0.6 mg/dl (35.69 µmol/l) – serum/plasma – Multi+
- Linearity:**
up to 40 mg/dl (2379.2 µmol/l) – serum/plasma – Biolis 30i
up to 56 mg/dl (3330.88 µmol/l) – urine – Biolis 30i
up to 36 mg/dl (2141.28 µmol/l) – serum/plasma – Multi+

For higher concentration, in serum or plasma, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 30 mg/dl for determinations in serum, ascorbate up to 50 mg/dl for determinations in urine, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision (Multi+)

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	5.47	0.16	2.8
level 2	9.75	0.29	3.0

Precision (Biolis 30i)

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	4.87	0.06	1.27
level 2	9.35	0.06	0.67
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	4.86	0.07	1.5
level 2	9.19	0.12	1.3

Method comparison

A comparison between uric acid values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 0.973x + 0.386 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between uric acid values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 34 plasma samples gave following results:

$$y = 0.9416x + 0.4715 \text{ mg/dl};$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

A comparison between uric acid values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:

$$y = 1.0126x - 0.1377 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between uric acid values determined at **Multi+(y)** and at **BECKMAN COULTER AU680** (x) using 21 serum samples gave following results:

$$y = 1.0227x + 0.1646$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

TRACEABILITY

URIC ACID STANDARD 5 is traceable to the SRM 1950/ 909C reference material.

WASTE MANAGEMENT

Please refer to local legal requirements.

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- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).