## CORMAY LDL DIRECT

	(PL) Nr kat.	CORMAY LDL DIRECT 120	CORMAY LDL DIRECT 500	CORMAY LDL DIRECT
Nazwa zestawu				
CORMAY LDL DIRECT mini	2-216			
CORMAY LDL DIRECT 30	2-191	1-Reagent 3 x 100 ml	3 x 300 ml	4 x 30 ml
CORMAY LDL DIRECT 60	2-192			
CORMAY LDL DIRECT 120	2-193	2-Reagent 1 x 100 ml	1 x 300 ml	4 x 10 ml
CORMAY LDL DIRECT 500	2-194			
CORMAY LDL DIRECT	2-180			

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### INTENDED USE

Diagnostic kit for determination of LDL-cholesterol concentration (direct method) used both for manual assay (Sample Start and Reagent Start method) and in several automatic analysers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions

### INTRODUCTION

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The relative protein and lipid determine the density of these lipoproteins and provide the basis on which to begin their classification. The classes are: chylomicron, very-low-density lipoprotein (VLDL), low-density-lipoprotein (LDL) and high-density lipoprotein (HDL).

LDL are synthesised in the liver by the action of various lipolytic enzymes on triglyceride rich VLDL. LDL-cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.

Accurate measurement of LDL-cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture.

### METHOD PRINCIPLE

The assay is a homogenous method for directly measuring LDL-cholesterol concentrations in serum or plasma, without the need for any off-line pretreatment or centrifugation steps. Liquid selective detergent method.

The method is in a two reagents format and depends on the properties of a unique detergent. This detergent (Reagent 1) solubilizes only the non LDL particles (HDL, VLDL, CM). The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction. A second detergent (Reagent 2) solubilizes the remaining LDL particles and chromogenic coupler allows for color formulation. The enzyme reaction with LDL-cholesterol in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

### REAGENTS

<b>Package</b>	<b>CORMAY LDL DIRECT mini</b>	<b>CORMAY LDL DIRECT 30</b>	<b>CORMAY LDL DIRECT 60</b>
1-Reagent	2 x 30 ml	3 x 30 ml	3 x 50 ml
2-Reagent	2 x 10 ml	1 x 30 ml	1 x 50 ml
	<b>CORMAY LDL DIRECT 120</b>	<b>CORMAY LDL DIRECT 500</b>	<b>CORMAY LDL DIRECT</b>
1-Reagent	3 x 100 ml	3 x 300 ml	4 x 30 ml
2-Reagent	1 x 100 ml	1 x 300 ml	4 x 10 ml

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Concentrations in the test

<b>1-Reagent</b>	
Buffer	< 1.0 %
Detergent 1	< 1500 U/l
Cholesterol esterase ( <i>Pseudomonas sp.</i> )	< 1500 U/l
Cholesterol oxidase ( <i>Cellulomonas sp.</i> )	< 1500 U/l
Peroxidase (horseradish)	< 1300 ppg U/l
4-aminoantipyrine	< 0.1 %
Ascorbic Acid Oxidase ( <i>Curcubita sp.</i> )	< 3000 U/l
Preservative	
<b>2-Reagent</b>	
Buffer	
Detergent 2	< 1.0 %
N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT)	< 1.0 mM
Preservative	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.

### SPECIMEN

Serum, sodium heparinized or EDTA plasma. Anticoagulants containing citrate should not be used. **Serum:** Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).

**Plasma:** Specimens may be collected in EDTA or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).

If not analysed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens need to be stored for more than 5 days, they may be preserved at -80°C. Samples may be frozen once.

Nevertheless it is recommended to perform the assay with freshly collected samples.

### ADDITIONAL EQUIPMENT

- automated analyser or photometer able to read at 630 nm;
- thermostat at 37°C;
- general laboratory equipment;

### PROCEDURE

The reagents are ready to use. Applications for automatic analysers are available on request.

### Manual procedure

wavelength	630 nm
temperature	37°C
cuvette	1 cm

Pipette into the cuvettes:

	test (T)	standard (S)
1-Reagent	900 µl	900 µl
Bring up to the temperature of determination (37 °C). Then add		
calibrator	-	9 µl
sample	9 µl	-

Mix well and after 5 min. of incubation at 37°C add:

2-Reagent	300 µl	300 µl
-----------	--------	--------

Mix well and after exactly 5 min. of incubation read the absorbance of standard A(S) and test A(T) against air or water.

### Calculation

$$\text{LDL cholesterol concentration} = \frac{A(T)}{A(S)} \times \text{calibrator concentration}$$

### REFERENCE VALUES<sup>7</sup>

	Adults	
NCEP* Classification:	mg/dl	mmol/l
Optional	< 100	< 2.59
Near optional	< 130	< 3.37
Bordeline high	130-159	3.37-4.12
High	160-189	4.14-4.90
Very high	≥ 190	≥ 4.92

\* National Cholesterol Education Program

As LDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2 (Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended.

The calibration curve should be prepared every 12, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using Multi+ for manual assay and automatic analyser Biolis 24i Premium. Results may vary if a different instrument is used.

#### Sensitivity:

5.5 mg/dl (0.14 mmol/l) – Multi+  
4.2 mg/dl (0.11 mmol/l) – Biolis 24i Premium

#### Linearity:

up to 316 mg/dl (8,18 mmol/l) – Multi+  
up to 700 mg/dl (18,13 mmol/l) - Biolis 24i Premium

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Triglycerides up to 1293 mg/dl, bilirubin conjugated up to 20 mg/dl, bilirubin total up to 20 mg/dl, haemoglobin up to 0.5 g/dl, ascorbic acid up to 500 mg/l and gamma-globulins up to 5000 mg/dl do not interfere with the test.

N-acetyl-p-benzoquinone imine (NAPQI), a metabolite of paracetamol (acetaminophen), may cause falsely low results for patients with a toxic level of paracetamol.

#### Precision (Multi+)

Repeatability (run to run) n = 10	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	62.72	1.22	1.94
level 2	158.77	0.69	0.44

#### Precision (Biolis 24i Premium)

Repeatability (run to run) n = 10	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	135.60	4.35	3.21
level 2	187.07	5.19	2.77
Reproducibility (day to day) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	135.06	3.67	2.71
level 2	187.24	5.62	3.00

#### Method comparison

A comparison between LDL cholesterol values determined at Multi+ (y) and at Hitachi 912 (x) using 48 samples gave following results:

$$y = 0.8968x + 11.543 \text{ mg/dl};$$

$$R = 0.982 \quad (R - \text{correlation coefficient})$$

A comparison between LDL cholesterol values determined at Biolis 24i Premium (y) and at Hitachi 912 (x) using 29 samples gave following results:

$$y = 1.0365x + 3.7224 \text{ mg/dl};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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