## BIL TOTAL MALLOY-EVELYN

Nazwa zestawu	(PL)	Nr kat.
Liquick Cor-BIL TOTAL MALLOY-EVELYN 30	2-343	
Liquick Cor-BIL TOTAL MALLOY-EVELYN 60	2-344	
Liquick Cor-BIL TOTAL MALLOY-EVELYN 120	2-345	

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### INTENDED USE

Diagnostic kit for determination of total bilirubin concentration intended to use for manual assay and in several automatic analyzers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Bilirubin is a yellow pigment – product of heme degradation. For clinical purposes, bilirubin is expressed as two fractions: conjugated and unconjugated. In hepatocytes bilirubin is enzymatically conjugated with glucuronic acid residues. This form is called direct or conjugated. Bilirubin without glucuronic acid modification is bound to albumin and is termed unconjugated or indirect. Indirect bilirubin is calculated

as the difference between total and direct bilirubin.

Hiperbilirubinemia is usually the result of jaundice (mechanical, hemolytic), Dubin-Jonson syndrome, Gilbert's syndrome, Crigler-Najjar syndrome, bile ducts disease.

### METHOD PRINCIPLE

Bilirubin and bilirubin glucuronate react with sulphodiazonium salt and form coloured derivative – azobilirubin. Bilirubin glucuronate is soluble in water and reacts directly, bilirubin associated with albumin must be previously hydrolysed with detergents. The colour intensity of formed azobilirubin is proportional to total bilirubin concentration in the sample.

### REAGENTS

#### Package

	Liquick Cor-BIL TOTAL MALLOY-EVELYN 30	Liquick Cor-BIL TOTAL MALLOY-EVELYN 60	Liquick Cor-BIL TOTAL MALLOY-EVELYN 120
1-REAGENT	5 x 24 ml	5 x 48 ml	5 x 96 ml
2-REAGENT	1 x 30 ml	1 x 60 ml	1 x 120 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 10 weeks on board the analyser at 2-10°C.

#### Working reagent preparation and stability

Assay can be performed with use of separate 1-REAGENT and 2-REAGENT reagents or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-REAGENT with 1 part of 2-REAGENT. Avoid foaming!

Stability of working reagent: 7 days at 2-8°C  
1 day at 15-25°C

Concentrations in the test	
sulphanilic acid	25.6 mmol/l
hydrochloric acid	40 mmol/l
sodium nitrite	1 mmol/l
detergent	49.6 mmol/l

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Ingredients:

1-REAGENT contains hydrochloric acid.

#### Danger

H314 Causes severe skin burns and eye damage.

EUH208 Contains sulfanilic acid. May produce an allergic reaction.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

### ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 530 nm (Hg 546, 550 nm);
- thermostat at 37°C;
- general laboratory equipment;

### SPECIMEN

Serum free from hemolysis.

Serum should be separated from red blood cells as soon as possible after blood collection.

Lipemic specimens may show falsely increased bilirubin concentration thus fasting specimen is recommended.

Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight.

Serum can be stored up to 3 days at 2-8°C or up to 3 months at -70°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Applications for them are available on request.

### Manual procedure

wavelength	530 nm (Hg 546, 550 nm)
temperature	37°C
cuvette	1 cm

### Sample Start method

Pipette into the cuvettes:

	test (T)	test blank (TB)	standard (S)	standard blank (SB)
1-REAGENT	-	1000 µl	-	1000 µl
working reagent	1000 µl	-	1000 µl	-

Bring up to the temperature of determination. Then add:

sample	50 µl	50 µl	-	-
calibrator	-	-	50 µl	50 µl

Mix well, incubate for 5 min. at 37°C. Then read absorbance of test (T) against test blank (TB) and standard (S) against standard blank (SB). The colour is stable for 30 min.

### Calculation

$$\text{total bilirubin concentration} = \frac{A(T)}{A(S)} \times \text{calibrator concentration}$$

### Reagent Start method

The determination can be also performed with use of separate 1-REAGENT and 2-REAGENT reagents.

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
1-REAGENT	1000 µl	1000 µl	1000 µl
calibrator	-	-	100 µl
sample	-	100 µl	-
distilled water	100 µl	-	-

Bring up to the temperature of determination. Then add:

2-REAGENT	250 µl	250 µl	250 µl
-----------	--------	--------	--------

Mix well and read the absorbance A1 of test (T) and standard (S) against reagent blank (RB). Then add:

2-REAGENT	250 µl	250 µl	250 µl
-----------	--------	--------	--------

Mix well and after 5 min. of incubation read the absorbance A2 of test (T) and standard (S) against reagent blank (RB). The intensity of colour is stable for 30 minutes. Calculate  $\Delta A (A2 - A1)$  for the test and standard.

### Calculation

$$\text{total bilirubin concentration} = \frac{\Delta A(T)}{\Delta A(S)} \times \text{calibrator concentration}$$

### REFERENCE VALUES<sup>3</sup>

serum (adults)	0.3 – 1.2 mg/dl 5 – 21 µmol/l
----------------	----------------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration is recommended the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177).

The calibration curve should be prepared every 10 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results characteristics have been obtained using automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 0.05 mg/dl (0.855 µmol/l)
- Linearity:** up to 25 mg/dl (428 µmol/l).

For higher bilirubin concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.08 g/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.09	0.01	0.84
level 2	5.57	0.03	0.59
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.10	0.01	0.95
level 2	5.55	0.06	1.09

### Method comparison

A comparison between bilirubin total values determined at **Biolis 24i Premium** (y) and at **ADVIA 1650** (x) using 46 samples gave following results:

$y = 1.2164x - 0.057$  mg/dl;

$R = 0.9995$  (R - correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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