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Data wydania: 04. 2022.

PRESTIGE 24i ALKALINE WASHING SOLUTION

Reagent name	Variant	(EN) Cat.No
PRESTIGE 24i ALKALINE WASHING SOLUTION	500 ml	3-104

INTENDED USE

PRESTIGE 24i ALKALINE WASHING SOLUTION is intended to use for cleaning the cuvettes of the following biochemistry analyzers: Biolis 24i Premium and Biolis 30i.

It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

PRESTIGE 24i ALKALINE WASHING SOLUTION is washing solution connected to the analyzer. Removal of biological material and reagents residues, especially protein-containing materials and microbial deposits from the cuvettes takes place due to the cleaning properties of the reagent. PRESTIGE 24i ALKALINE WASHING SOLUTION contains sodium hydroxide and an enzyme that dissolve and remove organic particles.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

sodium hydroxide	1.25 mol/l
SAVINASE	0.06 %
surfaktant	
preservative	

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.

Once opened, diluted reagent is stable for 30 days at room temperature.

WARNING AND NOTES

- Protect from direct sunlight!
- Do not use the reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- The product meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

PRESTIGE 24i ALKALINE WASHING SOLUTION contains sodium hydroxide and reaction mass of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

Danger

- H314 Causes severe skin burns and eye damage.
- H318 Causes serious eye damage.
- H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.
H290 May be corrosive to metals.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

P273 Avoid release to the environment.

P234 Keep only in original packaging.

ASSAY PROCEDURE

PRESTIGE 24i ALKALINE WASHING SOLUTION should be used as 2% solution. Before use, mix 2 parts of the washing solution with 98 parts of distilled water.

Use according to the analyzer's operation manual.

WASTE MANAGEMENT¹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements

- Reagents residues: 18 01 06*
- Empty packages: 15 01 10*
- Wastewater from the analyzer: 18 01 03*

INCIDENTS²

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (websites address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, a serious public health threat.

LITERATURE

- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

LIST OF CHANGES

Previous version: 05	Current version: 06
Sections updated: REAGENT STABILITY; WARNING AND NOTES	

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