ACCENT-200 ALKALINE WASHING SOLUTION

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INTENDED USE

ACCENT-200 ALKALINE WASHING SOLUTION is intended to use for cleaning the reagent and sample probes of the biochemistry analyzers: ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180 and BS-120 during routine analyzer operation. It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

ACCENT-200 ALKALINE WASHING SOLUTION is a washing solution placed on board the analyzer. Removal of biological material and reagents residues, especially protein-containing materials and microbial deposits from the reagent and sample probes takes place due to the cleaning properties of the reagent. ACCENT-200 ALKALINE WASHING SOLUTION contains sodium hydroxide and sodium hypochlorite that dissolve and remove organic particles.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

sodium hydroxide	10,3 mmol/l
sodium hypochlorite	<1%
carbonate buffer stabilizer	

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.

WARNING AND NOTES

- Protect from direct sunlight!
- Do not use the reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- The product meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

ACCENT-200 ALKALINE WASHING SOLUTION contains sodium hypochlorite.

Danger

- H314 Causes severe skin burns and eye damage.
- H318 Causes serious eye damage.
- H412 Harmful to aquatic life with long lasting effects.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

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