CONCENTRATED DETERGENT B

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Data wydania: 10. 2022.

CONCENTRATED DETERGENT B

Reagent name	Variant	Cat.No
CONCENTRATED DETERGENT B	12 x 1 L	3-117
CONCENTRATED DETERGENT B	1 x 1 L	3-118

INTENDED USE

CONCENTRATED DETERGENT B is intended to use for cleaning the cuvettes of the following biochemistry analyzers: ACCENT MC240 / BS-240Pro, ACCENT M320 / BS-360E, ACCENT 800 / BS-800M, BS-400, BS-480 and BS-800. It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

CONCENTRATED DETERGENT B is a washing solution connected to the analyzer. Removal of biological material and reagents residues, especially protein-containing materials and microbial deposits from the cuvettes takes place due to the cleaning properties of the reagent. CONCENTRATED DETERGENT B contains surfactants and sodium hydroxide that dissolve and remove organic particles.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

Thesit	1.25 %
citrate buffer	0.2 mol/l
sodium hydroxide surfactant	0.7 mol/l

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.

Once opened, concentrated reagent is stable for 30 days at room temperature.

Once opened, diluted reagent is stable for 30 days at room temperature.

WARNING AND NOTES

- Protect from direct sunlight!
- Do not use the reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- The product meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

CONCENTRATED DETERGENT B contains sodium hydroxide.

Danger

H314 Causes severe skin burns and eye damage.
 H318 Causes serious eye damage.
 H290 May be corrosive to metals.
 P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

P234 Keep only in original packaging.

ASSAY PROCEDURE

ANALYZERS: **BS-400, BS-480, BS-800, ACCENT 800 / BS-800M** Reagent is ready to use.

ANALYZERS: **ACCENT MC240/BS-240Pro, ACCENT M320/BS-360E:**

THE REAGENT MUST BE DILUTED BEFORE USE!
 USE AS 10% SOLUTION

Preparation of a diluted washing solution:

1. Pour the concentrated washing solution (1l) into the container for the diluted washing solution.
2. Add deionized water to the container of the diluted washing solution to obtain the appropriate dilution factor (1 liter of washing liquid : 9 liters of deionized water).
3. Replace the container cap and shake the container gently to thoroughly mix the liquids with each other.

Use according to the analyzer's operation manual.

WASTE MANAGEMENT ¹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements

- Reagents residues: 18 01 06*
- Empty packages: 15 01 10*
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ²

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (websites address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

1. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
2. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

LIST OF CHANGES

Previous version: 03	Current version: 04
<i>Sections updated: WARNING AND NOTES</i>	

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