PHOSPHORUS

Nazwa zestawu	(PL) Nr kat.
Liquick Cor-PHOSPHORUS 30	3-243
Liquick Cor-PHOSPHORUS 120	3-338
HC – PHOSPHORUS	4-543
OS – PHOSPHORUS	9-415
B50 - PHOSPHORUS	5-534

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INTENDED USE

PHOSPHORUS reagent is intended to determine quantitatively phosphorus levels in plasma, serum and urine. It is intended to monitor and as an aid to the diagnosis of clinical conditions associated with abnormal phosphorus levels. PHOSPHORUS reagent is intended to use on automatic analysers and manually.

2-STANDARD is intended to calibrate manual determinations of phosphorus, performed with Liquick Cor-PHOSPHORUS 30 reagent variant.

It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY 1-4

Determinations of phosphorus levels in plasma, serum and urine are used to monitor and as an aid to the diagnosis of parathyroid and kidney diseases or vitamin D imbalance.

Increased phosphorus levels may be associated with increased phosphate intake, cell lysis (leukemia), transcellular phosphate shift (acidosis).

Decreased phosphorus levels may occur in refeeding syndrome, decreased net intestinal phosphate absorption, lowered phosphate threshold in the kidneys.

METHOD PRINCIPLE 5,6

Direct phosphomolybdate reaction without deproteinization. Phosphate ions form with molybdate ions in acid solution proportional amounts of unreduced phosphomolybdate complex. The concentration of the complex formed is determined by measuring its absorbance.

REAGENTS

Package

	Liquick Cor-PHOSPHORUS 30	Liquick Cor-PHOSPHORUS 120	HC – PHOSPHORUS	OS - PHOSPHORUS S	B50 – PHOSPHORUS S
1-REAGENT	6 x 30 ml	6 x 100 ml	6 x 97 ml	3 x 53 ml	2 x 58.5 ml
2-STANDARD	1 x 2 ml	-			

2-STANDARD is a standard solution of phosphate ions with a concentration within the range 4.50–5.50 mg/dl (1.46–1.78 mmol/l). The exact concentration is printed on the label of the each vial.

CONCENTRATION OF THE ACTIVE INGREDIENTS

1-REAGENT	
ammonium molybdate	0.4 mmol/l
sulphuric acid	
hydrochloric acid	
detergents	
2-STANDARD	
potassium dihydrogen phosphate	5 mg/dl
hydrochloric acid	

REAGENT STABILITY

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (Biolis 30i).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. Disposable plastic ware is recommended for the test.
- The reagent is usable when its absorbance is less than 0.200 (read against distilled water, wavelength $\lambda=340$ nm, cuvette l = 1 cm, at temp. 37°C).
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-STANDARD meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains sulfuric acid (VI) and hydrochloric acid.

2-STANDARD contains hydrochloric acid.

Danger

 H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER or doctor.

EQUIPMENT REQUIRED BUT NOT PROVIDED

- automatic analyzer or photometer able to read at 340 nm;
- thermostat at 37°C;
- general laboratory equipment;

SPECIMEN 7,8,9

Serum, heparinized plasma free from hemolysis, 24-hours urine.

Serum is the preferred specimen. Level of inorganic phosphate in heparinized plasma is about 0.2 - 0.3 mg/dl (0.06 – 0.10 mmol/l) lower than in serum.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain several times higher phosphate concentration than normal serum.

Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Urine preparation: For collection of 24-hour urine specimens, 20-30 ml of 6M HCl should be added into the container to avoid phosphate precipitations. Dilute 1 part of acidified urine with 9 parts of distilled water before determination. Multiply the result by the dilution factor. 24-hours urine samples can be stored up to 2 days at 20-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE


Reagents are ready to use.

Carry-over effect may occur when the reagent is used on automated chemistry analyzers. It is based on interference between specific test combinations and results in either underestimating or inflating patient sample results. To minimize this effect, the preventive actions should be taken, such as ordering an extra washing program or performing tests in a separate order.

General procedure for automatic determinations:

R1: 300 μ l

Sample: 3 μ l

Main Wavelength: 340 nm

2nd wavelength: 700 nm

Method: End Point

Calibration Type: Linear

Direction: Increasing

R1 Volume: 300 μ l

Sample Volume: 3 μ l

Application guides for analysers are available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

Manual procedure

wavelength	340 nm
temperature	37°C
cuvette	1 cm

Pipette into the cuvette:

	reagent blank (RB)	test (T)	standard (S)
1-REAGENT	1000 μ l	1000 μ l	1000 μ l

Bring up to the temperature of determination. Then add:

2-STANDARD	-	-	10 μ l
sample	-	10 μ l	-

Mix well, incubate for 5 min. at the determination temperature. Read the absorbance of test A(T) and standard A(S) against reagent blank (RB). The absorbance is stable within 30 minutes.

Calculation

$$\text{phosphorus concentration} = \frac{A(T)}{A(S)} \times \text{standard concentration}$$

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum, CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of manual assay 2-STANDARD is recommended.

For calibration of the automatic analysers systems CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

Calibration stability depends on type of analyser used for analysis. The calibration curve should be prepared every 12 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES 7

serum / plasma	mg/dl	mmol/l
age:		
0 – 9 d	4.5 – 9.0	1.45 – 2.91
10 d - 2 y	4.0 – 6.5	1.29 – 2.10
3 – 9 y	3.2 – 5.8	1.03 – 1.87
10 – 15 y	3.3 – 5.4	1.07 – 1.74
16 – 59 y	2.4 – 4.4	0.78 – 1.42
60 - 89 y males	2.3 – 3.7	0.74 - 1.20
60 - 89 y females	2.8 – 4.0	0.90 – 1.26
>89 y females	2.5 - 4.2	0.81 - 1.36
>89 y males	2.2 - 3.9	0.71 - 1.26

24-hours urine	g/24h	mmol/24h
	0.4 – 1.3	12.9 – 42.0

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

Phosphorus concentration in 24-hours urine – calculation

phosphorus concentration in 24-hours urine [g/24h]	=	phosphorus concentration of 24-hours urine [mg/dl]	x	urine volume of 24-hours urine [dl/24h]	÷	1000
--	---	--	---	---	---	------

PERFORMANCE CHARACTERISTICS

The following results have been obtained using a Multi+ analyser for manual assay and/or an automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**
0.04 mg/dl (0.013 mmol/l) - Biolis 30i
- **LoD (Limit of Detection):**
0.07 mg/dl (0.023 mmol/l) - Biolis 30i
- **LoQ (Limit of Quantitation):**
0.4 mg/dl (0.129 mmol/l) - Multi +
0.13 mg/dl (0.042 mmol/l) - Biolis 30i
- **Linearity:**
up to 18 mg/dl (5.8 mmol/l) – Multi +
up to 20 mg/dl (6.46 mmol/l) - Biolis 30i

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- **Measurement range**
0.4 mg/dl (0.129mmol/l) - 18 mg/dl (5.8mmol/l) –Multi +
0.13 mg/dl (0.042 mmol/l) - 20 mg/dl (6.46 mmol/l) - Biolis 30i
- **Specificity / Interferences**
Haemoglobin up to 0.12 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

▪ Precision (Multi+)

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	3.24	0.11	3.4
level 2	7.18	0.20	2.8

▪ Precision (Biolis 30i)

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	3.36	0.08	2.4
level 2	7.31	0.11	1.6

Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	3.31	0.07	2.0
level 2	7.27	0.09	1.3

▪ Method comparison

A comparison between phosphorus values determined at **Multi+** (y) and at **BECKMAN COULTER AU680** (x) using 30 serum samples gave following results:
y = 1.03 x - 0.0731 mg/dl;
R = 0.995 (R – correlation coefficient)

A comparison between phosphorus values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 72 serum samples gave following results:
y = 0.989 x + 0.0358 mg/dl;
R = 0.995 (R – correlation coefficient)

A comparison between phosphorus values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 50 urine samples gave following results:
y = 0.9804 x - 0.7055 mg/dl;
R = 0.998 (R – correlation coefficient)

TRACEABILITY

2-STANDARD is traceable to the spectrophotometric reference method.

WASTE MANAGEMENT ¹⁰

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents: 18 01 06*
- Empty packages: 15 01 10*
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ¹¹

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cornay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

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LIST OF CHANGES

Previous version: 05	Current version: 06
Sections updated: <i>WARNING AND NOTES; WASTE MANAGEMENT</i>	

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