## CALCIUM

Nazwa zestawu	Nr kat.
Liquick Cor-CALCIUM 30	3-250
Liquick Cor-CALCIUM 60	3-251
Liquick Cor-CALCIUM 120	3-252
HC-CALCIUM	4-551
OS-CALCIUM	9-413

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### INTENDED USE

Diagnostic kit for determination of calcium concentration used both for manual assay (Sample Start and Reagent Start method) and in several automatic analysers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Calcium and phosphorus as a hydroxyapatite constitute mineral portion of bone. Calcium occurs also as divalent cations (free or bound with negatively charged proteins) which participate in blood coagulation, neuromuscular excitability, skeletal and cardiac muscle contractility and in multiple cellular functions. Calcium flux in organism is controlled by action of parathyroid hormone (PTH), vitamin D and calcitonin. Calcium serum level abnormalities are caused usually by parathyroid or thyroid disease, disorders of vitamin D metabolism or acute pancreatitis.

### METHOD PRINCIPLE

Calcium ions form a violet complex with o-cresolphthalein complexone in alkaline solution. The intensity of violet colour of this complex measured at 570-580 nm is proportional to the calcium concentration in the sample.

### REAGENTS

Package	Liquick Cor-CALCIUM 30	Liquick Cor-CALCIUM 60	Liquick Cor-CALCIUM 120
1-REAGENT	5 x 24 ml	5 x 48 ml	5 x 96 ml
2-REAGENT	1 x 30 ml	1 x 60 ml	1 x 120 ml
3-STANDARD	1 x 2ml	1 x 2 ml	-

	HC-CALCIUM	OS-CALCIUM
1-REAGENT	6 x 76 ml	2 x 56 ml
2-REAGENT	6 x 19,5 ml	2 x 18,5 ml

3-STANDARD is a standard solution of calcium ions with a concentration within the range 9.00 – 11.00 mg/dl (2.25 – 2.75 mmol/l). Exact concentration is printed on the label of the each vial.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 3 weeks on board the analyser at 2-10°C.

### Working reagent preparation and stability

Assay can be performed with use of separate 1-REAGENT and 2-REAGENT reagents or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-REAGENT with 1 part of 2-REAGENT. Avoid foaming.

Stability of working reagent:	7 days at 2-8°C
	2 days at 15-25°C

Concentrations in the test	(EN)
o-cresolphthalein complexone	0.06 mmol/l
8-quinolinol	8.6 mmol/l
hydrochloric acid	30 mmol/l
ethanolamine	377 mmol/l


### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The working reagent is slightly pink-coloured. This does not influence the results of the test.
- The reagents are usable when the absorbance of the working reagent is less than 0.250 (read against distilled water, wavelength  $\lambda=575$  nm, cuvette l=1 cm, at temp. 25°C).
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2 M HCl solution and then thoroughly rinsed with distilled water.
- 1-REAGENT, 2-REAGENT and 3-STANDARD meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT and 3-STANDARD contain hydrochloric acid.




### Danger

-  H314 Causes severe skin burns and eye damage.
-  P280 Wear protective gloves, protective clothing, eye protection or face protection.
-  P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
-  P310 Immediately call a POISON CENTER or doctor.

### Ingredients:

2-REAGENT contains 2-aminoethanol.

### Danger

-  H314 Causes severe skin burns and eye damage.
-  H335 May cause respiratory irritation.
-  P280 Wear protective gloves, protective clothing, eye protection or face protection.
-  P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
-  P310 Immediately call a POISON CENTER or doctor.
-  P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

### ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 575 nm (570-580 nm);
- thermostat at 37°C;
- general laboratory equipment;

### SPECIMEN

Serum, heparinized plasma free from hemolysis, 24-hours urine. Recommended anticoagulants: heparine lithium, sodium or ammonium salt.

Urine preparation: to prevent calcium salt precipitation in urine, specimens should be collected in HCl, 20-30 ml of 6 M for 24-h specimen. Any specimens collected without acid should be acidified using 20-30 ml of 6 M HCl, well mixed and allowed to stand for 1 h before aliquotting. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor. Serum and plasma can be stored up to 8 hours at 15-25°C or up to 1 day at 2-8°C. Samples frozen at -20°C can be stored up to 1 year. 24-hours urine samples should be kept at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Applications for analyzers are available on request.

### Manual procedure

wavelength	575 nm (570-580 nm)
temperature	20-25°C / 37°C
cuvette	1 cm

### Sample Start method

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
working reagent	1000 $\mu$ l	1000 $\mu$ l	1000 $\mu$ l

Bring up to the temperature of determination. Then add:

standard / calibrator	-	-	10 $\mu$ l
sample	-	10 $\mu$ l	-

Mix well, incubate for 10 min. at 20-25°C or 5 min. at 37°C. Read the absorbance of the test A(T) and standard A(S) against reagent blank (RB). The intensity of colour is stable for 30 min.

### Reagent Start method

The determination can be also performed with use of separate 1-REAGENT and 2-REAGENT reagents.

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
1-REAGENT	1000 $\mu$ l	1000 $\mu$ l	1000 $\mu$ l

Bring up to the temperature of determination. Then add:

standard / calibrator	-	-	10 $\mu$ l
sample	-	10 $\mu$ l	-

Mix well, incubate for 5 min. Then add:

2-REAGENT	250 $\mu$ l	250 $\mu$ l	250 $\mu$ l
-----------	-------------	-------------	-------------

Mix well; perform measurement as described for Sample Start method.

### Calculation

$$\text{calcium concentration} = \frac{A(T)}{A(S)} \times \text{standard / calibrator concentration}$$

### REFERENCE VALUES \*

serum, plasma	mg/dl	mmol/l
premature	6.2 – 11.0	1.55 – 2.75
adults	8.6 – 10.0	2.15 – 2.50
18 – 60 yr	8.8 – 10.2	2.20 – 2.55
60 – 90 yr	8.2 – 9.6	2.05 – 2.40
> 90 yr		
24-hours urine	mg/24h	mmol/24h
	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine with each batch of samples.

For the calibration of manual assay the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176), LEVEL 2 (Cat. No 5-175; 5-177) or CALCIUM STANDARD (Cat. No 5-132) is recommended.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and/or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 2 days, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 0.27 mg/dl (0.068 mmol/l).
- Linearity:** up to 15 mg/dl (3.75 mmol/l).

For higher calcium concentrations dilute the sample with 0.9% NaCl in the ratio of 1 to 1 and reassay. Multiply the result by 2.

### Specificity / Interferences

Haemoglobin up to 2.5 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	8.73	0.08	0.90
level 2	11.97	0.06	0.52
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	8.40	0.25	2.98
level 2	11.46	0.36	3.16

### Method comparison

A comparison between calcium values determined at Biolis 24i Premium (y) and at COBAS INTEGRA 400 (x) using 60 samples gave following results:

$$y = 0.9946x - 0.1547 \text{ mg/dl;}$$

$$R = 0.9962 \quad (R - \text{correlation coefficient})$$

### TRACEABILITY

CALCIUM STANDARD is traceable to the SRM 909B reference material.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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