PRESTIGE 24i LQ LIPASE

Nr kat. **4-185, 4-329** (PL)

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LITERATURA

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

Data wydania: 04.2021.

PRESTIGE 24i LQ LIPASE

Cat. No **4-185, 4-329** (EN)

INTENDED USE

Diagnostic kit for determination of lipase activity intended to use in automatic analysers: analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lipase is a digestive enzyme released into the intestine from the pancreas where it breaks down triglycerides into fatty acids and glycerol prior to absorption. Lipase measurements are used in the diagnosis and treatment of diseases of the pancreas such as acute pancreatitis, obstruction of the pancreatic duct and pancreatic tumours.

METHOD PRINCIPLE

The colorimetric method is based on a lipase specific degradation of a chromogenic substrate. The specific lipase substrate-DGGMR [1,2-o-dilauryl-racglycero-3-glutaric acid-(6'-methylresorufin) ester] is cleaved by the catalytic action of lipase to form 1,2-o-dilauryl-racglycerol and an unstable intermediate, glutaric acid-(6-methyl resorufin) ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The lipase activity in the specimen is proportional to the production of methylresorufin in the reaction and can be determined photometrically.

REAGENTS

Package	Cat. No 4-185 (24-TRAY)	Cat. No 4-329 (36-TRAY)
1-Reagent	4 x 38 ml	4 x 23 ml
2-Reagent	4 x 20 ml	4 x 12.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 7 weeks, Biolis 24i Premium – 7 weeks.

Concentrations in the test

1-Reagent

TAPS [N-Tris(hydroxymethyl)methyl-3-aminopropanesulfonic acid]	100 mM
sodium hydroxide	40 mM
sodium deoxycholate preservative	34 mM

2-Reagent

tartaric acid	9.5 mM
sodium hydroxide	19 mM
colipase	460 IU/ml
2-propanol	0.65 M
DGGMR [1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)-ester]	0.4 mM

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, heparinized plasma free from hemolysis. Sample may be stored for up to 5 days at 2-8°C or 24 hours at 20-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LIPASE – TG MONO, CHOL – LIPASE, LDL DIRECT – LIPASE, TG – LIPASE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LIPASE – CALCIUM ARSENAZO, CHOL – LIPASE, HDL DIRECT – LIPASE, LDL DIRECT – LIPASE, TG – LIPASE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁴

Normal range	13 – 60 U/l	0,22 – 1,00 µkat/l
--------------	-------------	--------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 7 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 0.4 U/l (0.007 µkat/l).

- LoD (Limit of Detection):** 0.7 U/l (0.012 µkat/l).
- LoQ (Limit of Quantitation):** 5 U/l (0.083 µkat/l).
- Linearity :** up to 627 U/l (10.45 µkat/l).

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	58.7	0.58	1.00
level 2	103.4	0.88	0.85
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	54.2	2.13	3.9
level 2	97.1	3.25	3.3

Method comparison

A comparison between lipase values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 62 serum samples gave following results:

$$y = 0.9951x + 2.5134 \text{ U/l};$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

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Item name	26	LIPA				
Auto Rerun SW						
ON						
Auto Rerun Range (Result)						
	ON	ON				
	Lower	Higher				
Serum	7	600				
Urine						
Plasma						
CSF						
Dialysis						
Other						
Auto Rerun Condition (Absorbance)						
Absorbance Range	Lower	Higher	OFF			
	Higher		OFF			
Auto Rerun Condition (Prozone)						
Prozone Range	OFF					

■ Prestige 24i Premium, Biolis 24i Premium

Item No.	26	Item Name	LIPA	Optical	
Data information					
Units	U/l				
Decimals	1				
Analysis					
Type	RATE method				
Main Wave Length	570 nm				
Sub Wave Length	700 nm				
Method	Colorimetric				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Correlation					
Y=	Slope	Intercept			
	1	0			

Item No.	26	Item Name	LIPA	Optical	
Aspiration					
Kind	Double				
Vol.	Kind	Vol.	Add	Units	
	Sample	3	5	µl	
	Reagent 1	200	10	µl	
	Reagent 2	100	10	µl	
Data Process					
Read	Main	Start	End		
	Sub	35	42		
Abs.Limit					
Low	-3		High	3	
Correction value					
Blank correction					
End Point Limit					
Linear Check (%)					
2					
80					
Prozone Check					
	First	Start	End	Limit (%)	
	Second			Low	
Blank value					
Water Blank					
Reaction Monitor					
0 Level Point	1				
Span	3				
Third mixing					
OFF					

Item No.	26	Item Name	LIPA	Optical			
Normal Range							
	Male		Female				
	Low	High	Low	High			
Serum	13	60	13	60			
Urine							
Plasma	13	60	13	60			
CSF							
Dialysis							
Other							
Panic Range							
	Male		Female				
	Low	High	Low	High			
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							

Item No.	26	Item Name	LIPA	Optical			
Auto Rerun SW							
ON							
Auto Rerun Range (Conc.)							
	First Dil	Low	High				
		Re	Value	Dil	Re	Value	Dil
Serum			6			1000	
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Absorbance)							
	Lower	Higher	OFF				
	Higher		OFF				
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
100:Dil2							

PRESTIGE 24i LQ LIPASE

▪ Biolis 30i

Item no	26	Item name	LIPASE	Specimen	SERUM / PLASMA	OPTICAL
Data information			Aspiration volume			
UNITS	U/L		TYPE	Double		
DECIMALS	1					
Analysis				SAMPLE	REAGENT 1	REAGENT 2
METHOD	RATE method		VOL. (µL)	3	200	100
Main Wave Length	570 nm		BOTTLE (ml)			
Sub Wave Length	700 nm		FIRST DIL.			
CORRELATION (Y= AX + B)			Data processing read			
A =	1		START	END		
B =	0		MAIN	37		53
Blank value			SUB			
	• WATER ° REAGENT		ABS LIMIT			
Calibration				-3	TO	3
TYPE	Linear 1		Collection value			
STABILITY			END POINT	2.5		
			LINEARITY CHECK (%)	80		
			Prozone check			
				START	END	LIMIT (%)
			FIRST			
			SECOND			
				MINIMUM ABS.		
			°HIGH	MEAN		
			•LOW	VARIATE		

Item No	26	Item Name	LIPASE	Specimen	SERUM / PLASMA	OPTICAL
Reference intervals			Auto rerun			
MALE		FEMALE			•ON	°OFF
LOW	HIGH	LOW	HIGH			
13	60	13	60	Auto rerun range (conc.)		
Panic range			Re	Value	Dil.	
MALE		FEMALE			5	627
LOW	HIGH	LOW	HIGH	Auto rerun condition (abs.)		
						DIL.
Decision limit			LOWER	°ON	•OFF	
	MALE	FEMALE	HIGH	°ON	•OFF	
			Auto rerun condition (prozone)			
			°ON	•OFF		
			SAMPLE VOL.			
Reaction check			Dilution			
	°ON	•OFF	•DIL 1	° DIL 2		
CHECK						
LOW						
HIGH						
VL CHECK		VH CHECK				
°ON	•OFF	°ON	•OFF			