## PRESTIGE 24i LQ GLUCOSE

Nr kat. 4-201, 4-401 (PL)

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#### LITERATURA

1. Padhi, S., Nayak, A. K., Behera, A. Type II diabetes mellitus: a review on recent drug based therapeutics. Biomedicine & Pharmacotherapy, 131, 110708 (2020).
2. Pisarczyk-Wiza D., Zozulińska-Ziółkiewicz D. Glikokortykosteroidy a zaburzenia metabolizmu glukozy. Diabetologia Kliniczna 4,3, 110-116 (2015).
3. Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
4. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
5. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
6. Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
7. 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland
8. Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical Chemistry 48, 3, 436-472 (2002).
9. Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clinical Chemistry 50, 1704-1706 (2004).

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Cat. No 4-201, 4-401 (EN)

#### INTENDED USE


GLUCOSE reagent is intended to determine quantitatively glucose level in plasma, serum and cerebrospinal fluid. It is intended to diagnosis, monitoring and as an aid in the diagnosis of clinical conditions connected with abnormal glucose level. GLUCOSE reagent is intended to use on automatic analysers Biolis 24i Premium and Biolis 30i. It is only for *in vitro* diagnostics, for healthcare professional users.

#### SUMMARY <sup>1,2</sup>

Determinations of glucose level in serum, plasma and cerebrospinal fluid are mainly used in diagnosis and treatment monitoring of diabetes mellitus. Measurement of glucose concentration is also used as an aid in diagnosis of clinical conditions associated with hypoglycemia and hyperglycemia. Increased glucose levels can be observed in transient hyperglycemia, acute stress reaction, Cushing's syndrome, glucocorticoid therapy, hyperthyroidism. Decreased concentrations of glucose can be associated with neonatal hypoglycemia or hypothyroidism.

#### METHOD PRINCIPLE <sup>3</sup>

Colorimetric, enzymatic method with glucose oxidase.

The colour intensity is proportional to the glucose concentration.

#### REAGENTS

Package	Cat. No 4-201 (24-TRAY)	Cat. No 4-401 (36-TRAY)
1-REAGENT	6 x 60 ml	10 x 25 ml
<b>Number of tests:</b>		
Biolis 24i Premium	1320	900
Biolis 30i	1330	880

#### CONCENTRATION OF THE ACTIVE INGREDIENTS

##### IN THE REAGENT

##### 1-GLUCOSE / 1-REAGENT

phenol	5 mmol/l
glucose oxidase (GOD)	333.3 µkat/l
peroxidase (POD)	38.33 µkat/l
4- aminoantipyrene (4-AA)	0,75 mmol/l
phosphate buffer	
metal ion chelator	
stabilizer	
preservatives	

#### REAGENT STABILITY

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Biolis 24i Premium, Biolis 30i – 12 weeks.

#### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- EUH208 Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
- EUH210 Safety data sheet available on request.

#### SPECIMEN <sup>4,5,6</sup>

EDTA or heparinized plasma in tubes containing sodium fluoride or sodium iodoacetate additive/ serum, free from hemolysis, cerebrospinal fluid.

**Plasma / Serum.** Serum and plasma specimens should be separated from cells within 30 minutes after collection.

Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level.

Serum and plasma can be stored up to 8 hours at 25°C or 72 hours at 4 °C.

Plasma is the specimen recommended for the glucose determination in the blood.

**Cerebrospinal fluid.** Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes .

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

#### PROCEDURE

1-REAGENT is ready to use.

1-REAGENT put on basic position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- **Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE – CREA ENZYMATYC, GLUCOSE – PHOSPHORUS, GLUCOSE – TG, GLUCOSE – TG MONO, CK–MB – GLUCOSE. To avoid this effect follow the recommendations contained in the instruction:

- **Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE - PHOSPHORUS. To avoid this effect follow the recommendations contained in the instruction: 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER

Application guide for Prestige 24i analyser is available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

## QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionized water should be used to perform a blank test.

The calibration curve should be prepared every 12 weeks (Biolis 24i Premium) or every 10 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

## REFERENCE VALUES <sup>4,7,8,9</sup>

	mg/dl	mmol/l
plasma, serum <sup>7,8,9</sup>	70 – 99	3.9 – 5.5
cerebrospinal fluid <sup>4</sup>	40 – 70	2.2 – 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

## PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**  
0.0 mg/dl (0.0 mmol/l) – BS-400
- **LoD (Limit of Detection):**  
0.3 mg/dl (0.017 mmol/l) – BS-400

- **LoQ (Limit of Quantitation):**

6.0 mg/dl (0.33 mmol/l) – Biolis 30i

- **Linearity:**

up to 600 mg/dl (33.3 mmol/l) – Biolis 30i

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- **Measurement range:**

6.0 mg/dl (0.33 mmol/l) - 600 mg/dl (33.3 mmol/l)

- **Specificity / Interferences**

**a) samples containing low glucose concentration**

Haemoglobin up to 0.63 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

**b) samples containing high glucose concentration**

Haemoglobin up to 2.50 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

- **Precision (Biolis30i)**

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	86.6	0.45	0.5
level 2	293.6	2.38	0.8
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	87.8	2.23	2.5
level 2	291.6	5.32	1.8

- **Method comparison**

A comparison between glucose values determined at **Biolis 30i** (y) and at **SIEMENS ADVIA 1800** (x) using 59 plasma samples gave following results:

$$y = 1.0133 x + 0.1086 \text{ mg/dl;} \\ R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **SIEMENS ADVIA 1800** (x) using 56 serum samples gave following results:

$$y = 1.0559 x - 1.8604 \text{ mg/dl;} \\ R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 30 PMR samples gave following results:

$$y = 1.02 x + 0.4699 \text{ mg/dl;} \\ R = 0.992 \quad (R - \text{correlation coefficient})$$

## WASTE MANAGEMENT <sup>10</sup>

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03\*

## INCIDENTS <sup>11</sup>

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cornay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

## LITERATURE

1. Padhi, S., Nayak, A. K., Behera, A. Type II diabetes mellitus: a review on recent drug based therapeutics. Biomedicine & Pharmacotherapy, 131, 110708 (2020).
2. Pisarczyk-Wiza D., Zozulińska-Ziółkiewicz D. Glikokortykosteroidy a zaburzenia metabolizmu glukozy. Diabetologia Kliniczna 4,3, 110-116 (2015).
3. Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
4. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
5. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
6. Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
7. 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland
8. Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical Chemistry 48, 3, 436-472 (2002).
9. Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clinical Chemistry 50, 1704-1706 (2004).
10. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
11. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

## LIST OF CHANGES

Previous version: 03	Current version: 04
Sections updated: <i>WARNING AND NOTES</i>	

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- Padhi, S., Nayak, A. K., Behera, A. Type II diabetes mellitus: a review on recent drug based therapeutics. Biomedicine & Pharmacotherapy, 131, 110708 (2020).
- Pisarczyk-Wiza D., Zozulińska-Ziółkiewicz D. Glikokortykosteroidy a zaburzenia metabolizmu glukozy. Diabetologia Kliniczna 4,3, 110-116 (2015).
- Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
- 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical Chemistry 48, 3, 436-472 (2002).
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clinical Chemistry 50, 1704-1706 (2004).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Item No.	7	Item Name	GLUC	Optical	
<b>Aspiration</b>					
Kind	Single				
Vol.					
	Kind	Vol.	Add	Units	
	Sample	3	5	µl	
	Reagent 1	250	10	µl	
	Reagent 2				
<b>Blank value</b>					
Water Blank					
<b>Reaction Monitor</b>					
0 Level Point	1				
Span	3				
<b>Third mixing</b>					
OFF					
<b>Data Process</b>					
Read	Start	End			
	Main	46	47		
	Sub	9	10		
Abs.Limit	Low	-0.1	~	High	2
<b>Correction value</b>					
Blank correction					
End Point Limit					
Linear Check (%)					
<b>Prozone Check</b>					
	Start	End	Limit (%)		
	First				
	Second		Low		

Item No.	7	Item Name	GLUC	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	74	106	74	106
Urine				
Plasma	74	106	74	106
CSF	40	70	40	70
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	7	Item Name	GLUC	Optical			
<b>Auto Rerun SW</b>							
ON							
<b>Auto Rerun Range (Conc.)</b>							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum			6			500	
Urine							
Plasma							
CSF							
Dialysis							
Other							
<b>Auto Rerun Condition (Absorbance)</b>							
	Lower	OFF					
	Higher	OFF					
<b>Auto Rerun Condition (Prozone)</b>							
OFF							
<b>Dilution</b>							
100:DiI2							

• **Biolis 30i**

Item no	7	Item name	GLUC	Specimen	SERUM/ PLASMA/CSF	OPTICAL												
<b>Data information</b>				<b>Aspiration volume</b>														
UNITS		mg/dL		TYPE														
DECIMALS		1		Single														
<b>Analysis</b>				<table border="1"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>3</td> <td>250</td> <td></td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>				SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	3	250		BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	3	250																
BOTTLE (ml)																		
METHOD				END method														
Main Wave Length		505 nm		FIRST DIL.														
Sub Wave Length		700 nm																
<b>CORRELATION (Y= AX + B)</b>				<b>Data processing read</b>														
A =		1		START		END												
B =		0		MAIN		46												
				SUB		47												
<b>Blank value</b>				ABS LIMIT														
• WATER		° REAGENT		-0.1 TO 2														
<b>Calibration</b>				<b>Collection value</b>														
TYPE		Linear 2		END POINT		2.5												
STABILITY				LINEARITY CHECK (%)		0												
<b>Prozone check</b>				<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> <td>LIMIT (%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> </table>				START	END	LIMIT (%)	FIRST				SECOND			
	START	END	LIMIT (%)															
FIRST																		
SECOND																		
°HIGH				MINIMUM ABS.														
•LOW				MEAN														
				VARIATE														

Item No	7	Item Name	GLUC	Specimen	SERUM/PLASMA/CSF	OPTICAL
<b>Reference intervals</b>						
MALE			FEMALE			
LOW	HIGH	LOW	HIGH			
70	99	70	99			
<b>Panic range</b>						
MALE			FEMALE			
LOW	HIGH	LOW	HIGH			
<b>Decision limit</b>						
MALE			FEMALE			
<b>Reaction check</b>						
°ON		•OFF				
CHECK						
LOW						
HIGH						
<b>VL CHECK</b>			<b>VH CHECK</b>			
°ON	•OFF	°ON	•OFF			

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