## PRESTIGE 24i LQ CK-MB

Nr kat. 4-227, 4-427 (PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of CK-MB fraction activity used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased CK-MB serum level is a strong marker of myocardial infarction.

### METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC) with use of antibodies against CK-M fraction. Specific antibodies against CK-M inhibit the complete CK-MM activity (which is the main part of total CK activity) and the CK-M subunit of CK-MB. Only CK-B activity is measured.

The rate of absorbance changes at  $\lambda=340$  nm is directly proportional to half of CK-MB activity (B subunit activity).

### REAGENTS

#### Package

	Cat. No 4-227 (24-TRAY)	Cat. No 4-427 (36-TRAY)
1-Reagent	2 x 40 ml	2 x 22.5 ml
2-Reagent	2 x 10.5 ml	2 x 6 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks.

### Concentrations in the test

1-Reagent		
imidazole buffer		100 mmol/l
glucose		20 mmol/l
N-acetylcysteine		20 mmol/l
magnesium acetate		10 mmol/l
EDTA		2 mmol/l
NADP		2 mmol/l
ADP		2 mmol/l
AMP		5 mmol/l
HK		> 2.5 U/ml
polyclonal antibodies against CK-M; inhibiting capacity		8000 U/l

### 2-Reagent

diadenosinepentaphosphate	10 μmol/l
glucose-6-phosphate-dehydrogenase (G6P-DH)	> 1.5 U/ml
creatine phosphate preservatives	30 mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- Do not use reagents past the expiry date.
- Do not interchange caps among reagents.
- Results CK-MB can be falsely high in case of prostate, kidney, ovary, breast and bladder cancer when isoenzyme CK-BB appears in the blood.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-Reagent contains imidazole.

### Danger

- H360 May damage fertility or the unborn child.
- P201 Obtain special instructions before use.
- P202 Do not handle until all safety precautions have been read and understood.
- P308+P313 IF exposed or concerned: Get medical advice/attention.
- P405 Store locked up.
- P501 Dispose of the contents/containers in accordance with the current legislation on waste treatment.

### SPECIMEN

Serum, free from hemolysis. CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light. Specimens can be stored up to 4-8 hours at 15-25°C or 1-2 days at 2-8°C or 1 month at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

- 1-Reagent and 2-Reagent are ready to use.
- 1-Reagent put on basic position in reagent tray.
- 2-Reagent put on start position in reagent tray.
- For reagent blank deionized water is recommended.

### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CK – MB – CREATININE, CK – MB – GLUCOSE, CREATININE – CK - MB. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CK-MB – CREATININE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

### REFERENCE VALUES <sup>9</sup>

serum	37°C	
adults	up to 24 U/l	up to 0.401 μkat/l

The probability that cardiac infarction has occurred is high when CK-MB and total CK activities are above normal values and CK-MB activity is between 6 and 25% of the total CK activity. It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY CK-MB CONTROL N (Cat. No 5-183) and CORMAY CK-MB CONTROL P (Cat. No 5-184).

For the calibration the CORMAY CK-MB CALIBRATOR (Cat. No 5-182) is recommended.

The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 0.69 U/l (0.012 μkat/l)
- LoD (Limit of Detection):** 1.33 U/l (0.022 μkat/l)
- LoQ (Limit of Quantitation):** 5 U/l (0.083 μkat/l)
- Linearity:** up to 2800 U/l (46,67 μkat/l).

For higher activity, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

### Specificity / Interferences

Haemoglobin interfere even in small amounts, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	38.55	0.35	0.90
level 2	132.63	0.56	0.42
Reproducibility (day to day) n = 80	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	38.4	1.90	5.0
level 2	135.2	4.45	3.3

### Method comparison

A comparison between CK-MB values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 63 serum samples gave following results:

$$y = 0.9953x + 0.5669 \text{ U/l};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	20	CKMB
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	10	2150
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
<b>Prozone Range</b>		
		OFF

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	20	Item Name	CKMB	Optical
<b>Data information</b>				
Units	U/L			
Decimals	2			
<b>Calibration</b>				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
<b>Analysis</b>				
Type	RATE method			
Main Wave Length	340nm			
Sub Wave Length	405nm			
Method	IFCC			
<b>Correlation</b>				
Y=	Slope	Intercept		
	1	0		

Item No.	20	Item Name	CKMB	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		10	5	µl
Reagent 1		200	10	µl
Reagent 2		40	10	µl
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	38	52	
Abs.Limit Low High				
	-3		3	
<b>Correction value</b>				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	90			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				

Item No.	20	Item Name	CKMB	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	24	0	24
Urine				
Plasma	0	24	0	24
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	20	Item Name	CKMB	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low Value	Dil	High Value
Serum		6		2100
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	OFF		
	Higher	OFF		
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:DiI2				

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