## PRESTIGE 24i LQ PHOSPHORUS

Nr kat. 4-243, 4-443 (PL)

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(EN)

#### INTENDED USE

PHOSPHORUS reagent is intended to determine quantitatively phosphorus levels in plasma, serum and urine. It is intended to monitor and as an aid to the diagnosis of clinical conditions associated with abnormal phosphorus levels. PHOSPHORUS reagent is intended to use on automatic analysers Biolis 24i Premium and Biolis 30i. It is only for *in vitro* diagnostics, for healthcare professional users.

#### SUMMARY<sup>1-4</sup>

Determinations of phosphorus levels in plasma, serum and urine are used to monitor and as an aid to the diagnosis of parathyroid and kidney diseases or vitamin D imbalance. Increased phosphorus levels may be associated with increased phosphate intake, cell lysis (leukemia), transcellular phosphate shift (acidosis). Decreased phosphorus levels may occur in refeeding syndrome, decreased net intestinal phosphate absorption, lowered phosphate threshold in the kidneys.

#### METHOD PRINCIPLE<sup>5,6</sup>

Direct phosphomolybdate reaction without deproteinization. Phosphate ions form with molybdate ions in acid solution proportional amounts of unreduced phosphomolybdate complex. The concentration of the complex formed is determined by measuring its absorbance.

#### REAGENTS

##### Package

	Cat. No 4-243 (24-TRAY)	Cat. No 4-443 (36-TRAY)
1-REAGENT	6 x 59 ml	12 x 27.5 ml

##### Number of tests:

Biolis 24i Premium	1090	1000
Biolis 30i	1090	990

#### CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

##### 1-REAGENT

ammonium molybdate	0,4 mmol/l
sulphuric acid	
hydrochloric acid	
detergents	

#### REAGENT STABILITY

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Biolis 24i Premium, Biolis 30i – 12 weeks.

#### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Ingredients:

1-REAGENT contains sulfuric acid (VI) and hydrochloric acid.

#### Danger

H314 Causes severe skin burns and eye damage.  
 P280 Wear protective gloves, protective clothing, eye protection or face protection.  
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P310 Immediately call a POISON CENTER or doctor.

#### SPECIMEN<sup>7,8,9</sup>

Serum, heparinized plasma free from hemolysis, 24-hours urine.

Serum is the preferred specimen. Level of inorganic phosphate in heparinized plasma is about 0.2 - 0.3 mg/dl (0.06 – 0.10 mmol/l) lower than in serum.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain several times higher phosphate concentration than normal serum.

Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

**Urine preparation:** For collection of 24-hour urine specimens, 20-30 ml of 6M HCl should be added into the container to avoid phosphate precipitations. Dilute 1 part of acidified urine with 9 parts of distilled water before determination. Multiply the result by the dilution factor.

24-hours urine samples can be stored up to 2 days at 20-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

#### ASSAY PROCEDURE

1-REAGENT is ready to use.

1-REAGENT put on basic position in reagent tray.

For reagent blank deionized water is recommended.

### Actions required:

- **Biolis 24i Premium:** When performing assays in analyser, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE – PHOSPHORUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- **Biolis 30i:** When performing assays in analyser, there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO - PHOSPHORUS, GLUCOSE - PHOSPHORUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

Application guide for Prestige 24i analyser is available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionized water should be used to perform a blank test.

The calibration curve should be prepared every 12 weeks (Biolis 24i Premium, Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

### REFERENCE VALUES <sup>7</sup>

serum / plasma	mg/dl	mmol/l
age:		
0 – 9 d	4.5 – 9.0	1.45 – 2.91
10 d - 2 y	4.0 – 6.5	1.29 – 2.10
3 – 9 y	3.2 – 5.8	1.03 – 1.87
10 – 15 y	3.3 – 5.4	1.07 – 1.74
16 – 59 y	2.4 – 4.4	0.78 – 1.42
60 - 89 y males	2.3 – 3.7	0.74 – 1.20
60 - 89 y females	2.8 – 4.0	0.90 – 1.26
>89 y females	2.5 - 4.2	0.81 - 1.36
>89 y males	2.2 - 3.9	0.71 - 1.26

24-hours urine	g/24h	mmol/24h
	0.4 – 1.3	12.9 – 42.0

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result

### Phosphorus concentration in 24-hours urine – calculation

phosphorus concentration in 24-hours urine [g/24h]	=	phosphorus concentration in sample of 24-hours urine [mg/dl]	x	urine volume of 24-hours urine [dl/24h]	÷	1000
--	---	--	---	---	---	------

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**  
0.04 mg/dl (0.013 mmol/l)
- **LoD (Limit of Detection):**  
0.07 mg/dl (0.023 mmol/l)
- **LoQ (Limit of Quantitation):**  
0.13 mg/dl (0.042 mmol/l)
- **Linearity:**  
up to 20 mg/dl (6.46 mmol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- **Measurement range**  
0.13 mg/dl (0.042 mmol/l) - 20 mg/dl (6.46 mmol/l)

### Specificity / Interferences

Haemoglobin up to 0.12 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	3.36	0.08	2.4
level 2	7.31	0.11	1.6
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	3.31	0.07	2.0
level 2	7.27	0.09	1.3

### Method comparison

A comparison between phosphorus values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 72 serum samples gave following results:

$$y = 0.989x + 0.0358 \text{ mg/dl};$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 50 urine samples gave following results:

$$y = 0.9804x - 0.7055 \text{ mg/dl};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT <sup>10</sup>

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents: 18 01 06\*
- Empty packages: 15 01 10\*
- Wastewater from the analyzer: 18 01 03\*

### INCIDENTS <sup>11</sup>

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

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### LIST OF CHANGES

Previous version: 03	Current version: 04
Sections updated: <i>WARNING AND NOTES</i>	

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Item No.	21	Item Name	PHOS	Optical
<b>Aspiration</b>				
Kind	Single			
Vol.	Kind	Vol.	Add	Units
	Sample	3	5	µl
	Reagent 1	300	10	µl
<b>Blank value</b>				
<b>Water Blank</b>				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
<b>OFF</b>				
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	9	10	
Abs.Limit	Low	-0.1	High	2
<b>Correction value</b>				
Blank correction				
End Point Limit				
Linear Check (%)				
<b>Prozone Check</b>				
	First	Start	End	Limit (%)
	Second			Low

Item No.	21	Item Name	PHOS	Optical
<b>Normal Range *</b>				
	Male		Female	
	Low	High	Low	High
Serum	2.4	4.4	2.4	4.4
Urine				
Plasma	2.4	4.4	2.4	4.4
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	21	Item Name	PHOS	Optical
<b>Auto Rerun SW</b>				
<b>ON</b>				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low Value	Dil	High Value
	Re	0.13		20
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	OFF		
	Higher	OFF		
<b>Auto Rerun Condition (Prozone)</b>				
<b>OFF</b>				
<b>Dilution</b>				
<b>100:Dil2</b>				

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Item No	21	Item Name	PHOS	Specimen	SERUM / PLASMA / URINE	OPTICAL
<b>Reference intervals *</b>						
MALE			FEMALE			
LOW	HIGH	LOW	HIGH			
2.4	4.4	2.4	4.4			
<b>Panic range</b>						
MALE			FEMALE			
LOW	HIGH	LOW	HIGH			
<b>Decision limit</b>						
MALE			FEMALE			
<b>Reaction check</b>						
°ON			°OFF			
CHECK						
LOW						
HIGH						
<b>VL CHECK</b>			<b>VH CHECK</b>			
°ON			°OFF			
°OFF			°ON			
			°OFF			

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