PRESTIGE 24i LQ TG

Nr kat. 4-253, 4-453 (PL)

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Cat. No **4-253, 4-453** (EN)

INTENDED USE

Diagnostic kit for determination of triglycerides concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.





The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

The colour intensity is proportional to the triglycerides concentration.

REAGENTS

Package

	Cat. No 4-253 (24-TRAY)	Cat. No 4-453 (36-TRAY)
1-Reagent	6 x 40 ml	8 x 23 ml
2-Reagent	6 x 12.5 ml	8 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 10 weeks.

Concentrations in the test

buffer PIPES (pH 7.0)	40 mmol/l
4-aminoantipyrine (4-AA)	0.4 mmol/l
ATP	1.5 mmol/l
Mg ²⁺	1.6 mmol/l
ADPS	0.6 mmol/l
glycerol kinase (GK)	> 66.67 µkat/l
glycerol-3-phosphate oxidase (GPO)	> 60.00 µkat/l
peroxidase (POD)	> 20.00 µkat/l
lipoprotein lipase (LPL)	> 16.67 µkat/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis.

Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 3 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: TG – LIPASE, TG – LIPASE II GEN, ALAT – TG, ALAT (II GEN) – TG, ALP – TG, GLUCOSE – TG, HDL DIRECT – TG, LDL DIRECT – TG, TOTAL PROTEIN – TG. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: TG – LIPASE, ALAT (II GEN) – TG, LDL DIRECT – TG, TOTAL PROTEIN – TG. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁷

serum, plasma	< 150 mg/dl < 1.7 mmol/l
---------------	-----------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control, it is recommended to use, with each batch of samples the CORMAY SERUM HN (Cat. No 5-172), CORMAY SERUM HP (Cat. No 5-173) and CORMAY LIPID CONTROL LEVEL 1 (Cat. No 5-179), CORMAY LIPID CONTROL LEVEL 2 (Cat. No 5-180).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 6 weeks (Prestige 24i) or every 10 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.5 mg/dl (0.006 mmol/l)
- LoD (Limit of Detection):**
1.1 mg/dl (0.012 mmol/l)
- LoQ (Limit of Quantitation):**
3.5 mg/dl (0.04 mmol/l)
- Linearity:**
up to 2000 mg/dl (22.6 mmol/l)

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.50 g/dl, bilirubin up to 20 mg/dl and ascorbate up to 62 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	84.3	1.17	1.39
level 2	188.0	2.66	1.42
Reproducibility (day to day) n = 80	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	177.8	5.96	3.4
level 2	90.1	3.56	4.0

Method comparison

A comparison between triglycerides values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 69 serum samples gave following results:

$$y = 1.017x - 1.4956 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	11	TRIG			
Auto Rerun SW					
ON					
Auto Rerun Range (Result)					
	ON	ON			
	Lower	Higher			
Serum	6	1500			
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Absorbance)					
Absorbance Range		Lower	OFF		
		Higher	OFF		
Prozone Range		OFF			

• Prestige 24i Premium, Biolis 24i Premium

Item No.	11	Item Name	TRIG	Optical		
Data information						
Units	mg/dl					
Decimals	1					
Calibration						
Type	Linear2					
Std sample conc.						
Blank	0	#1	*	#2	*	
#3		#4		#5		
#6						
Analysis						
Type	END method					
Main Wave Length	505nm					
Sub Wave Length	700nm					
Method	GPO-POD					
Correlation						
Y=	Slope		Intercept			
	1		0			

Item No.	11	Item Name	TRIG	Optical		
Aspiration						
Kind	Double					
Data Process						
Read		Start	End			
	Main	51	52			
	Sub	30	31			
Abs.Limit						
Low	-0.1		High	2		
Blank value						
Water Blank						
Reaction Monitor						
0 Level Point	1					
Span	3					
Third mixing						
OFF						
Correction value						
Blank correction						
End Point Limit						
2						
Linear Check (%)						
Prozone Check						
	Start	End	Limit (%)			
First						
Second			Low			

Item No.	11	Item Name	TRIG	Optical		
Normal Range						
	Male		Female			
	Low	High	Low	High		
Serum	0	150	0	150		
Urine						
Plasma						
CSF						
Dialysis						
Other						
Panic Range						
	Male		Female			
	Low	High	Low	High		
Serum						
Urine						
Plasma						
CSF						
Dialysis						
Other						

Item No.	11	Item Name	TRIG	Optical		
Auto Rerun SW						
ON						
Auto Rerun Range (Conc.)						
	First Dil	Low	High			
		Value	Dil	Value	Dil	
Serum		11.5		2000		
Urine						
Plasma						
CSF						
Dialysis						
Other						
Auto Rerun Condition (Absorbance)						
Lower		OFF				
Higher		OFF				
Auto Rerun Condition (Prozone)						
OFF						
Dilution						
100:Dil2						

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