## PRESTIGE 24i LQ AMYLASE

Nr kat. **4-255, 4-455** (PL)

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Cat. No **4-255, 4-455** (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -amylase activity intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\alpha$ -Amylase is a digestive enzyme secreted by salivary glands and pancreas. Low level of amylase is also found in skeletal muscle, adipose tissue and fallopian tubes.  $\alpha$ - Amylase is measured generally in pancreas diseases. Elevation of amylase activity is observed also due to inflammation of abdominal cavity or salivary glands.

### METHOD PRINCIPLE

2-Chloro-4-nitrophenyl- $\alpha$ -maltotriose (CNP-G3) is a direct substrate for determination of  $\alpha$ -amylase activity, which does not require the presence of ancillary enzymes.

10 CNP-G3  $\xrightarrow{\alpha\text{-amylase}}$  9 CNP + CNP-G2 + 9 maltotriose + glucose

The rate of 2-chloro-4-nitrophenol formation can be monitored at 405 nm and is proportional to the  $\alpha$ -amylase activity.

### REAGENTS

#### Package

|           |                      |                      |
|-----------|----------------------|----------------------|
|           | <b>Cat. No 4-255</b> | <b>Cat. No 4-455</b> |
|           | <b>(24-TRAY)</b>     | <b>(36-TRAY)</b>     |
| 1-REAGENT | 4 x 60 ml            | 10 x 25 ml           |

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. Stability on board of the analyzer at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 12 weeks, Biolis 30i - 12 weeks.

### Concentrations in the test

|   |               |
|---|---------------|
| MES buffer                                    | < 120 mmol/l  |
| calcium acetate                               | < 7 mmol/l    |
| potassium hydroxide                           | < 40 mmol/l   |
| potassium thiocyanate                         | < 1100 mmol/l |
| 2-chloro-4-nitrophenyl- $\alpha$ -maltotriose | < 2 mmol/l    |
| preservative, stabilizer                      |               |

### Warnings and notes

- Protect from direct sunlight!
- Prevent the reagent from microbiological contamination and from saliva and sweat  $\alpha$ -amylase! Saliva and sweat contain  $\alpha$ -amylase. Do not pipette by the mouth, avoid skin contact with reagent, specimens, tips, cuvettes. Ensure to use automatic pipettes and laboratory gloves.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum or plasma collected on heparin, free from hemolysis, urine.

Do not use anticoagulants: EDTA, citrates and oxalates as they inhibit amylase activity.

Serum / plasma can be stored for 7 days at 15-25°C or for one month at 2-8°C.<sup>7</sup>

Urine can be stored for 2 days at 15-25°C or for 10 days at 2-8°C.<sup>9</sup> Amylase is very unstable in acid urine. Adjust pH to approximately 7.0 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT is ready to use.

1-REAGENT put on basic position in reagent tray.

For reagent blank deionized water is recommended.

### REFERENCE VALUES <sup>8</sup>

| serum / plasma | U/l      | $\mu$ kat/l  |
|----------------|----------|--------------|
|                | 20 – 104 | 0.34 – 1.77  |
| urine          | U/l      | $\mu$ kat/l  |
|                | 32 – 641 | 0.54 – 10.90 |

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analyzers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i, Biolis 24i Premium) or 12 weeks (Biolis 30i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzer Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
1.1 U/l (0.018  $\mu$ kat/l)
- LoD (Limit of Detection):**  
1.7 U/l (0.028  $\mu$ kat/l)

- LoQ (Limit of Quantitation):**  
4.5 U/l (0.075  $\mu$ kat/l)

- Linearity:**  
up to 2300 U/l (38.3  $\mu$ kat/l)

### Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

|   | Mean [U/l] | SD [U/l] | CV [%] |
|---|------------|----------|--------|
| <b>Repeatability</b><br>(run to run) n = 20   |            |          |        |
| level 1                                       | 70.4       | 0.97     | 1.38   |
| level 2                                       | 433.7      | 3.96     | 0.91   |
| <b>Reproducibility</b><br>(day to day) n = 80 |            |          |        |
| level 1                                       | 69.2       | 2.88     | 4.2    |
| level 2                                       | 440.0      | 14.70    | 3.3    |

### Method comparison

A comparison between amylase values determined at **Biolis 30i** (y) and at **BS-800** (x) using 60 serum samples gave following results:

$y = 1.0142x - 5.8464$  U/l;

$R = 0.999$  (R – correlation coefficient)

A comparison between amylase values determined at **Biolis 30i** (y) and at **BS-800** (x) using 62 urine samples gave following results:

$y = 1.0214x - 7.6945$  U/l;

$R = 1.000$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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|                      |          |                  |                  |
|----------------------|----------|------------------|------------------|
| Item name            | 14       | AMYL             |                  |
| <b>Aspiration</b>    |          |                  |                  |
| Kind                 | Single   |                  |                  |
| Vol.                 |          |                  |                  |
| Kind                 | Vol.     | Add              | Units            |
| Sample               | 5        | 5                | µl               |
| Reagent1             | 250      | 10               | µl               |
| Reagent2             |          |                  |                  |
| <b>Data Process</b>  |          |                  |                  |
| Read                 | Start    | End              | Absorbance Limit |
| Main                 | 20       | 30               | Low -0.100       |
| Sub                  |          |                  | High 2.000       |
| <b>Factor</b>        |          |                  |                  |
| Blank correction     |          | Endpoint Limit   | 2.100            |
|                      |          | Linear Check (%) | 90               |
| <b>Dilution</b>      |          |                  |                  |
| Diluent              | 100:Di12 |                  |                  |
| <b>Prozone Check</b> |          |                  |                  |
|                      | Start    | End              | Limit (%)        |
| First                |          |                  |                  |
| Second               |          |                  | Low              |
| Third                |          |                  | Low              |
| <b>Monitor</b>       |          |                  |                  |
| 0 Level Point        | 1        |                  |                  |
| Span                 | 3.000    |                  |                  |
| <b>Third mixing</b>  |          |                  |                  |
| OFF                  |          |                  |                  |

|  |        |        |
|--|--------|--------|
| Item name                                | 14     | AMYL   |
| <b>Auto Rerun SW</b>                     |        |        |
| ON                                       |        |        |
| <b>Auto Rerun Range (Result)</b>         |        |        |
|  | ON     | ON     |
|  | Lower  | Higher |
| Serum                                    | 5.3    | 1500   |
| Urine                                    |        |        |
| Plasma                                   |        |        |
| CSF                                      |        |        |
| Dialysis                                 |        |        |
| Other                                    |        |        |
| <b>Auto Rerun Condition (Absorbance)</b> |        |        |
| Absorbance Range                         |        |        |
|  | Lower  | OFF    |
|  | Higher | OFF    |
| <b>Prozone Range</b>                     |        |        |
| OFF                                      |        |        |

### ■ Prestige 24i Premium, Biolis 24i Premium

|                         |             |           |      |         |   |
|-------------------------|-------------|-----------|------|---------|---|
| Item No.                | 14          | Item Name | AMYL | Optical |   |
| <b>Data information</b> |             |           |      |         |   |
| Units                   | U/L         |           |      |         |   |
| Decimals                | 1           |           |      |         |   |
| <b>Analysis</b>         |             |           |      |         |   |
| Type                    | RATE method |           |      |         |   |
| Main Wave Length        | 405nm       |           |      |         |   |
| Sub Wave Length         | 660nm       |           |      |         |   |
| Method                  | CNP-G3      |           |      |         |   |
| <b>Calibration</b>      |             |           |      |         |   |
| Type                    | Linear2     |           |      |         |   |
| Std sample conc.        |             |           |      |         |   |
| Blank                   | 0           | #1        | *    | #2      | * |
| #3                      |             | #4        |      | #5      |   |
| #6                      |             |           |      |         |   |
| <b>Correlation</b>      |             |           |      |         |   |
| Y=                      | Slope       | Intercept |      |         |   |
|                         | 1           | 0         |      |         |   |

|                         |        |           |           |         |
|-------------------------|--------|-----------|-----------|---------|
| Item No.                | 14     | Item Name | AMYL      | Optical |
| <b>Aspiration</b>       |        |           |           |         |
| Kind                    | Single |           |           |         |
| Vol.                    |        |           |           |         |
| Kind                    | Vol.   | Add       | Units     |         |
| Sample                  | 3      | 5         | µl        |         |
| Reagent 1               | 200    | 10        | µl        |         |
| <b>Data Process</b>     |        |           |           |         |
| Read                    | Main   | Start     | End       |         |
|                         | Sub    | 20        | 30        |         |
| Abs.Limit Low High      |        |           |           |         |
|                         | -0.1   | ~         | 2         |         |
| <b>Blank value</b>      |        |           |           |         |
| Water Blank             |        |           |           |         |
| <b>Reaction Monitor</b> |        |           |           |         |
| 0 Level Point           | 1      |           |           |         |
| Span                    | 3      |           |           |         |
| <b>Correction value</b> |        |           |           |         |
| Blank correction        |        |           |           |         |
| End Point Limit 2.1     |        |           |           |         |
| Linear Check (%) 10     |        |           |           |         |
| <b>Prozone Check</b>    |        |           |           |         |
|                         | Start  | End       | Limit (%) |         |
| First                   |        |           |           |         |
| Second                  |        |           |           | Low     |

|                     |      |           |        |         |
|---------------------|------|-----------|--------|---------|
| Item No.            | 14   | Item Name | AMYL   | Optical |
| <b>Normal Range</b> |      |           |        |         |
|                     | Male |           | Female |         |
|                     | Low  | High      | Low    | High    |
| Serum               | 20   | 104       | 20     | 104     |
| Urine               | 32   | 641       | 32     | 641     |
| Plasma              | 20   | 104       | 20     | 104     |
| CSF                 |      |           |        |         |
| Dialysis            |      |           |        |         |
| Other               |      |           |        |         |
| <b>Panic Range</b>  |      |           |        |         |
|                     | Male |           | Female |         |
|                     | Low  | High      | Low    | High    |
| Serum               |      |           |        |         |
| Urine               |      |           |        |         |
| Plasma              |      |           |        |         |
| CSF                 |      |           |        |         |
| Dialysis            |      |           |        |         |
| Other               |      |           |        |         |

|  |           |           |      |            |     |
|--|-----------|-----------|------|------------|-----|
| Item No.                                 | 14        | Item Name | AMYL | Optical    |     |
| <b>Auto Rerun SW</b>                     |           |           |      |            |     |
| ON                                       |           |           |      |            |     |
| <b>Auto Rerun Range (Conc.)</b>          |           |           |      |            |     |
|  | First Dil | Low Value | Dil  | High Value | Dil |
| Serum                                    |           | 2.5       |      | 1500       |     |
| Urine                                    |           |           |      |            |     |
| Plasma                                   |           |           |      |            |     |
| CSF                                      |           |           |      |            |     |
| Dialysis                                 |           |           |      |            |     |
| Other                                    |           |           |      |            |     |
| <b>Auto Rerun Condition (Absorbance)</b> |           |           |      |            |     |
| ON                                       |           |           |      |            |     |
|  | Lower     | OFF       |      |            |     |
|  | Higher    | OFF       |      |            |     |
| <b>Auto Rerun Condition (Prozone)</b>    |           |           |      |            |     |
| OFF                                      |           |           |      |            |     |
| <b>Dilution</b>                          |           |           |      |            |     |
| 100:Di12                                 |           |           |      |            |     |

## PRESTIGE 24i LQ AMYLASE

**Biolis 30i**

|                                |             |           |                          |           |                    |         |
|--------------------------------|-------------|-----------|--------------------------|-----------|--------------------|---------|
| Item no                        | 14          | Item name | AMYL                     | Specimen  | SERUM/PLASMA/URINE | OPTICAL |
| <b>Data information</b>        |             |           |                          |           |                    |         |
| UNITS                          | U/L         |           | <b>Aspiration volume</b> |           |                    |         |
| DECIMALS                       | 1           |           | TYPE                     |           |                    |         |
| <b>Analysis</b>                |             |           | Single                   |           |                    |         |
| METHOD                         | RATE method |           | SAMPLE                   | REAGENT 1 | REAGENT 2          |         |
| Main Wave Length               | 405 nm      |           | VOL. (µL)                | 3         | 200                |         |
| Sub Wave Length                | 660 nm      |           | BOTTLE (ml)              |           |                    |         |
| <b>CORRELATION (Y= AX + B)</b> |             |           | FIRST DIL.               |           |                    |         |
| A =                            | 1           |           | Data processing read     |           |                    |         |
| B =                            | 0           |           | START                    |           | END                |         |
| <b>Blank value</b>             |             |           | MAIN                     |           |                    |         |
| ° WATER ° REAGENT              |             |           | SUB                      |           |                    |         |
| <b>Calibration</b>             |             |           | ABS LIMIT                |           |                    |         |
| TYPE                           | Linear 2    |           | -0.1                     |           | TO 2               |         |
| STABILITY                      |             |           | <b>Collection value</b>  |           |                    |         |
|                                |             |           | END POINT                |           | 2.5                |         |
|                                |             |           | LINEARITY CHECK (%)      |           | 90                 |         |
|                                |             |           | <b>Prozone check</b>     |           |                    |         |
|                                |             |           | START                    | END       | LIMIT (%)          |         |
|                                |             |           | FIRST                    |           |                    |         |
|                                |             |           | SECOND                   |           |                    |         |
|                                |             |           | °HIGH                    |           | MINIMUM ABS.       |         |
|                                |             |           | °LOW                     |           | MEAN               |         |
|                                |             |           |                          |           | VARIATE            |         |

|                            |      |                 |      |                                       |                    |         |
|----------------------------|------|-----------------|------|---------------------------------------|--------------------|---------|
| Item No                    | 14   | Item Name       | AMYL | Specimen                              | SERUM/PLASMA/URINE | OPTICAL |
| <b>Reference intervals</b> |      |                 |      |                                       |                    |         |
| MALE                       |      | FEMALE          |      | <b>Auto rerun</b>                     |                    |         |
| LOW                        | HIGH | LOW             | HIGH | °ON                                   |                    | °OFF    |
| 20                         | 104  | 20              | 104  | <b>Auto rerun range (conc.)</b>       |                    |         |
| <b>Panic range</b>         |      |                 |      | Re                                    | Value              | Dil.    |
| LOW                        | HIGH | LOW             | HIGH | 4.5                                   | 2300               |         |
|                            |      |                 |      | <b>Auto rerun condition (abs.)</b>    |                    |         |
|                            |      |                 |      |                                       |                    | DIL.    |
|                            |      |                 |      | LOWER                                 | °ON                | °OFF    |
|                            |      |                 |      | HIGH                                  | °ON                | °OFF    |
|                            |      |                 |      | <b>Auto rerun condition (prozone)</b> |                    |         |
|                            |      |                 |      | °ON                                   |                    | °OFF    |
|                            |      |                 |      | SAMPLE VOL.                           |                    |         |
|                            |      |                 |      | <b>Dilution</b>                       |                    |         |
|                            |      |                 |      | °DIL 1                                |                    | °DIL 2  |
| <b>Reaction check</b>      |      |                 |      |                                       |                    |         |
|                            |      |                 |      | °ON                                   |                    | °OFF    |
|                            |      |                 |      | CHECK                                 |                    |         |
|                            |      |                 |      | LOW                                   |                    |         |
|                            |      |                 |      | HIGH                                  |                    |         |
| <b>VL CHECK</b>            |      | <b>VH CHECK</b> |      |                                       |                    |         |
| °ON                        | °OFF | °ON             | °OFF |                                       |                    |         |

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