## PRESTIGE 24i CHOLINESTERASE

Nr kat. **4-256, 4-456** (PL)

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Cat. No **4-256, 4-456** (EN)

### INTENDED USE

Diagnostic kit for determination of cholinesterase activity intended to use in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

There are two cholinesterases (CHE and ACHE) differing in substrate specificity, tissue of origin and biological role. The term cholinesterase (ACHE), also known as acetylcholine acetylhydrolase, is found in erythrocytes, in the lungs, spleen and in grey matter of the brain. The pseudo-cholinesterase (CHE), also referred to as acylcholine acylhydrolase, is found in serum, the liver, pancreas, heart and in the white matter of the brain. The assay of serum cholinesterase (CHE) is useful to diagnose: liver disorders, hepatitis, cirrhosis, carcinoma with metastasis, sensitivity to succinylcholine administration and pesticide poisoning. Levels decrease in all of the diseases above.

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The method uses butyrylthiocholine as the specific substrate for cholinesterase (CHE). Cholinesterase catalyses the hydrolysis of butyrylthiocholine substrate forming butyrate and thiocholine, in presence of potassium hexacyanoferrate (III). Thiocholine reduces potassium hexacyanoferrate (III) (yellow colour) to potassium hexacyanoferrate (II) (colourless). The decrease in absorbance is directly proportional to the CHE activity in the sample.

### REAGENTS

#### Package

	Cat. No 4-256 (24-TRAY)	Cat. No 4-456 (36-TRAY)
1-Reagent	2 x 40 ml	2 x 23 ml
2-Reagent	2 x 12.5 ml	2 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks..

### Concentrations in the test

<b>1-Reagent</b>	
pyrophosphate buffer, pH 7.7	65 mmol/l
hexacyanoferrate (III)	2 mmol/l
<b>2-Reagent</b>	
Good's Buffer pH 4.0	20 mmol/l
butyrylthiocholine iodide	65 mmol/l

### Warnings and notes

- Protect from direct sunlight, evaporation and avoid contamination!
- Do not freeze reagents
- Immediately after use, recap the bottles and store at 2-8°C.
- Do not use after expiry date.
- Do not interchange caps.
- Reagents should be mixed before use by gentle inverting the bottle several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of reagent instability.

### SPECIMEN

Fresh serum free from haemolysis, plasma (EDTA, heparin) not hemolyzed. Serum / plasma should be separated from red blood cells as soon as possible after blood collection. Do not use citrate, borate, oxalate and fluoride as an anticoagulant because it inhibits cholinesterase activity.

It is recommended to follow CLSI procedures regarding specimen collecting and handling.

Sample may be stored for up to 15 days at 2-8°C or 12 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CHOLINESTERASE – BIL TOTAL. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.

### REFERENCE VALUES <sup>4,5</sup>

serum / plasma	37°C	
female	4000 – 12600 U/l	67 – 210 µkat/l
male	5100 – 11700 U/l	85 – 195 µkat/l

In infants up to 6 months of age, cholinesterase activity is 40% to 50% higher than in adults. In young adult (< 35 years) women, the enzyme activity is approximately 64% to 74% of that in adult males. The activity decreases during pregnancy.

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers: Prestige 24i, Biolis 24i Premium the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number.

For the calibration of Biolis 30i automatic analyser the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 29 U/l (0.48 µkat/l)
- LoD (Limit of Detection):** 53 U/l (0.88 µkat/l)
- LoQ (Limit of Quantitation):** 69 U/l (1.15 µkat/l)
- Linearity:** up to 16200 U/l (270 µkat/l)

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	4274	35.84	0.84
level 2	7793	81.31	1.04
Reproducibility (day to day) n = 80	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	4161.7	70.29	1.7
level 2	7734.2	125.47	1.6

#### Method comparison

A comparison between cholinesterase values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 67 serum samples gave following results:

$$y = 0.9073 x + 903.01 \text{ U/l};$$

$$R = 0.986 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	58	Item Name	CHE	Optical
<b>Data information</b>				
Units	U/L			
Decimals	0			
<b>Calibration</b>				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
<b>Analysis</b>				
Type	RATE method			
Main Wave Length	405 nm			
Sub Wave Length	660 nm			
Method	DGKC			
<b>Correlation</b>				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	58	Item Name	CHE	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	3	5	µl	
Reagent 1	200	10	µl	
Reagent 2	50	10	µl	
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	38	48	
<b>Abs.Limit</b>				
Low	0.6			High
				2.5
<b>Correction value</b>				
Blank correction	1			
End Point Limit	2.5			
Linear Check (%)	15			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				

Item No.	58	Item Name	CHE	Optical
<b>Normal Range</b>				
	Male	Female		
	Low	High	Low	High
Serum	5100	11700	4000	12600
Urine				
Plasma	5100	11700	4000	12600
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male	Female		
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	58	Item Name	CHE	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	OFF		
	Higher	OFF		
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum			223	21000
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:DiI2				

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