

CORMAY IgE CALIBRATORS

CALIBRATORS FOR DETERMINATION OF IgE LEVELS



Kit name	Kit size	Cat. No
CORMAY IgE CALIBRATORS	5 x 2 ml	4-280

INTRODUCTION

IgE is an immunoglobulin with a molecular weight of approximately 190 kD normally present in the blood in trace amounts. Continual production of IgE antibodies in response to common naturally occurring allergens, however, often results in elevated serum levels and in the development of such clinically important Type I allergic reactions as asthma, hay fever, dermatitis and food allergies. Elevated IgE levels are also seen in parasitic (helminth) diseases, IgE myeloma, and in hepatitis. The measurement of IgE in human serum is thus considered to be useful in the diagnosis, treatment, assessment of disease progression, or postoperative prognosis for such conditions.

CORMAY IgE CALIBRATORS kit is prepared by diluting IgE with buffer solution containing 1 w/v% bovine serum albumin at various concentrations. Concentrations are shown on each vial and in this kit insert.

CORMAY IgE CALIBRATORS kit is dedicated for calibration of IgE determinations on automatic analysers.

PACKAGE

Cap colour	IgE level	Size
CREAM	50 IU/ml	1 x 2 ml
GREEN-YELLOW	100 IU/ml	1 x 2 ml
BLUE	200 IU/ml	1 x 2 ml
PINK	500 IU/ml	1 x 2 ml
DARK RED	1000 IU/ml	1 x 2 ml

The calibrators are ready to use.

STABILITY AND STORAGE

Unopened calibrator is stable at 2-10°C until the expiration date printed on the vial. Protect from light.

WARNINGS AND NOTES

- Product for in vitro diagnostic use only.
- As the colours of caps vary according to IgE levels in the vial, care should be taken not to interchange the caps.
- The calibrators should not be mixed with ones in the other vial.
- The calibrator should be used according to this kit insert.
- The calibrator is based on human serum that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Product contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.

TRACEABILITY

IgE calibrators are prepared referring to the WHO standardization.

WASTE MANAGEMENT

Please refer to local legal requirements.

Date of issue: 02.2012.

MANUFACTURER

PZ CORMAY S.A.
ul. Wiosenna 22,
05-092 Łomianki, POLAND
tel.: +48 (0) 22 751 79 10
fax: +48 (0) 22 751 79 14
<http://www.pzcormay.pl>