## CORMAY IMMUNO-CONTROL I CORMAY IMMUNO-CONTROL II

(PL)

### CORMAY IMMUNO- CONTROL I

Nr kat. 4-288 2 x 3 ml

### CORMAY IMMUNO- CONTROL II

Nr kat. 4-290 2 x 3 ml

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Cat. No 4-290 2 x 3 ml

#### INTENDED USE

CORMAY IMMUNO-CONTROL I (high level-H, low level L) and CORMAY IMMUNO-CONTROL II (high level-H, low level L) sera are lyophilised, pooled human sera which are as a multianalyte control sera when measuring various serum proteins by automatic analysers with turbidimetric capabilities. Assigned values for both the high and low controls vary from lot to lot and are given in values sheet.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

#### RECONSTITUTION

Reconstitute the contents of one vial with 3 ml of distilled water. Allow to stand at room temperature for 20 minutes, gently swirling occasionally to facilitate dissolution. Once reconstituted treat as human serum specimen.

#### STABILITY AND STORAGE

1. Lyophilised serum remains stable when stored at 2-10°C until expiry date given on the product label. Protect from light and contamination.
2. Stability after reconstitution, at 2-8°C: 4 weeks (C-reactive protein, anti-streptolysin O, rheumatoid factor, ferritin, total IgE, myoglobin).

#### NOTES

1. Once reconstituted control sera are ready for use.
2. Remove the cap carefully in order not to lose the powder in the bottle.
3. The controls are based on human serum that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as thought capable of transmitting infectious disease.

#### WASTE MANAGEMENT

Please refer to local legal requirements.

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