PRESTIGE 24i COMPLEMENT C4

Nr kat. 4-310, 4-336 (PL)

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Cat. No **4-310, 4-336** (EN)

INTENDED USE

Diagnostic kit for determination of complement C4 concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Complement is a group of 20 immunologically distinct proteins present in blood and tissues. They are able to interact with antigen-antibody complexes, with each other and with cell membranes, in a complex way intended at destroying viruses and bacteria. They are synthesised in liver and are present in serum as functionally inactive molecules. They are activated by antigen-antibody complexes. C4 complement is a α -glycoprotein of 3 subunits. It is an acute phase reactant whose levels are increased during the acute phase. Low levels are found in immune complex diseases and in inherited angioedema, while C3 complement levels are normal.

METHOD PRINCIPLE

The complement C4 present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to complement C4 concentration in the sample.

REAGENTS

Package

	Cat. No 4-310 (24-TRAY)	Cat. No 4-336 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-Reagent: anti-human complement C4 antiserum, inorganic salt, buffer, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma.

Serum should be separated from red blood cells as soon as possible after blood collection. If the test cannot be done immediately, the sample should be stored at -70°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ²

adults	0.1 – 0.4 g/l
--------	---------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzer Hitachi and Biolis 30i. Results may vary if a different instrument is used.

- LoQ (limit of quantitation):** 0.002 g/l

- Linearity:** up to 0.825 g/l.

- Specificity / Interferences**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl triglycerides up to 155.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.147	0.001	0.77
level 2	0.321	0.003	0.76
level 3	0.462	0.003	0.64

Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.137	0.005	3.48
level 2	0.321	0.008	2.52
level 3	0.472	0.010	2.08

- Method comparison**

A comparison between complement C4 values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 64 serum samples gave following results:

$y = 1.048x - 0.0237$ g/l;

$R = 0.997$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	35	C4
Auto Rerun SW		
OFF		
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦІЯ: Prestige 24i Premium, Biolis 24i Premium

Item No.	35	Item Name	C4	Optical
Data information				
Units	g/l			
Decimals	2			
Calibration				
Type	Logit2			
Std sample conc.				
Blank	0	#1	*	#2 *
#3	*	#4	*	#5 *
#6				
Analysis				
Type	END method			
Main Wave Length	340 nm			
Sub Wave Length	700 nm			
Method	Immuno			
Correlation				
	Slope	Intercept		
Y=	1	X+	0	

Item No.	35	Item Name	C4	Optical
Aspiration				
Kind	Double			
Data Process				
Read	Main	Start	End	
	Sub	30	31	
Vol.				
	Kind	Vol.	Add	Units
Sample		7	5	µl
Reagent 1		250	10	µl
Reagent 2		50	10	µl
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Correction value				
Blank correction	0.8371			
End Point Limit	2			
Linear Check (%)	0			
Reaction Monitor				
0 Level Point	1			
Span	3			
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Third mixing				
ON				

Item No.	35	Item Name	C4	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.1	0.4	0.1	0.4
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	35	Item Name	C4	Optical			
Auto Rerun SW							
OFF							
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Absorbance)							
OFF							
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
99:Dil1							

PRESTIGE 24i COMPLEMENT C4

• **Biolis 30i**

Item no	35	Item name	C4	Specimen	SERUM / PLASMA	OPTICAL												
Data information																		
UNITS	g/L		Aspiration volume															
DECIMALS	2		TYPE Double															
Analysis																		
METHOD	END method		<table border="1"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>7</td> <td>250</td> <td>50</td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>					SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	7	250	50	BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	7	250	50															
BOTTLE (ml)																		
Main Wave Length	340 nm		FIRST DIL.															
Sub Wave Length	700 nm																	
CORRELATION (Y= AX + B)																		
A =	1		Data processing read															
B =	0		<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> </tr> <tr> <td>MAIN</td> <td>51</td> <td>52</td> </tr> <tr> <td>SUB</td> <td>30</td> <td>31</td> </tr> </table>					START	END	MAIN	51	52	SUB	30	31			
	START	END																
MAIN	51	52																
SUB	30	31																
Blank value																		
• WATER ° REAGENT																		
Calibration																		
TYPE	Logit 2		ABS LIMIT															
STABILITY			-3 TO 3															
Collection value																		
END POINT			2.5															
LINEARITY CHECK (%)			0															
Prozone check																		
		START	END	LIMIT (%)														
FIRST																		
SECOND																		
°HIGH				MINIMUM ABS.														
•LOW				MEAN														
				VARIATE														
Item No	35	Item Name	C4	Specimen	SERUM / PLASMA	OPTICAL												
Reference intervals																		
MALE				FEMALE														
LOW	HIGH	LOW	HIGH															
0.1	0.4	0.1	0.4															
Panic range																		
MALE				FEMALE														
LOW	HIGH	LOW	HIGH															
Decision limit																		
MALE				FEMALE														
Reaction check																		
		°ON	•OFF															
CHECK																		
LOW																		
HIGH																		
VL CHECK																		
°ON	•OFF		VH CHECK															
			°ON •OFF															

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