## PRESTIGE 24i HAPTOGLOBIN

Nr kat. 4-311, 4-337

(PL)

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Cat. No **4-311, 4-337** (EN)

### INTENDED USE

Diagnostic kit for determination of haptoglobin concentration used in automatic analysers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Haptoglobin is an acute phase protein whose primary function consists in binding free haemoglobin in serum. The complex is removed within minutes by the reticulo-endothelial system where its components are metabolised to free aminoacids and iron. Haptoglobin consequently plays a major role in preventing the loss of haemoglobin in urine and the consequent iron loss from the iron pool. The haptoglobin levels are increased during the acute phase and in such conditions as burns or nephrotic syndrome. Haptoglobin levels are abnormally high in intravascular hemolysis and when haemoglobin turnover is increased such as during haemolytic anaemia, transfusion reactions and malaria.

### METHOD PRINCIPLE

The haptoglobin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at  $\lambda=340$  nm is proportional to haptoglobin concentration in the sample.

### REAGENTS

#### Package

	Cat. No 4-311 (24-TRAY)	Cat. No 4-337 (36-TRAY)
1-REAGENT	1 x 40 ml	2 x 23 ml
2-REAGENT	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

### Reagent components

**1-Reagent:** imidazole buffer, polymer, inorganic salt, buffer, preservative.

**2-Reagent:** anti-human haptoglobin antiserum, inorganic salt, buffer, preservative.

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.

- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum. Nonhemolyzed, fresh serum is recommended. Samples may be stored up to 2 weeks at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
 1-Reagent put on basic position in reagent tray.  
 2-Reagent put on start position in reagent tray.  
 For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES

adults	0.26 – 1.85 g/l
newborns	0.05 – 0.48 g/l
> 60 year M	0.35 – 1.64 g/l
> 60 year F	0.40 – 1.75 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems, the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers Hitachi or Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LOQ (Limit of Quantitation):** 0.002 g/l
- Linearity:** up to 3.86 g/l.
- Interferences:**  
 Hemoglobin up to 0.06 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 1669.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

### Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.49	0.007	1.44
level 2	1.31	0.008	0.63
level 3	2.14	0.017	0.77
Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.62	0.018	2.84
level 2	1.48	0.037	2.48
level 3	2.27	0.050	2.21

### Method comparison

A comparison between haptoglobin values determined at **Biolis 30i** (y) and at **BS-800** (x) using 61 serum samples gave following results:

$$y = 1.0131 x + 0.0082 \text{ g/l};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	36	HAPTO
<b>Auto Rerun SW</b>		
OFF		
<b>Auto Rerun Range (Result)</b>		
	Lower	Higher
	OFF	OFF
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range		
	Lower	OFF
	Higher	OFF
<b>Auto Rerun Condition (Prozone)</b>		
Prozone Range		
OFF		

#### • Prestige 24i Premium, Biolis 24i Premium

Item No.	36	Item Name	HAPTO	Optical
<b>Data information</b>				
Units	g/l			
Decimals	2			
<b>Analysis</b>				
Type	END method			
Main Wave Length	340nm			
Sub Wave Length	700nm			
Method	Immunoturbidimetry			
<b>Calibration</b>				
Type	Logit2			
Std sample conc.				
Blank	0	#1 *	#2 *	#3 *
#3 *	#4 *	#5 *	#6 *	
#6				
<b>Correlation</b>				
Y=	Slope	Intercept		
	1	0		

Item No.	36	Item Name	HAPTO	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		2.5	5	µl
Reagent 1		250	10	µl
Reagent 2		50	10	µl
<b>Data Process</b>				
Read				
	Main	Start	End	
	Sub	30	31	
<b>Abs.Limit</b>				
Low	-3	High	3	
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
ON				
<b>Correction value</b>				
Blank correction				0.8350
End Point Limit				2
Linear Check (%)				0
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	36	Item Name	HAPTO	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0.26	1.85	0.26	1.85
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	36	Item Name	HAPTO	Optical			
<b>Auto Rerun SW</b>							
OFF							
<b>Auto Rerun Range (Conc.)</b>							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
<b>Auto Rerun Condition (Absorbance)</b>							
OFF							
<b>Auto Rerun Condition (Prozone)</b>							
OFF							
<b>Dilution</b>							
99:Dil1							

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