## PRESTIGE 24i ALPHA 1-ANTITRYPSIN

Nr kat. 4-315, 4-341 (PL)

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The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\alpha$  1-antitrypsin is an acute phase protein showing anti-protease activity. The main function of this protein consists in neutralising lysosomal elastase released upon phagocytosis by polymorphonuclear leukocytes. Inherited deficiency of the protein is associated with lung and liver diseases. Low levels are encountered in neonatal respiratory distress syndrome and in severe protein losing disorders. Increased levels are more common, particularly during the acute phase.

### METHOD PRINCIPLE

The  $\alpha$  1-antitrypsin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured is proportional to  $\alpha$  1-antitrypsin concentration in the sample.

### REAGENTS

#### Package

	Cat. No 4-315 (24-TRAY)	Cat. No 4-341 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

### Reagent components

**1-Reagent:** TRIS buffer, polymer, inorganic salt, preservative.

**2-Reagent:** anti-human  $\alpha$  1-antitrypsin antiserum, buffer, inorganic salt, preservative.

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product

### SPECIMEN <sup>3</sup>

Serum or plasma (EDTA).

Sample may be stored up to 24 hours at 20-25°C, up to 7 days at 4-8°C or up to 3 months at -70°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Application guide for Biolis 30i analyser is available on request.

### REFERENCE VALUES <sup>3</sup>

	mg/dl	g/l
newborns	145 - 270	1.45 - 2.70
adults	78 - 200	0.78 - 2.00
> 60 year	115 - 200	1.15 - 2.00

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi. Results may vary if a different instrument or a manual procedure is used.

- LOQ (Limit of Quantitation):** 0.685 mg/dl

- Linearity:** up to 350 mg/dl.

- Interferences:**

Hemoglobin up to 0.32 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 1669.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	58.3	0.65	1.12
level 2	117.7	1.01	0.85
level 3	176.7	1.97	1.11
Reproducibility (day to day) n = 66	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	59.3	0.96	1.61
level 2	117.9	1.48	1.25
level 3	177.2	2.82	1.59

### Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using at least 30 samples gave following results:

$$y = 1.05x - 0.02 \text{ mg/dl};$$

$$R = 0.994 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	40	ATT
<b>Auto Rerun SW</b>		
OFF		
<b>Auto Rerun Range (Result)</b>		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
<b>Prozone Range</b>		
		OFF

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	40	Item Name	ATT	Optical
<b>Data information</b>				
Units	g/l			
Decimals	2			
<b>Analysis</b>				
Type	END method			
Main Wave Length	340nm			
Sub Wave Length	700nm			
Method	Immuno			
<b>Calibration</b>				
Type	Logit2			
Std sample conc.				
Blank	0	#1 *	#2 *	#3 *
#3	*	#4 *	#5 *	#6 *
#6				
<b>Correlation</b>				
Y=	Slope	Intercept		
	1	0		

Item No.	40	Item Name	ATT	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	3	5	µl	
Reagent 1	250	10	µl	
Reagent 2	50	10	µl	
<b>Data Process</b>				
Read	Start	End		
Main	44	45		
Sub	30	31		
<b>Abs.Limit</b>				
Low	-3	High	3	
<b>Correction value</b>				
Blank correction	0.8350			
End Point Limit	2			
Linear Check (%)				
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
ON				

Item No.	40	Item Name	ATT	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0.78	2.00	0.78	2.00
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	40	Item Name	ATT	Optical
<b>Auto Rerun SW</b>				
OFF				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	Higher	OFF	
			OFF	
<b>Auto Rerun Condition (Prozone)</b>				
	OFF			
<b>Dilution</b>				
	99:Dil1			

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