## PRESTIGE 24i LQ LIPASE

Nr kat. **4-185, 4-329** (PL)

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### LITERATURA

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

Data wydania: 04.2021.

## PRESTIGE 24i LQ LIPASE

Cat. No **4-185, 4-329** (EN)

### INTENDED USE

Diagnostic kit for determination of lipase activity intended to use in automatic analysers: analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Lipase is a digestive enzyme released into the intestine from the pancreas where it breaks down triglycerides into fatty acids and glycerol prior to absorption. Lipase measurements are used in the diagnosis and treatment of diseases of the pancreas such as acute pancreatitis, obstruction of the pancreatic duct and pancreatic tumours.

### METHOD PRINCIPLE

The colorimetric method is based on a lipase specific degradation of a chromogenic substrate. The specific lipase substrate-DGGMR [1,2-o-dilauryl-racglycero-3-glutaric acid-(6'-methylresorufin) ester] is cleaved by the catalytic action of lipase to form 1,2-o-dilauryl-racglycerol and an unstable intermediate, glutaric acid-(6-methyl resorufin) ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The lipase activity in the specimen is proportional to the production of methylresorufin in the reaction and can be determined photometrically.

### REAGENTS

Package	Cat. No 4-185 (24-TRAY)	Cat. No 4-329 (36-TRAY)
1-Reagent	4 x 38 ml	4 x 23 ml
2-Reagent	4 x 20 ml	4 x 12.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 7 weeks, Biolis 24i Premium – 7 weeks.

### Concentrations in the test

#### 1-Reagent

TAPS [N-Tris(hydroxymethyl)methyl-3-aminopropanesulfonic acid]	100 mM
sodium hydroxide	40 mM
sodium deoxycholate preservative	34 mM

#### 2-Reagent

tartaric acid	9.5 mM
sodium hydroxide	19 mM
colipase	460 IU/ml
2-propanol	0.65 M
DGGMR [1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)-ester]	0.4 mM

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, heparinized plasma free from hemolysis. Sample may be stored for up to 5 days at 2-8°C or 24 hours at 20-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LIPASE – TG MONO, CHOL – LIPASE, LDL DIRECT – LIPASE, TG - LIPASE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LIPASE – CALCIUM ARSENAZO, CHOL – LIPASE, HDL DIRECT – LIPASE, LDL DIRECT – LIPASE, TG - LIPASE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

Normal range	13 – 60 U/l	0,22 – 1,00 µkat/l
--------------	-------------	--------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 7 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 0.4 U/l (0.007 µkat/l).

- LoD (Limit of Detection):** 0.7 U/l (0.012 µkat/l).
- LoQ (Limit of Quantitation):** 5 U/l (0.083 µkat/l).
- Linearity :** up to 627 U/l (10.45 µkat/l).

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

#### Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	58.7	0.58	1.00
level 2	103.4	0.88	0.85
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	54.2	2.13	3.9
level 2	97.1	3.25	3.3

#### Method comparison

A comparison between lipase values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 62 serum samples gave following results:

$$y = 0.9951x + 2.5134 \text{ U/l};$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

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Item name	26	LIPA				
<b>Auto Rerun SW</b>						
ON						
<b>Auto Rerun Range (Result)</b>						
	ON	ON				
	Lower	Higher				
Serum	7	600				
Urine						
Plasma						
CSF						
Dialysis						
Other						
<b>Auto Rerun Condition (Absorbance)</b>						
Absorbance Range	Lower	Higher	OFF			
			OFF			
<b>Auto Rerun Condition (Prozone)</b>						
Prozone Range	OFF					

### ■ Prestige 24i Premium, Biolis 24i Premium

Item No.	26	Item Name	LIPA	Optical	
<b>Data information</b>					
Units	U/l				
Decimals	1				
<b>Calibration</b>					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
<b>Analysis</b>					
Type	RATE method				
Main Wave Length	570 nm				
Sub Wave Length	700 nm				
Method	Colorimetric				
<b>Correlation</b>					
Y=	Slope	Intercept			
	1	0			

Item No.	26	Item Name	LIPA	Optical	
<b>Aspiration</b>					
Kind	Double				
Vol.	Kind	Vol.	Add	Units	
	Sample	3	5	µl	
	Reagent 1	200	10	µl	
	Reagent 2	100	10	µl	
<b>Data Process</b>					
Read	Main	Start	End		
	Sub	35	42		
<b>Abs.Limit</b>					
Low	-3	High	3		
<b>Correction value</b>					
Blank correction					
End Point Limit					
Linear Check (%)					
2					
80					
<b>Prozone Check</b>					
First	Start	End	Limit (%)		
Second			Low		

Item No.	26	Item Name	LIPA	Optical			
<b>Normal Range</b>							
	Male		Female				
	Low	High	Low	High			
Serum	13	60	13	60			
Urine							
Plasma	13	60	13	60			
CSF							
Dialysis							
Other							
<b>Panic Range</b>							
	Male		Female				
	Low	High	Low	High			
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							

Item No.	26	Item Name	LIPA	Optical			
<b>Auto Rerun SW</b>							
ON							
<b>Auto Rerun Range (Conc.)</b>							
	First Dil	Low	High				
		Re	Value	Dil	Re	Value	Dil
Serum			6			1000	
Urine							
Plasma							
CSF							
Dialysis							
Other							
<b>Auto Rerun Condition (Absorbance)</b>							
	Lower	Higher	OFF				
			OFF				
<b>Auto Rerun Condition (Prozone)</b>							
	OFF						
<b>Dilution</b>							
	100:Dil2						

## PRESTIGE 24i LQ LIPASE

▪ **Biolis 30i**

Item no	26	Item name	LIPASE	Specimen	SERUM / PLASMA	OPTICAL
<b>Data information</b>				<b>Aspiration volume</b>		
UNITS		U/L		TYPE		
DECIMALS		1		Double		
<b>Analysis</b>				SAMPLE REAGENT 1 REAGENT 2		
METHOD		RATE method		VOL. (µL)		
Main Wave Length		570 nm		3 200 100		
Sub Wave Length		700 nm		BOTTLE (ml)		
<b>CORRELATION (Y= AX + B)</b>				FIRST DIL.		
A =		1		Data processing read		
B =		0		START END		
<b>Blank value</b>				MAIN 37 53		
• WATER ° REAGENT				SUB		
<b>Calibration</b>				ABS LIMIT		
TYPE		Linear 1		-3 TO 3		
STABILITY				Collection value		
				END POINT 2.5		
				LINEARITY CHECK (%) 80		
				Prozone check		
				START END LIMIT (%)		
				FIRST		
				SECOND		
				MINIMUM ABS.		
				°HIGH MEAN		
				•LOW VARIATE		

Item No	26	Item Name	LIPASE	Specimen	SERUM / PLASMA	OPTICAL
<b>Reference intervals</b>						
MALE		FEMALE		Auto rerun		
LOW	HIGH	LOW	HIGH	•ON °OFF		
13	60	13	60	Auto rerun range (conc.)		
<b>Panic range</b>						
MALE		FEMALE		Re	Value	Dil.
LOW	HIGH	LOW	HIGH	5	627	
<b>Decision limit</b>						
MALE		FEMALE		Auto rerun condition (abs.)		
				DIL.		
				LOWER	°ON	•OFF
				HIGH	°ON	•OFF
<b>Reaction check</b>						
				Auto rerun condition (prozone)		
				°ON •OFF		
				SAMPLE VOL.		
<b>Dilution</b>						
				•DIL 1 ° DIL 2		
<b>VL CHECK</b>						
				°ON •OFF		
<b>VH CHECK</b>						
				°ON •OFF		