PRESTIGE 24i LQ BILE ACIDS

Nr kat. **4-345, 4-495** (PL)

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Cat. No 4-345, 4-495 (EN)

INTENDED USE

Diagnostic kit for determination of total bile acids concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.


INTRODUCTION

Bile acids are the main product of degradation of endogenous cholesterol formed in the liver. Total bile acids are metabolized in the liver and are a valuable indicator of normal or abnormal liver function. Serum total bile acids are increased in patients with viral hepatitis, liver cirrhosis and liver cancer.

METHOD PRINCIPLE

Enzymatic method with 3- α -hydroxysteroid dehydrogenase (3- α -HSD).

Bile acids under the influence of 3-hydroxysteroid dehydrogenase (3- α -HSD) in the presence of thio-NAD are converted to 3-ketosteroids and thio-NADH. The reaction is reversible and 3- α -HSD can convert 3-ketosteroids and NADH to bile acids and NAD.

The rate of thio-NADH formation can be monitored at 405 nm and is proportional to the bile acids activity.

REAGENTS

Package

	Cat. No 4-345 (24-TRAY)	Cat. No 4-495 (36-TRAY)
1-Reagent	2 x 16.5 ml	2 x 15.5 ml
2-Reagent	2 x 6.5 ml	2 x 6 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 7 weeks, Biolis 24i Premium – 7 weeks.

Concentrations in the test

1-Reagent	
Thio-NAD	> 0.1 mmol
Buffer	

2-Reagent

3- α -HSD	> 2 kU/l
NADH	> 0.1 mmol
Buffer	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Avoid contact with skin and mucous membranes.
- Reagents from different lots must not be interchanged.
- Yellow or yellow-brown color of the reagent does not affect the reagents performance.
- Samples from patients treated with ursodeoxycholic acid (UDCA) are not suitable for the determination of total bile acid concentrations.

SPECIMEN

Serum.

Total bile acids concentration is increased after meals, therefore samples should be collected under fasting conditions. Serum samples are stable for 7 days at 4°C or for 3 month at -20 °C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionised water is recommended.

REFERENCE VALUES ³

serum	2.5 – 6.8 μ mol/l (1.25 – 3.4 μ g/ml)
-------	---

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY BILE ACIDS CONTROLS (Cat. No 5-149).

For the calibration of automatic analysers systems the CORMAY BILE ACIDS CALIBRATOR (Cat. No 3-125) is recommended.

The calibration curve should be prepared every 7 weeks (Prestige 24i, Biolis 24i Premium) with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i or Hitachi 717. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.91 μ mol/l (0.46 μ g/ml)

LoD (Limit of Detection):

1.30 μ mol/l (0.65 μ g/ml)

LoQ (Limit of Quantitation):

4 μ mol/l (2 μ g/ml)

Linearity:

up to 146 μ mol/l (73 μ g/ml)

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 50 mg/dl, ascorbic acid up to 50 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [μ mol/l]	SD [μ mol/l]	CV [%]
level 1	27.68	0.42	1.52
level 2	100.61	1.24	1.23
Reproducibility (day to day) n = 80	Mean [μ mol/l]	SD [μ mol/l]	CV [%]
level 1	27.8	1.11	4.0
level 2	103.2	3.71	3.6

Method comparison

A comparison between bile acids values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$y = 0.9468x + 0.3185 \mu$ mol/l;

$R = 1.000$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATORY / APPLICATION for / АДАПТАЦИЯ для:

• **Prestige 24i, Biolis 24i**

Item name	72	Bile Acids
Data information		
Units	$\mu\text{mol/l}$	
Decimals	2	
Analysis		
Type	RATE	
Main W.Length1	405	
Sub W.Length2	700	
Method	Enzymatic	
Calibration		
Type	Linear	
Standard		
#1	*	#4
#2		#5
#3		#6
Normal Range		
	Male	Female
	Low High	Low High
Serum	2.5 6.8	2.5 6.8
Urine		
Plasma		
CSF		
Dialysis		
Other		
Corr		
Y=	Slope	Inter
	1.000	0.000

Item name	72	Bile Acids
Aspiration		
Kind	Double	
Vol.		
Kind	Vol.	Add
Sample	3	5
Reagent1	270	10
Reagent2	90	10
Data Process		
Read	Start	End
Main	35	50
Sub		
Absorbance Limit		
Low	-3.000	
High	3.000	
Factor		
Blank correction	Endpoint Limit	2.000
	Linear Check (%)	0
Dilution		
Diluent	100:DiI2	
Prozone Check		
First	Start	End
Second		Low
Third		Low
Monitor		
0 Level Point	1	
Span	3.000	
Third Mix.		
Third Mix.	OFF	
R1 Blank		
R1 Blank	Water-B	

Item name	72	Bile Acids
Auto Rerun SW		
Auto Rerun SW	OFF	
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum	4.9	170
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	72	Item Name	Bile Acids	Optical
Data information				
Units	$\mu\text{mol/l}$			
Decimals	2			
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	RATE			
Main Wave Length	405 nm			
Sub Wave Length	700 nm			
Method	Enzymatic			
Correlation				
Slope	Intercept			
Y=	1	X+	0	

Item No.	72	Item Name	Bile Acids	Optical
Aspiration				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	2	5	μl	
Reagent 1	250	10	μl	
Reagent 2	50	10	μl	
Data Process				
Read	Main	Start	End	
	Sub	35	50	
Abs.Limit				
Low	High			
-3	3			
Correction value				
Blank correction				
End Point Limit	2			
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				

Item No.	72	Item Name	Bile Acids	Optical
Normal Range				
	Male	Female		
	Low High	Low	High	
Serum	2.5 6.8	2.5	6.8	
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male	Female		
	Low High	Low	High	
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	72	Item Name	Bile Acids	Optical	
Auto Rerun SW					
Auto Rerun SW	OFF				
Auto Rerun Range (Conc.)					
	First Dil	Low	High		
		Re Value	Dil	Re Value	Dil
Serum		2.9		150	
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Absorbance)					
	Lower	OFF			
	Higher	OFF			
Auto Rerun Condition (Prozone)					
	OFF				
Dilution					
	100:DiI2				

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