PRESTIGE 24i LQ ASAT

Nr kat. **4-214, 4-414** (PL)

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LITERATURA

1. Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
2. Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
3. Bergmeyer H.U., IFCC methods for the measurement of catalytic concentrations of enzymes: Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2), Clin. Chem. Acta 105, 147-172 (1980).
4. Bergmeyer HU, Hørdler M, Rej R. International Federation of Clinical Chemistry (IFCC) Scientific Committee, Analytical Section: approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2). J Clin Chem Clin Biochem. 1986 Jul;24(7):481-95. PMID: 3734711.
5. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 911 (1989).
6. NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition. NCCLS document H18-A3 (ISBN 1-56238-555-0). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA (2004).
7. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 413, 1800 (2018).
8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
9. Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2257, (2006).

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
Cat. No **4-214, 4-414** (EN)

INTENDED USE

ASAT reagent is intended to determine quantitatively aspartate aminotransferase (ASAT) activity in serum. It is used for monitoring and as an aid to diagnosis of clinical conditions associated with abnormal aspartate aminotransferase activity. ASAT reagent is intended to use on automatic analysers Biolis 24i Premium and Biolis 30i. It is only for *in vitro* diagnostics, for healthcare professional users.

INTRODUCTION^{1,2}

Determinations of elevated aspartate aminotransferase activity (ASAT, AST, GOT) in serum are mainly used to aid in the diagnosis and monitor of liver diseases e.g. hepatitis, hepatic cirrhosis, drug-induced liver injury, hepatic metastasis, mononucleosis. Also skeletal muscle diseases such as muscle trauma, surgery, burns, muscular dystrophy, heat stroke and other diseases such as acute hemolytic anemia, acute pancreatitis and patients with higher risk of COVID-19 adverse outcome. Decreased levels of ASAT can be observed in acute renal disease, beriberi, diabetic ketoacidosis, pregnancy, chronic renal dialysis.

METHOD PRINCIPLE^{3,4}

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

L-aspartate + 2-oxoglutarate $\xleftarrow{\text{ASAT}}$ oxalacetate + L-glutamate

oxalacetate + NADH + H⁺ $\xleftarrow{\text{MDH}}$ malate + NAD⁺

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to aspartate aminotransferase activity.

REAGENTS

	Cat. No 4-214 (24-TRAY)	Cat. No 4-414 (36-TRAY)
1-REAGENT	6 x 40 ml	8 x 23 ml
2-REAGENT	6 x 12.5 ml	8 x 7.5 ml
Number of tests		
Biolis 24i Premium	1070	810
Biolis 30i	1080	810

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT	
L-aspartate	300 mmol/l
LDH	1.95 U/ml
MDH	1.15 U/ml
buffer Tris	
pH adjuster	

stabilizers	
preservative	
2-REAGENT	
2-oxoglutarate	62 mmol/l
NADH	1.4 mmol/l
buffer	
pH adjuster	
preservatives	

REAGENT STABILITY

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C (Biolis 24i Premium, Biolis 30i).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

 H315 Causes skin irritation.
H319 Causes serious eye irritation.
P280 Wear protective gloves, protective clothing, eye protection or face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN^{5,6,7,8}

Serum free from hemolysis.
Hemolysis should be avoided, since ASAT activity in erythrocytes is 15 times higher than in normal serum.
ASAT activity remains stable in specimen up to 4 days at 15-25°C, up to 7 days at 2-8°C or up to 3 months at -20°C.
Nevertheless it is recommended to perform the assay with freshly collected samples!
Follow tube manufacturers' instructions carefully when using collection tubes.
Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

PROCEDURE

Reagents are ready to use.
1-REAGENT put on basic position in reagent tray.
2-REAGENT put on start position in reagent tray.
For reagent blank deionized water is recommended.

Actions required:

- **Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: ASAT II GEN – LDH. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

Application guide for Prestige 24i analyser is available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks (Biolis 24i Premium, Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ⁹

serum	37°C	
female	up to 31 U/l	up to 0.53 µkat/l
male	up to 35 U/l	up to 0.60 µkat/l

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**
1.8 U/l (0.03 µkat/l) – BS-400

- **LoD (Limit of Detection):**

3.0 U/l (0.05 µkat/l) – BS-400

- **LoQ (Limit of Quantitation):**

9.0 U/l (0.15 µkat/l) – Biolis 30i

- **Linearity:**

up to 770 U/l (12.8 µkat/l) – Biolis 30i

- **Measurement range**

9,0 U/l (0,15 µkat/l) - 770 U/l (12,8 µkat/l) – Biolis 30i

- **Specificity / Interferences**

Ascorbic acid up to 62 mg/l, triglycerides up to 1000 mg/dl, bilirubin up to 20 mg/dl, haemoglobin up to 0.31 g/dl (in samples containing low activity of ASAT) and up to 2.5 g/dl (in samples containing high activity of ASAT) do not interfere with the test.

- **Precision (Biolis 30i)**

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	44.6	1.56	3.5
level 2	200.4	1.33	0.7
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	40.6	1.36	3.4
level 2	210.4	1.99	0.9

- **Method comparison**

A comparison between ASAT values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 59 serum samples gave following results:

$$y = 1.0272 x + 1.1659 \text{ U/l}$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT ¹⁰

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ¹¹

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

1. Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
2. Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
3. Bergmeyer H.U., IFCC methods for the measurement of catalytic concentrations of enzymes: Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2), Clin. Chem. Acta 105, 147-172 (1980).
4. Bergmeyer HU, Hørder M, Rej R. International Federation of Clinical Chemistry (IFCC) Scientific Committee, Analytical Section: approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2). J Clin Chem Clin Biochem. 1986 Jul;24(7):481-95. PMID: 3734711.
5. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 911 (1989). MAT BIOLOG
6. NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition. NCCLS document H18-A3 (ISBN 1-56238-555-0). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA (2004).
7. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 413, 1800 (2018).
8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
9. Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2257, (2006).
10. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
11. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 05	Current version: 06
<i>Sections updated: WARNING AND NOTES.</i>	

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7. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 413, 1800 (2018).
8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
9. Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2257, (2006).
10. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
11. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Item No.	13	Item Name	ASAT II GEN	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
	Sample	18	5	µl
	Reagent 1	200	10	µl
	Reagent 2	50	10	µl
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Data Process				
Read		Start	End	
		Main	36	
		Sub		
Abs.Limit Low ~ High				
		-3	3	
Correction value				
Blank correction				
End Point Limit 2				
Linear Check (%) 80				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	13	Item Name	ASAT II GEN	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0	35	0	31
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	13	Item Name	ASAT II GEN	Optical			
Auto Rerun SW							
ON							
Auto Rerun Condition (Absorbance)							
		Lower	OFF				
		Higher	OFF				
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum			7			750	
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
100:DiI2							

▪ **Biolis 30i**

Item no	13	Item name	ASAT II GEN	Specimen	SERUM	OPTICAL												
Data information																		
UNITS	U/L		Aspiration volume															
DECIMALS	1		TYPE Double															
Analysis			<table border="1"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>18</td> <td>200</td> <td>50</td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>					SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	18	200	50	BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	18	200	50															
BOTTLE (ml)																		
METHOD	RATE method		FIRST DIL.															
Main Wave Length	340 nm																	
Sub Wave Length	450 nm																	
CORRELATION (Y= AX + B)																		
A =	1		Data processing read															
B =	0		<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> </tr> <tr> <td>MAIN</td> <td>38</td> <td>46</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table>					START	END	MAIN	38	46	SUB					
	START	END																
MAIN	38	46																
SUB																		
Blank value																		
• WATER ° REAGENT																		
Calibration																		
TYPE	Linear 2		ABS LIMIT															
STABILITY			-3 TO 3															
Collection value																		
END POINT			2.5															
LINEARITY CHECK (%)			80															
Prozone check																		
			START END LIMIT (%)															
FIRST																		
SECOND																		
			MINIMUM ABS.															
°HIGH			MEAN															
•LOW			VARIATE															

Item No	13	Item Name	ASAT II GEN	Specimen	SERUM	OPTICAL
Reference intervals						
MALE		FEMALE		Auto rerun		
LOW	HIGH	LOW	HIGH	°ON °OFF		
0	35	0	31	Auto rerun range (conc.)		
				Re	Value	Dil.
					9,0	770
Panic range						
MALE		FEMALE		Auto rerun condition (abs.)		
LOW	HIGH	LOW	HIGH	DIL.		
				LOWER	°ON	•OFF
				HIGH	°ON	•OFF
Decision limit						
MALE		FEMALE		Auto rerun condition (prozone)		
				°ON •OFF		
				SAMPLE VOL.		
				Dilution		
				•DIL 1 ° DIL 2		
Reaction check						
		°ON •OFF		CHECK		
				LOW		
				HIGH		
VL CHECK						
°ON		•OFF		VH CHECK		
				°ON •OFF		

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