PRESTIGE 24i LQ TOTAL PROTEIN

Nr kat. 4-236, 4-436 (PL)

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
Cat. No **4-236, 4-436** (EN)

INTENDED USE

TOTAL PROTEIN reagent is intended to determine quantitatively total protein level in plasma and serum. It is intended to monitoring and as an aid in the diagnosis of clinical conditions associated with abnormal protein concentration (hypo-, hyperproteinemia). TOTAL PROTEIN reagent is intended to use on automatic analysers Biolis 24i Premium and Biolis 30i.

It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY¹⁻⁵

Determinations of total protein level in serum and plasma are used to diagnosis and monitoring of clinical conditions associated with abnormal protein concentration. Increased total protein level can be observed in multiple myeloma, Waldenstrom macroglobulinemia, dehydration, hyperimmunoglobulinemia caused by polyclonal or monoclonal gammopathies, acute phase response. Decreased values can be observed in malnutrition, malabsorption, liver diseases (cirrhosis, fulminant failure), kidney diseases, protein-losing gastroenteropathies and acute burns.

METHOD PRINCIPLE⁴

The method is based on the biuret reaction. Protein forms the coloured complex with cupric ions in alkaline medium. The colour intensity is proportional to the protein concentration.

REAGENTS

Package

	Cat. No 4-236 (24-TRAY)	Cat. No 4-436 (36-TRAY)
1-REAGENT	6 x 60 ml	10 x 25 ml

Number of tests

Biolis 24i Premium	1320	900
Biolis 30i	1330	880

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT

sodium-potassium tartrate	30 mmol/l
potassium iodide	30 mmol/l
copper sulfate	12 mmol/l
sodium hydroxide	

REAGENT STABILITY

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagent is stable for 8 weeks (Biolis 24i Premium) or 11 weeks (Biolis 30i) on board the analyser at 2-10°C.

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains sodium hydroxide.

Danger

H314 Causes severe skin burns and eye damage.
H412 Harmful to aquatic life with long lasting effects.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

P273 Avoid release to the environment.

SPECIMEN^{2, 6, 7}

Serum or plasma.

Nonhemolyzed, nonlipaemic, fresh serum is recommended. Recommended anticoagulants: EDTA, heparine. Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

PROCEDURE

1-REAGENT is ready to use.

1-REAGENT put on basic position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: TOTAL PROTEIN – FERRUM, TOTAL PROTEIN - TG. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: TOTAL PROTEIN – FERRUM, TOTAL PROTEIN – TG. To avoid this effect follow

the recommendations contained in the instruction 51_03_24_009_BIOLIS 30i_CARRYOVER.

Application guide for Prestige 24i analyser is available on request. The application guide should be verified and validated by the user prior to testing patient samples.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analyser Biolis 24i Premium, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) is recommended.

For the calibration of automatic analyser Biolis 30i, the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 4.5 weeks (Biolis 24i Premium) or 2 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ^{4,8}

serum	g/dl	g/l
premature 1 day	3.6 – 6.0	36 – 60
children 1day – 2 weeks	5.3 – 8.3	53 – 83
15 days – 12 months	4.4 – 7.1	44 – 71
1 – 3 years	6.1 – 7.5	61 – 75
3 – 5 years	6.3 – 8.1	63 – 81
6 – 19 years	6.8 – 8.2	68 – 82
adults 20 – 29 years	6.5 – 8.3	65 – 83
30 – 79 years	6.5 – 7.8	65 – 78

Results obtained from plasma analysis might be slightly elevated due to fibrinogen presence:

Sample	Protein concentration increase from serum to plasma
Blood donors	+ 0.25 g/dl
Nonhospitalized patients	+ 0.36 g/dl
Hospitalized patients	+ 0.46 g/dl
Hospitalized patients with CRP<50 mg/dl	+ 0.66 g/dl

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical

history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**
0.03 g/dl (0.3 g/l)
- **LoD (Limit of Detection):**
0.05 g/dl (0.5 g/l)
- **LoQ (Limit of Quantitation):**
0.15 g/dl (1.5 g/l)
- **Linearity:**
up to 19 g/dl (190 g/l)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- **Measurement range:**
0.15 g/dl (1.5 g/l) - 19 g/dl (190 g/l)

Specificity / Interferences

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 950 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	4.54	0.03	0.59
level 2	6.88	0.06	0.87
Reproducibility (day to day) n = 80	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	6.7	0.10	1.5
level 2	4.3	0.07	1.7

Method comparison

A comparison between total protein values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 63 serum samples gave following results:

$$y = 0.9766x + 0.258 \text{ g/dl};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT ⁹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents:
1-REAGENT - 18 01 06*
- Empty packages:
1-REAGENT - 15 01 10*
- Wastewater from the analyzer:
1-REAGENT - 18 01 03*

INCIDENTS ¹⁰

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (websites address: incidents@cornay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

1. Szczeklik A., Gajewski P., ed. Interna Szczeklika. Mały podręcznik 2018/2019 10th ed., Kraków, Medycyna Praktyczna, (2018)
2. Alan H.B., Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B Saunders Company, 916 (2006)
3. Chowdhury, M., Akhter, N., Haque, M., Aziz, R., Nahar, N., Serum Total Protein and Albumin Levels in Different Grades of Protein Energy Malnutrition, *Journal of Bangladesh Society of Physiologist*, 3, 58-60 (2009).
4. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 396, 1780 (2018)
5. McPherson R.A., Pincus M.R., ed. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed., ELSEVIER, Philadelphia, (2011).
6. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
7. NCCLS - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition, NCCLS document H18-A3 (ISBN 1-56238-555-0), NCCLS, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, USA (2004)
8. Thomas L.:Clinical Laboratory Diagnostics. 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 644-647 (1998).
9. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 04	Current version: 05
Sections updated: WARNING AND NOTES ;	

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- McPherson R.A., Pincus M.R., ed. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed., ELSEVIER, Philadelphia, (2011).
- WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
- NCCLS - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition, NCCLS document H18-A3 (ISBN 1-56238-555-0), NCCLS, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, USA (2004)
- Thomas L.: Clinical Laboratory Diagnostics. 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 644-647 (1998).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Item No.	2	Item Name	TP	Optical
Aspiration				
Kind	Single			
Vol.				
	Kind	Vol.	Add	Units
	Sample	6	5	µl
	Reagent 1	250	10	µl
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Data Process				
Read	Start	End		
	Main	46	47	
	Sub	9	10	
Abs.Limit	Low	High		
	-0.1	2		
Correction value				
Blank correction				
End Point Limit 2				
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First			Low	
Second				

Item No.	2	Item Name	TP	Optical
Normal Range *				
	Male		Female	
	Low	High	Low	High
Serum	6.5	7.8	6.5	7.8
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	2	Item Name	TP	Optical			
Auto Rerun SW							
ON							
Auto Rerun Condition (Absorbance)							
	Lower	OFF					
	Higher	OFF					
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum			0.13			16.50	
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
100:Di2							

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• **Biolis 30i**

Item no	2	Item name	TP	Specimen	SERUM/ PLASMA	OPTICAL
Data information						
UNITS	g/dL		Aspiration volume			
DECIMALS	2		TYPE			
			Single			
			VOL. (µL)	SAMPLE	REAGENT 1	REAGENT 2
			BOTTLE (ml)	6	250	
			FIRST DIL.			
Analysis						
METHOD			END method			
Main Wave Length			546 nm			
Sub Wave Length			700 nm			
CORRELATION (Y= AX + B)						
A =			1		START	
B =			0		END	
Blank value						
			• WATER ° REAGENT			
Calibration						
TYPE			Linear 1			
STABILITY						
Data processing read						
MAIN			46		47	
SUB						
Abs LIMIT						
-0.1			TO		2	
Collection value						
END POINT			2.5			
LINEARITY CHECK (%)			0			
Prozone check						
FIRST			START		LIMIT (%)	
SECOND						
			MINIMUM ABS.			
°HIGH			MEAN			
•LOW			VARIATE			

Item No	2	Item Name	TP	Specimen	SERUM/PLASMA	OPTICAL
Reference intervals *						
MALE		FEMALE				
LOW	HIGH	LOW	HIGH			
6.5	7.8	6.5	7.8			
Panic range						
MALE		FEMALE				
LOW	HIGH	LOW	HIGH			
Decision limit						
MALE		FEMALE				
Reaction check						
			°ON		•OFF	
CHECK						
LOW						
HIGH						
VL CHECK						
°ON			•OFF		VH CHECK	
			°ON		•OFF	
Auto rerun						
			•ON		°OFF	
Auto rerun range (conc.)						
Re	Value	Dil.	Re	Value	Dil.	
	0.15			19		
Auto rerun condition (abs.)						
LOWER			°ON		•OFF	
HIGH			°ON		•OFF	
Auto rerun condition (prozone)						
°ON			•OFF			
SAMPLE VOL.						
Dilution						
•DIL 1			° DIL 2			

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