PRESTIGE 24i CREA ENZYMATYC

Nr kat. 4-237, 4-437

(PL)

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LITERATURA

- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhard V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

Data wydania: 04.2021.

PRESTIGE 24i CREA ENZYMATIC

Cat. No 4-237 and 4-437

(EN)

2-Reagent

Good's buffer	≤ 5%
creatininase	≤ 1%
peroxidase	≤ 5%
4-ammoantipyrine (4-AA)	≤ 0.01%
stabilizers and preservatives	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum and urine.

Urine preparation: before analysis urine sample should be diluted with 0,9% NaCl 2-10 times. Multiply the result by dilution factor.

Urine can be stored up to 1 day at 20-25°C, 4 days at 2-8°C. For longer storage samples should be frozen at -20°C. Serum can be stored up to 1 day at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CREA ENZYMATIC – URINE PROTEINS, GLUCOSE - CREA ENZYMATIC, URINE PROTEINS – CREA ENZYMATIC. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CREA ENZYMATIC – URINE PROTEINS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ^{2,3}

serum / plasma	mg/dl	µmol/l
newborns	0.3 – 1.0	26.5 – 88.4
infants	0.2 – 0.4	17.7 – 35.4
children	0.2 – 0.8	17.7 – 70.7
female	0.5 – 1.0	44.2 – 88.4
male	0.7 – 1.2	61.9 – 106.1
urine (morning)	mg/dl	mmol/l
female	29 – 226	2.56 – 20.0
male	40 – 278	3.54 – 24.6

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers: Prestige 24i and Biolis 24i Premium, the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analysers Biolis 30i, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 2 weeks (Prestige 24i) or every 4 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.01 mg/dl (0.88 µmol/l)
- LoD (Limit of Detection):**
0.03 mg/dl (2.65 µmol/l)
- LoQ (Limit of Quantitation):**
0.1 mg/dl (8.84 µmol/l)
- Linearity:**
up to 27 mg/dl (2386.8 µmol/l)

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbic acid up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and creatine up to 20 mg/dl both at serum and urine do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.6	0.01	0.42
level 2	5.2	0.03	0.52
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.7	0.04	2.5
level 2	5.2	0.09	1.8

Method comparison

A comparison between creatinine values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 61 serum samples gave following results:

$$y = 1.051x - 0.0236 \text{ mg/dl;} \\ R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between creatinine values determined at **Biolis 30i** (y) and at **BS-800** (x) using 59 urine samples gave following results:

$$y = 1.0331x + 0.0699 \text{ mg/dl;} \\ R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liapps N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

Date of issue: 04.2021.

Components and concentrations

1-Reagent

Good's buffer	≤ 5%
creatininase	≤ 5%
N-ethyl-N-(3-sulfopropyl)-3-methylaniline (ESPMT)	≤ 5%
sarcosine oxidase	≤ 0.01%
ascorbate oxidase	≤ 1%
detergents, stabilizers and preservatives	

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1. Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
2. Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
3. Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
4. Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
5. Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
6. Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

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■ Prestige 24i Premium, Biolis 24i Premium

Item No.	69	Item Name	Crea Enz	Optical
Data information				
Units	mg/dl			
Decimals	2			
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	END method			
Main Wave Length	546 nm			
Sub Wave Length	800 nm			
Method				
Correlation				
Slope	Intercept			
Y=	1	X+	0	
Item No.	69	Item Name	Crea Enz	Optical
Aspiration				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	8	5	µl	
Reagent 1	180	10	µl	
Reagent 2	60	10	µl	
Data Process				
Read	Main	Start	End	
	Sub	28	30	
Abs.Limit	Low	High		
	-0.1	2		
Correction value				
Blank correction	0.7581			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low
Item No.	69	Item Name	Crea Enz	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.7	1.2	0.5	1.0
Urine	40	278	29	226
Plasma	0.7	1.2	0.5	1.0
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Item No.	69	Item Name	Crea Enz	Optical
Auto Rerun SW				
ON				
Auto Rerun Condition (Absorbance)				
	Lower	Higher	OFF	
			OFF	
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			0.08	24
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Prozone)				
	OFF			
Dilution				
	100:Dil2			

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