PRESTIGE 24i LQ BIL TOTAL

Nr kat. 4-245, 4-445 (PL)

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LITERATURA

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2):116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood E.R., WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

Data wydania: 04. 2021.

PRESTIGE 24i LQ BIL TOTAL

Cat. No **4-245, 4-445** (EN)

INTENDED USE

Diagnostic kit for determination of total bilirubin concentration intended to use in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Bilirubin is a yellow pigment – product of heme degradation. For clinical purposes, bilirubin is expressed as two fractions: conjugated and unconjugated. In hepatocytes bilirubin is enzymatically conjugated with glucuronic acid residues. This form is called direct or conjugated. Bilirubin without glucuronic acid modification is bound to albumin and is termed unconjugated or indirect. Indirect bilirubin is calculated as the difference between total and direct bilirubin.

Serum bilirubin measurement is widely used as a screening test for liver functions. Hiperbilirubinemia is usually the result of jaundice (mechanical, hemolytic), Dubin-Jonson syndrome, Gilbert's syndrome, Crigler-Najjar syndrome, bile ducts disease.

METHOD PRINCIPLE

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and vanadate in an acidic solution, total bilirubin (both conjugated – direct, and unconjugated bilirubin) is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the total bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

REAGENTS

Package

	Cat. No 4-245 (24-TRAY)	Cat. No 4-445 (36-TRAY)
1-Reagent	8 x 41.5 ml	11 x 23.0 ml
2-Reagent	8 x 12.5 ml	11 x 7.0 ml

The reagents when stored at 10-25°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 7 weeks, Biolis 24i Premium - 7 weeks.

Concentrations in the test

1-Reagent

citrate buffer (pH 2.8) 90 mmol/l
detergent

2-Reagent

phosphate buffer (pH 7.0) 4.6 mmol/l
sodium metavanadate 3.0 mmol/l

Warnings and notes

- Do not freeze reagents.
- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not interchange caps.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously

with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum free from hemolysis.

Serum should be separated from red blood cells as soon as possible after blood collection. Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended.

It is recommended to follow CLSI procedures regarding specimen collecting and handling.

Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark at 2-8°C, at the most 3 days.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser there is a probability of **cross-contamination** affecting the tests results: CHOLINESTERASE – BIL TOTAL, CREATININE BIL TOTAL. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

REFERENCE VALUES⁴

Age	mg/dl	μmol/l
0 - 1 days	< 8	< 137
1 - 2 days	< 12	< 205
3 - 5 days	< 16	< 274
5 days - 60 years	0.3 - 1.2	5 - 21
60 - 90 years	0.2 - 1.1	3 - 19
> 90 years	0.2 - 0.9	3 - 15

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems are recommended CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177).

The calibration curve should be prepared every 7 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.04 mg/dl (0.68 μmol/l)
- LoD (Limit of Detection):**
0.05 mg/dl (0.86 μmol/l)
- LoQ (Limit of Quantitation):**
0.15 mg/dl (2.57 μmol/l)
- Linearity:**
up to 88 mg/dl (1504.8 μmol/l)
- Specificity / Interferences**
Haemoglobin up to 0.25 g/dl, ascorbic acid up to 500 mg/l and triglycerides up to 250 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	0.96	0.01	0.82
level 2	4.30	0.05	1.22
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	0.9	0.03	3.3
level 2	4.5	0.07	1.7

Method comparison

A comparison between total bilirubin values determined at **Biolis 30i** (y) and at **SIEMENS ADVIA 1800** (x) using 46 serum samples gave following results:

$$y = 1.0083x - 0.0471 \text{ mg/dl;}$$

$$R = 0.992$$

(R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2);116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis C.A. and Ashwood E.R., WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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Item name	3	Bil T		
Auto Rerun SW				
ON				
Auto Rerun Range (Result)				
	ON	ON		
	Lower	Higher		
Serum	0.1	70		
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
Absorbance Range				
	Lower	Higher	OFF	
			OFF	
Auto Rerun Condition (Prozone)				
Prozone Range				
OFF				

▪ Prestige 24i Premium, Biolis 24i Premium

Item No.	3	Item Name	Bil T	Optical	
Data information					
Units	mg/dl				
Decimals	2				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Analysis					
Type	END method				
Main Wave Length	450 nm				
Sub Wave Length	546 nm				
Method	Vanadate				
Correlation					
	Slope	Intercept			
Y=	1	X+	0		

Item No.	3	Item Name	Bil T	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		12	5	µl
Reagent 1		280	10	µl
Reagent 2		70	10	µl
Data Process				
Read	Start	End		
Main	51	52		
Sub	29	31		
Abs.Limit				
Low	-3			High
				3
Correction value				
Blank correction				0.8066
End Point Limit				2
Linear Check (%)				90
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				

Item No.	3	Item Name	Bil T	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.3	1.2	0.3	1.2
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	3	Item Name	Bil T	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			0.2	59
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
Lower				
Higher				
OFF				
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:DiI2				

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