PRESTIGE 24i LQ CALCIUM

Nr kat. 4-251, 4-451 (PL)

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Cat. No 4-251, 4-451 (EN)

INTENDED USE

Diagnostic kit for determination of calcium concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Calcium and phosphorus as a hydroxyapatite constitute mineral portion of bone. Calcium occurs also as divalent cations (free or bound with negatively charged proteins) which participate in blood coagulation, neuromuscular excitability, skeletal and cardiac muscle contractility and in multiple cellular functions. Calcium flux in organism is controlled by action of parathyroid hormone (PTH), vitamin D and calcitonin. Calcium serum level abnormalities are caused usually by parathyroid or thyroid disease, disorders of vitamin D metabolism or acute pancreatitis.

METHOD PRINCIPLE

Calcium ions form a violet complex with o-cresolphthalein complexone in alkaline solution. The intensity of violet colour of this complex measured at 570-580 nm is proportional to the calcium concentration in the sample.

REAGENTS

Package

| | Cat. No 4-251 (24-TRAY) | Cat. No 4-451 (36-TRAY) |
|-----------|----------------------------|----------------------------|
| 1-Reagent | 4 x 40 ml | 8 x 23 ml |
| 2-Reagent | 4 x 12.5 ml | 8 x 7.5 ml |

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 4 weeks, Biolis 24i Premium - 3 weeks.

Concentrations in the test

| | |
|------------------------------|-------------|
| o-cresolphthalein complexone | 0.06 mmol/l |
| 8-quinolinol | 8.6 mmol/l |
| hydrochloric acid | 30 mmol/l |
| ethanolamine | 377 mmol/l |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2 M HCl solution and then thoroughly rinsed with distilled water.
- 1-Reagent and 2-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent contains hydrochloric acid.

Danger

H314 Causes severe skin burns and eye damage.
 H318 Causes serious eye damage.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water (or shower).

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

Ingredients:

2-Reagent contains 2-aminoethanol.

Danger

H314 Causes severe skin burns and eye damage.
 H318 Causes serious eye damage.

H335 May cause respiratory irritation.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water (or shower).

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

SPECIMEN

Serum, heparinized plasma free from hemolysis, 24-hours urine. Recommended anticoagulants: heparine lithium, sodium or ammonium salt.

Urine preparation: to prevent calcium salt precipitation in urine, specimens should be collected in HCl, 20-30 ml of 6 M for 24-h specimen. Any specimens collected without acid should be acidified using 20-30 ml of 6 M HCl, well mixed and allowed to stand for 1 h before aliquotting. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor.

Serum and plasma can be stored up to 8 hours at 15-25°C or up to 1 day at 2-8°C. Samples frozen at -20°C can be stored up to 1 year.

24-hours urine samples should be kept at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionised water is recommended.

Actions required:

▪ **Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LIPASE II GEN - CALCIUM, RF - CALCIUM. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

▪ **Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: RF - CALCIUM. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁸

| serum, plasma | mg/dl | mmol/l |
|----------------|------------|-------------|
| premature | 6.2 – 11.0 | 1.55 – 2.75 |
| adults | | |
| 18 – 60 yr | 8.6 – 10.0 | 2.15 – 2.50 |
| 60 – 90 yr | 8.8 – 10.2 | 2.20 – 2.55 |
| > 90 yr | 8.2 – 9.6 | 2.05 – 2.40 |
| 24-hours urine | mg/24h | mmol/24h |
| | 100 – 300 | 2.5 – 7.5 |

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems Prestige 24i and Biolis 24i Premium the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analyser Biolis 30i the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) is recommended.

The calibration curve should be prepared every 3 days (Prestige 24i) or every 2 days (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):** 0.07 mg/dl (0.018 mmol/l)
- **LoD (Limit of Detection):** 0.15 mg/dl (0.038 mmol/l)
- **LoQ (Limit of Quantitation):** 0.22 mg/dl (0.055 mmol/l)
- **Linearity:** up to 17 mg/dl (4.25 mmol/l)

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

| Repeatability (run to run) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|--|-----------------|---------------|-----------|
| level 1 | 8.82 | 0.05 | 0.56 |
| level 2 | 11.69 | 0.07 | 0.62 |
| Reproducibility (day to day) n = 80 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
| level 1 | 8.7 | 0.29 | 3.4 |
| level 2 | 11.6 | 0.36 | 3.1 |

Method comparison

A comparison between calcium values determined at **Biolis 30i** (y) and at **BS-400** (x) using 60 serum samples gave following results:

$$y = 1.0778x - 0.4895 \text{ mg/dl};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between calcium values determined at **Biolis 30i** (y) and at **BS-400** (x) using 60 urine samples gave following results:

$$y = 0.9797x + 0.1159 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | |
|--|-------|--------|
| Item name | 23 | CALC |
| Auto Rerun SW | | |
| ON | | |
| Auto Rerun Condition (Absorbance) | | |
| Absorbance Range | Lower | Higher |
| | OFF | OFF |
| Auto Rerun Range (Result) | | |
| | ON | ON |
| | Lower | Higher |
| Serum | 1.11 | 15 |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| Other | | |
| Prozone Range | | |
| | OFF | |

• **Prestige 24i Premium, Biolis 24i Premium**

| | | | | |
|-------------------------|------------|-----------|------|---------|
| Item No. | 23 | Item Name | CALC | Optical |
| Data information | | | | |
| Units | mg/dl | | | |
| Decimals | 2 | | | |
| Calibration | | | | |
| Type | Linear1 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 | * | #2 |
| #3 | | #4 | | #5 |
| #6 | | | | |
| Analysis | | | | |
| Type | END method | | | |
| Main Wave Length | 570 nm | | | |
| Sub Wave Length | 700 nm | | | |
| Method | o-cresolo. | | | |
| Correlation | | | | |
| Y= | Slope | Intercept | | |
| | 1 | 0 | | |

| | | | | |
|-------------------------|--------|-----------|------------------|---------|
| Item No. | 23 | Item Name | CALC | Optical |
| Aspiration | | | | |
| Kind | Double | | | |
| Vol. | | | | |
| | Kind | Vol. | Add | Units |
| Sample | 5 | 5 | 5 | µl |
| Reagent 1 | 160 | 10 | 10 | µl |
| Reagent 2 | 40 | 10 | 10 | µl |
| Data Process | | | | |
| Read | Start | End | Abs.Limit | |
| Main | 35 | 36 | Low | -0.1 |
| Sub | 30 | 31 | High | 2 |
| Correction value | | | | |
| Blank correction | | | | |
| End Point Limit | | | | |
| 2 | | | | |
| Linear Check (%) | | | | |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | Low | |
| Blank value | | | | |
| Water Blank | | | | |
| Reaction Monitor | | | | |
| 0 Level Point | | | | |
| 1 | | | | |
| Span | | | | |
| 3 | | | | |
| Third mixing | | | | |
| OFF | | | | |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 23 | Item Name | CALC | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | 8.6 | 10.0 | 8.6 | 10.0 |
| Urine | | | | |
| Plasma | 8.6 | 10.0 | 8.6 | 10.0 |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | | |
|--|-----------|-----------|------|------------|-----|
| Item No. | 23 | Item Name | CALC | Optical | |
| Auto Rerun SW | | | | | |
| ON | | | | | |
| Auto Rerun Condition (Absorbance) | | | | | |
| | Lower | Higher | | | |
| | OFF | OFF | | | |
| Auto Rerun Range (Conc.) | | | | | |
| | First Dil | Low Value | Dil | High Value | Dil |
| Serum | | 0.27 | | 15 | |
| Urine | | | | | |
| Plasma | | | | | |
| CSF | | | | | |
| Dialysis | | | | | |
| Other | | | | | |
| Auto Rerun Condition (Prozone) | | | | | |
| OFF | | | | | |
| Dilution | | | | | |
| 100:Dil2 | | | | | |

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