PRESTIGE 24i LQ FERRUM

Cat. No 4-258, 4-458

(EN)

Danger

H314 Causes severe skin burns and eye damage.
EUH208 Contains thiourea. May produce an allergic reaction (1 Reagent).
EUH208 Contains 1-[1,3-Bis (hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-bis

(hydroxymethyl) urea. May produce an allergic reaction (2-Reagent).

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

SPECIMEN¹³

Serum free from hemolysis, collected in plastic tubes.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum can be stored up to 7 days at 15-25°C or up to 3 weeks at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- **Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: TOTAL PROTEIN – FERRUM. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

- **Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: FERRUM – dTIBC, FERRUM – URINE PROTEINS, FERRUM – FERRITIN, TOTAL PROTEIN – FERRUM. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES^{8, 10}

serum	µg/dl	µmol/l
newborns	100 – 250	17.9 – 44.8
infants	40 – 100	7.2 – 17.9
children	50 – 120	9.0 – 21.5
adult females	50 – 170	9.0 – 30.4
adult males	65 – 175	11.6 – 31.3

Samples should be taken in the morning from patients in a fasting state, since iron values decrease by 30% during the course of the day.

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems Prestige 24i and Biolis 24i Premium the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number.

For the calibration of the automatic analyser Biolis 30i the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 11 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):** 0.7 µg/dl (0.13 µmol/l)
- **LoD (Limit of Detection):** 2.3 µg/dl (0.41 µmol/l)
- **LoQ (Limit of Quantitation):** 4 µg/dl (0.72 µmol/l)
- **Linearity:** up to 1120 µg/dl (200.48 µmol/l)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and copper up to 500 µg/dl do not interfere with the test. Haemoglobin interferes even in small amount with the determination.

Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[µg/dl]	[µg/dl]	[%]
level 1	65.4	0.70	1.06
level 2	227.2	1.48	0.65
Reproducibility (day to day) n = 80	Mean	SD	CV
	[µg/dl]	[µg/dl]	[%]
level 1	234.6	5.92	2.5
level 2	66.9	2.16	3.2

Method comparison

A comparison between ferrum values determined at Biolis 30i (y) and at ADVIA SIEMENS 1800 (x) using 46 serum samples gave following results:

$$y = 1.0077 x + 4.4208 \text{ µg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Stookey L.L.: Anal. Chem. 42/7, 779-781 (1970).
2. Williams H.L., Johnson D.J., Haut M.J.: Clin. Chem. 23/2, 237-240 (1977).
3. Duffy J.R., Gaudin J.: Clin. Biochem. 10/3, 122-123 (1977).
4. Ceriotti F., Ceriotti G: Clin. Chem. 26/2, 327-331 (1980).
5. Valcour A., Krzymowski G., Onorowski M., Bowers G.N. Jr., McComb R.B.: Clin Chem. 36/10, 1789-1792 (1990).
6. Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, 2062 (1994).
7. Tietz N.W., Textbook of Clinical Chemistry, Philadelphia, PA: WB Saunders, 3:24, (1990).
8. Kaplan L.A., Pesce A.J., ed. Clinical chemistry, theory, analysis and correlation, 3rd ed. St Louis, MO: Mosby, 714 (1996).
9. Dembińska-Kieć A., Naskalski J.W.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Volumes, 24-25, (1998).
10. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 634, (2006).
11. Tietz NW, Rinker AD, Morrison SR. Clin Chem. 40(4):546-51 (1994).
12. Br J Haematol. 75(4):615-6 (1990).
13. Ehret W., Heil W., Schmitt Y., Töpfer G., Wisser H., Zawta B., et al. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 Rev. 2, p. 36

Date of issue: 04.2021.

The reagents are stable up to expiry date printed on the package, if stored at 2-8°C. Stability on board of the analyser at 2-10°C: Prestige 24i – 11 weeks, Biolis 24i Premium – 11 weeks.

Concentrations in the reagent

1-Reagent

citric acid (pH 1.9) ≤ 240 mmol/l
thiourea ≤ 108 mmol/l
detergent ≤ 7 %

2-Reagent

sodium ascorbate ≤ 150 mmol/l
3-(2-pyridyl)-5,6-bis(2-[5-furyl sulfonic acid])-
1,2,4-triazine sodium salt (ferrozine) ≤ 6 mmol/l
preservative
stabilizer

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent contains polyethylene glycol monoalkyl ether.

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Item name	24	FERRUM II GEN
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	7.0	1000
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	Higher
		OFF
		OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

• Prestige 24i Premium, Biolis 24i Premium

Item No.	24	Item Name	FERRUM II GEN	Optical
Data information				
Units	µg/dl			
Decimals	1			
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	END method			
Main Wave Length	570 nm			
Sub Wave Length	700 nm			
Method	Ferrozine			
Correlation				
Slope	Y=		Intercept	X+
	1			0

Item No.	24	Item Name	FERRUM II GEN	Optical
Aspiration				
Kind				
Vol.				
	Kind	Vol.	Add	Units
Sample	20	5		µl
Reagent 1	250	10		µl
Reagent 1	50	10		µl
Data Process				
Read	Main	Start	End	
	Sub	28	30	
Abs.Limit Low ~ High				
	-3			3
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	24	Item Name	FERRUM II GEN	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	65	175	50	170
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	24	Item Name	FERRUM II GEN	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low Value	Dil	High Value
Serum		3.6		1000
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	Higher		
			OFF	
			OFF	
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:Di2				

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