PRESTIGE 24i LACTATE

Nr kat. **4-386, 4-466** (PL)

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LITERATURA

- Tietz Textbook of Clinical Chemistry (Edited by Burtis CA and Ashwood ER Eds): Third Edition WB Saunders Company 787-8, (1999).
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 650-652, (2006).

Data wydania: 04.2021.

PRESTIGE 24i LACTATE

Cat. No **4-386, 4-466**

(EN)

INTENDED USE

Diagnostic kit for determination of lactate concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lactate is produced in Cori cycle, by anaerobic conversion of glucose, mainly in skeletal muscle. Its determination, frequently done together with pyruvate, is useful in discovering lactic acidosis due to i.a. reduced tissue oxygenation, enzymatic deficiencies, diabetes mellitus, liver and kidneys diseases.

METHOD PRINCIPLE

Lactate is oxidized by lactate oxidase to pyruvate and hydrogen peroxide, which, in presence of peroxidase (POD), reacts with 4-aminoantipyrine and phenol forming a compound, which colour intensity is proportional to the concentration of lactate in the examined sample.

REAGENTS

Package

	Cat. No 4-386 (24-TRAY)	Cat. No 4-466 (36-TRAY)
1-Reagent	2 x 19,5 ml	2 x 17,5 ml

Unopened reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. Stability on board of the analyser at 2-10°C: Prestige 24i - 12 weeks, Biolis 24i Premium - 12 weeks.

Concentrations in the test

Tris buffer (pH 7.5)	≥ 50 mmol/l
lactate oxidase	≥ 0.2 kU/l
peroxidase	≥ 2 kU/l
4-aminoantipyrine	≥ 0.4 mmol/l

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent!
- Do not use after expiry date.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lactate concentration rapidly increases during physical activities. Normal levels are reached again after usually 30 minutes but it may vary according to individuals.
- Draw blood with lowest venous stasis as possible (max. 30 seconds) from fasting and completely resting patient and avoid using a tourniquet.

SPECIMEN

Plasma. Avoid haemolysis. Collect samples in tubes containing sodium fluoride and potassium oxalate. Keep samples on ice. Centrifuge within 15 minutes after collection and separate from cells. Analyze promptly. Note whether sample is venous or arterial. It is recommended to follow NCCLS procedures regarding specimen collecting and handling. Lactate in plasma is stable up to 8 hours at room temperature or up to 14 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.
1-Reagent put on basic position in reagent tray.
For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LACTATE – ALAT (II GEN). To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LACTATE – CRP ULTRA. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ²

plasma (venous)	4.5 – 19.8 mg/dl	0.5 – 2.2 mmol/l
plasma (arterial)	4.5 – 14.4 mg/dl	0.5 – 1.6 mmol/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.07 mg/dl (0.0078 mmol/l)

- LoD (Limit of Detection):**
0.17 mg/dl (0.0189 mmol/l)

- LoQ (Limit of Quantitation):**
1.40 mg/dl (0.1554 mmol/l)

- Linearity:**
up to 139 mg/dl (15.43 mmol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

- In plasma samples containing approximately 12 mg lactate/dl, there is no interference up to: 0.23 g/dl haemoglobin, 8 mg/dl bilirubin, 337 mg/dl triglycerides, 15.5 mg/l ascorbic acid.
- In plasma samples containing approximately 40 mg lactate/dl, there is no interference up to: 1.25 g/dl haemoglobin, 10 mg/dl bilirubin, 1000 mg/dl triglycerides, 62 mg/l ascorbic acid.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	8.8	0.08	0.92
level 2	42.4	0.30	0.70
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	9.1	0.37	4.1
level 2	40.4	0.93	2.3

Method comparison

A comparison between lactate values determined at **Biolis 30i** (y) and at **BS-400** (x) using 58 plasma samples gave following results:

$$y = 1.035 x - 0.3969 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tietz Textbook of Clinical Chemistry (Edited by Burtis CA and Ashwood ER Eds): Third Edition WB Saunders Company 787-8, (1999).
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 650-652, (2006).

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Item name	66	LACT
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum		
Urine		
Plasma	2.40	90
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Prozone Range		
		OFF

▪ Prestige 24i Premium, Biolis 24i Premium

Item No.	66	Item Name	LACT	Optical	
Data information					
Units	mg/dl				
Decimals	2				
Analysis					
Type	END method				
Main Wave Length	546nm				
Sub Wave Length	700nm				
Method	Colorimetric				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Correlation					
Y=	Slope	X+	Intercept		
	1		0		

Item No.	66	Item Name	LACT	Optical
Aspiration				
Kind	Single			
Data Process				
Read	Main	Start	End	
	Sub	9	10	
Abs.Limit				
Low	-3			High
				3
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	66	Item Name	LACT	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma	4.5	19.8	4.5	19.8
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	66	Item Name	LACT	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma		1.4	90	
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
	OFF			
Dilution				
	100:DiI2			

PRESTIGE 24i LACTATE

• **Biolis 30i**

Item no	66	Item name	LACT	Specimen	PLASMA	OPTICAL
Data information			Aspiration volume			
UNITS	mg/dL			TYPE	Single	
DECIMALS	2				SAMPLE	REAGENT 1
Analysis						
METHOD	END method			VOL. (µL)	3	300
Main Wave Length	546 nm			BOTTLE (ml)		
Sub Wave Length	700 nm			FIRST DIL.		
CORRELATION (Y= AX + B)			Data processing read			
A =	1			START	END	
B =	0			MAIN	52	53
Blank value						
• WATER ° REAGENT						
Calibration						
TYPE	Linear 2			ABS LIMIT	-3	3
STABILITY				TO		
			Collection value			
			END POINT 2.5			
			LINEARITY CHECK (%) 0			
			Prozone check			
			START END LIMIT (%)			
			FIRST			
			SECOND			
			MINIMUM ABS.			
			°HIGH MEAN			
			•LOW VARIATE			

Item No	66	Item Name	LACT	Specimen	PLASMA	OPTICAL
Reference intervals			Auto rerun			
MALE		FEMALE		•ON		°OFF
LOW	HIGH	LOW	HIGH			
4.5	19.8	4.5	19.8			
Panic range			Auto rerun range (conc.)			
MALE		FEMALE		Re	Value	Dil.
LOW	HIGH	LOW	HIGH		1.4	139
Decision limit			Auto rerun condition (abs.)			
MALE		FEMALE				DIL.
LOW	HIGH	LOW	HIGH	LOWER	°ON	•OFF
			HIGH °ON •OFF			
			Auto rerun condition (prozone)			
			°ON •OFF			
			SAMPLE VOL.			
			Dilution			
			•DIL 1 ° DIL 2			
Reaction check						
°ON		•OFF				
CHECK						
LOW						
HIGH						
VL CHECK		VH CHECK				
°ON	•OFF	°ON	•OFF			

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