## PRESTIGE 24i LQ AMYLASE EPS

Nr kat. 4-176, 4-476

(PL)

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#### LITERATURA

- Pesce, A.J., Kaplan, L.A.: "Methods in Clinical Chemistry", Mosby Ed. (1987).
- Burtis C.A., Ashwood E.R.: "Tietz Textbook of Clinical Chemistry", W.B. Saunders Company Ed. (3rd edition, 1999).
- Lorentz K. Approved recommendation on IFCC methods for the measurement of catalytic concentration of enzymes. Part 9. IFCC method for alpha-amylase (1,4-alpha-D-glucan 4-glycanohydrolase, EC 3.2.1.1). International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Committee on Enzymes. Clin Chem Lab Med. 1998 Mar;36(3):185-203.
- Junge W, Wortmann W, Wilke B, Waldenström J, Kurlle-Weittenhiller A, Finke J, Klein G. Development and evaluation of assays for the determination of total and pancreatic amylase at 37 degrees C according to the principle recommended by the IFCC. Clin Biochem. 2001 May;34:607-15.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 100-104, (2006).
- Hohenwallner W, Hagele EO, Scholer A et al. Ber Oster Ges Klin Chem 1983;6:101-112.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 46-8 (1995).

**Data wydania:** 04. 2021

## PRESTIGE 24i LQ AMYLASE EPS

Cat. No. 4-176, 4-476

(EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -amylase activity used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION



$\alpha$ -Amylases are hydrolytic enzymes which hydrolyze 1,4- $\alpha$ -glucosidic bond in starch and other similar polysaccharides to maltose and other oligosaccharides. Several types of amylases can be distinguished, depending on the organ they are originating from.  $\alpha$ -amylase is the most commonly measured in the diagnosis of acute pancreatitis, when its activity in serum is very high. Elevation of  $\alpha$ -amylase activity in serum is also accompanied by increased excretion of enzyme in urine which can last longer than in the blood. Because of that activity in  $\alpha$ -amylase in urine is used as a indicator of acute pancreatitis. Hyperamylasemia occurs also in chronic pancreatitis, failures of kidneys, lungs, diseases of the salivary glands, cerebral traumas, surgical interventions and macroamylasemia. To confirm pancreatic specificity it is recommended to determine also other pancreas specific enzyme like lipase.

### METHOD PRINCIPLE

Enzymatic colorimetric method, with EPS substrate, in accordance to recommendations of IFCC – International Federation of Clinical Chemistry and Laboratory Medicine (modified IFCC method).

$\alpha$ -Amylase catalyzes hydrolysis of substrate 4,6-ethylidene-(G7)-p-nitrophenyl-(G1)- $\alpha$ -D-maltoheptaoside (EPS, Ethylidene Protected Substrate). Ethylidene group prevents the substrate from breaking down because of exo-enzymes activity, therefore in absence of  $\alpha$ -amylase no increase of absorbance is observed.

$\alpha$ -Amylase hydrolyses the substrate into smaller fragments which are acted upon by  $\alpha$ -glucosidase, causing the ultimate release of chromophore p-nitrophenol (pNP) and glucose.

Increase of absorbance related to formation of p-nitrophenol is proportional to the  $\alpha$ -amylase activity in sample and is measured spectrophotometrically at 405 nm wavelength.

### REAGENTS

#### Package

	Cat. No 4-176 (24-TRAY)	Cat. No 4-476 (36-TRAY)
1-Reagent	2 x 24 ml	2 x 23 ml
2-Reagent	2 x 8 ml	2 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks.

### Concentrations in the test

HEPES buffer, pH 7.2	52.5 mmol/l
sodium chloride	87 mmol/l
magnesium chloride	12.6 mmol/l
calcium chloride	0.075 mmol/l
$\alpha$ -glucosidase	$\geq 4\text{ kU/l}$
4,6-ethylidene G7pNP (EPS)	$> 4\text{ mmol/l}$
stabilizers and preservatives	

### Warnings and notes

- Prevent the reagents from microbiological contamination and from saliva and sweat  $\alpha$ -amylase.
- Protect from direct sunlight.
- The reagents must be clear, do not use if turbid.
- A slight yellow colour of 2-Reagent does not influence the reagent performance.
- 1-Reagent and 2-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H - isothiazol-3-one, mixture (3:1).

### Warning

H317 - May cause an allergic skin reaction.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P333 + P313 If skin irritation or rash occurs: Get medical advice.  
P363 Wash contaminated clothing before reuse.

### SPECIMEN

Serum or plasma collected on heparin, free from hemolysis, urine.

Do not use anticoagulants: EDTA, citrates and oxalates as they inhibit amylase activity.

Serum / plasma can be stored for 7 days at 15-25°C or for one month at 2-8°C.<sup>7</sup>

Urine can be stored for 2 days at 15-25°C or for 10 days at 2-8°C.<sup>6</sup> Amylase is very unstable in acid urine. Adjust pH to approximately 7.0 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

### REFERENCE VALUES<sup>5</sup>

serum / plasma	28 – 100 U/l	0.47 – 1.7 $\mu\text{kat/l}$
urine	$\leq 460\text{ U/l}$	$\leq 7.7\text{ }\mu\text{kat/l}$

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum, CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers: Prestige 24i, Biolis 24i Premium the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of Biolis 30i automatic analyser the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 0.1 U/l (0.00167  $\mu\text{kat/l}$ )
- LoD (Limit of Detection):** 0.6 U/l (0.01  $\mu\text{kat/l}$ )
- LoQ (Limit of Quantitation):** 3.0 U/l (0.05  $\mu\text{kat/l}$ )
- Linearity:** up to 2300 U/l (38.3  $\mu\text{kat/l}$ )

For higher activity, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.156 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l, triglycerides up to 1250 mg/dl and glucose up to 2000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	75.1	0.68	0.90
level 2	411.9	2.20	0.53
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	76.4	0.87	1.1
level 2	409.3	4.67	1.1

### Method comparison

A comparison between amylase values determined at **Biolis 30i** (y) and at **SIEMENS ADVIA 1800** (x) using 56 serum samples gave following results:

$y = 0.9946x + 0.4623\text{ U/l}$ ;

$R = 1.000$  (R – correlation coefficient)

A comparison between amylase values determined at **Biolis 30i** (y) and at **BS-800** (x) using 62 urine samples gave following results:

$y = 1.035x - 4.2576\text{ U/l}$ ;

$R = 1.000$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Pesce, A.J., Kaplan, L.A.: “Methods in Clinical Chemistry”, Mosby Ed. (1987).
- Burtis C.A., Ashwood E.R.: “Tietz Textbook of Clinical Chemistry”, W.B. Saunders Company Ed. (3<sup>rd</sup> edition, 1999).
- Lorentz K. Approved recommendation on IFCC methods for the measurement of catalytic concentration of enzymes. Part 9. IFCC method for alpha-amylase (1,4-alpha-D-glucan 4-gluconohydrolase, EC 3.2.1.1). International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Committee on Enzymes. Clin Chem Lab Med. 1998 Mar;36(3):185-203.
- Junge W, Wortmann W, Wilke B, Waldenström J, Kurlle-Weittenhiller A, Finke J, Klein G. Development and evaluation of assays for the determination of total and pancreatic amylase at 37 degrees C according to the principle recommended by the IFCC. Clin Biochem. 2001 May;34:607-15.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 100-104, (2006).
- Hohenwallner W, Hagele EO, Scholer A et al. Ber Oster Ges Klin Chem 1983;6:101-112.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 46-8 (1995).

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Item name	61	AMY EPS
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Condition (Absorbance)</b>		
	Lower	OFF
	Higher	OFF
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	4.3	1700
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Prozone Range</b>		
		OFF

• Prestige 24i Premium, Biolis 24i Premium

Item No.	61	Item Name	AMY EPS	Optical
<b>Data information</b>				
Units	U/l			
Decimals	1			
<b>Calibration</b>				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
<b>Analysis</b>				
Type	RATE			
Main Wave Length	405 nm			
Sub Wave Length	700 nm			
Method	IFCC mod.			
<b>Correlation</b>				
Y=	Slope	X+	Intercept	
	1	0		

Item No.	61	Item Name	AMY EPS	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		7	5	µl
Reagent 1		200	10	µl
Reagent 2		50	10	µl
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				
<b>Data Process</b>				
Read		Start	End	
	Main	41	52	
	Sub			
<b>Abs.Limit</b>				
Low	-3	High	3	
<b>Correction value</b>				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	80			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	61	Item Name	AMY EPS	Optical
<b>Normal Range</b>				
	Male	Female		
	Low High	Low High		
Serum	28 100	28 100		
Urine	0 460	0 460		
Plasma	28 100	28 100		
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male	Female		
	Low High	Low High		
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	61	Item Name	AMY EPS	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	OFF		
	Higher	OFF		
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re Value Dil	Re Value Dil	
Serum		1.1	2000	
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:Di12				

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