



PRESTIGE 24i ASO

DIAGNOSTIC KIT FOR DETERMINATION OF ANTI-STREPTOLYSIN O LEVELS

INTRODUCTION

Most people infected with hemolytic streptococcus produce anti-streptolysin O (ASO), antibodies against streptolysin O (SLO), an exotoxin of Streptococcus. Measuring the level of ASO is effective for diagnosing, judging the progress of medical treatment, and assessing recovery from diseases caused by hemolytic streptococcus such as rheumatic fever, acute glomerulonephritis, scarlatina and tonsillitis.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ASO in a sample and SLO which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ASO in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-270 (24-TRAY)	Cat. No 4-489 (36-TRAY)
1-Reagent	2 x 14 ml	3 x 10 ml
2-Reagent	2 x 20 ml	3 x 13 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks. Protect from light and avoid contamination!

Concentrations in the test

suspension of latex particles sensitized with SLO (pH 8.2) 0.17 w/v%
glycine buffer solution (pH 8.3)

Warnings and notes

- Product for in vitro diagnostic use only.
- Reagent bottles should be shaken before use by gently inverting several times.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- The reagents contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid). After blood has clotted thoroughly, the sample is centrifuged and the serum is separated from blood cells and fibrins. If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

These reagents may be used in automatic analysers Prestige 24i, Biolis 24i, Sapphire 400 and Prestige 24i Premium, Biolis 24i Premium, Sapphire 400 Premium.

- 1-Reagent and 2-Reagent are ready to use.
 - 1-Reagent put on basic position in reagent tray.
 - 2-Reagent put on start position in reagent tray.
- For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ³

serum, plasma	< 160 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples. For the calibration of automatic analysers systems the CORMAY ASO CALIBRATOR kit (Cat. No 4-278) is recommended. The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the automatic analysers Biolis 24i Premium and Hitachi 917. Results may vary if a different instrument is used.

- Sensitivity:** 38 IU/ml.
- Linearity:** up to 1100 IU/ml.
For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.
- Specificity / Interferences**
Haemoglobin up to 0.5 g/dl, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	46.8	1.07	2.29
level 2	80.2	1.31	1.63
level 3	221.7	2.62	1.18

Reproducibility (day to day) n = 12	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	48.8	2.83	5.81
level 2	78.8	2.60	3.30
level 3	219.8	5.24	2.38

Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 50 samples gave following results:

$$y = 1.11 x - 44 \text{ IU/ml};$$

$$R = 0.945 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Galvin J. P. et al.: Particle enhanced photometric immunoassay systems., Clin. Lab. Assays (Pap. Annu. Clin. Lab. Assays Conf.), 4th, 73 (1983).
- Singer J. M. et al.: The latex fixation test. I. Application to the serologic diagnosis of rheumatoid arthritis, Amer. J. Med., 21, 888 (1956).
- Shojiro Kano: antistreptolysin O (ASO), Nippon Rinsho, 57, 108 (1999).

APPLICATION for Prestige 24i, Biolis 24i and Sapphire 400

Item name	51	ASO
Data information		
Units	IU/ml	
Decimals	0	
Analysis		
Type	END	
Main W.Length1	570	
Sub W.Length2	800	
Method	Immuno	
Calibration		
Type	Linear	
Standard		
#1	*	#4
#2		#5
#3		#6
Normal Range		
	Male	
	Low	High
	Female	
	Low	High
Serum	0	160
Urine		
Plasma	0	160
CSF		
Dialysis		
Other		
Corr		
Y=	Slope	X+ Inter
	1.000	0.000

Item name	51	ASO
Aspiration		
Kind	Double	
Data Process		
Read	Start	End
Main	52	54
Sub	37	38
Absorbance Limit		
Low	-3.000	
High	3.000	
Factor	Endpoint Limit 2.000	
Blank correction	Linear Check (%) 0	
Dilution		
Diluent	99:Dil1	
Monitor		
0 Level Point	1	
Span	3.000	
Prozone Check		
First	Start	End
Second		
Third		Low

Item name	51	ASO
Auto Rerun SW		
OFF		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
	Lower	OFF
	Higher	OFF
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
OFF		

APPLICATION for Prestige 24i Premium, Biolis 24i Premium and Sapphire 400 Premium

Item No.	51	Item Name	ASO	Optical
Data information				
Units	IU/ml			
Decimals	0			
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	END method			
Main Wave Length	570 nm			
Sub Wave Length	800 nm			
Method	Immuno			
Correlation				
Slope		Intercept		
Y=	1	X+	0	

Item No.	51	Item Name	ASO	Optical
Aspiration				
Kind	Double			
Data Process				
Read	Start	End		
Main	50	52		
Sub	37	38		
Vol.				
Kind	Vol.	Add	Units	
Sample	3	5	µl	
Reagent 1	160	10	µl	
Reagent 2	110	10	µl	
Abs.Limit				
Low	-3	High 3		
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				
Correction value				
Blank correction		1		
End Point Limit		2		
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Third			Low	

Item No.	51	Item Name	ASO	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0	160	0	160
Urine				
Plasma	0	160	0	160
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	51	Item Name	ASO	Optical
Auto Rerun SW				
ON				
Auto Rerun Condition (Absorbance)				
Lower		OFF		
Higher		OFF		
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
	Re	Value	Dil	Re
Serum		38		1100
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
99:Dil1				

Date of issue: 01. 2013.

MANUFACTURER

PZ CORMAY S.A.
 22 Wiosenna Street,
 05-092 Łomianki, POLAND
 tel.: +48 (0) 22 751 79 10
 fax: +48 (0) 22 751 79 14
<http://www.cormay.pl>

01/13/01/13