IgA

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INTENDED USE

Diagnostic kit for determination of IgA concentration used in automatic analysers according to their user manual.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity. Immunoglobulin A (IgA) is the major Ig found in secretions, playing a major role in the protection of the respiratory, genitourinary, and gastrointestinal tracts against infection.

METHOD PRINCIPLE

The IgA present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to IgA concentration in the sample.

REAGENTS

Package	CORMAY IgA	OS- IgA	B50- IgA
1-Reagent	1 x 48.5 ml	1 x 44.5 ml	1 x 58 ml
2-Reagent	1 x 10 ml	1 x 10.5 ml	1 x 14 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human IgA antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found

to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

ADDITIONAL EQUIPMENT

- automated clinical chemistry analyser capable of accommodating two-reagent assays;
- general laboratory equipment;

SPECIMEN 4

Serum.

Specimen without lipemia or hemolysis is recommended.

Specimen can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Applications for analysers are available on request.

REFERENCE VALUES 4

adults	0.40 – 3.50 g/l
children (1 year – 12 years)	0.15 – 2.50 g/l
children (1 month – 12 months)	0.20 – 0.90 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. Calibrators and 0.9% NaCl should be used for calibration. The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Hitachi or Biolis 24i Premium. Results may vary if a different instrument is used.

- Measurement range:** 0.02 g/l to 8 g/l.

- Specificity / Interferences**

Haemoglobin up to 0.19 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 155.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 0.8 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.30	0.02	1.13
level 2	2.60	0.07	2.67
level 3	3.90	0.17	4.29

Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.95	0.08	8.75
level 2	2.27	0.09	4.17
level 3	3.54	0.10	3.00

- Method comparison**

A comparison between IgA values determined at **Biolis 24i Premium** (y) and at **ADVIA 1650** (x) using 22 samples gave following results:

$$y = 0.9035x + 0.0964 \text{ g/l};$$

$$R = 0.9960 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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