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CORMAY URINE CONTROL LEVEL 1 CORMAY URINE CONTROL LEVEL 2

CORMAY URINE CONTROL LEVEL 1
Cat. No 5-161 3 x 10 ml
CORMAY URINE CONTROL LEVEL 2
Cat. No 5-162 3 x 10 ml

INTENDED USE

CORMAY URINE CONTROL LEVEL 1 and LEVEL 2 was produced on the basis of human urine with added chemicals, preservatives and stabilizers to inhibit microbial growth. This controls are designed for routine quality control urine of organic and inorganic components in urine, which exact target values are given in the enclosed values sheet. Constituent concentrations are available at 2 levels.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INSTRUCTION FOR USE

1. Remove the controls from refrigerator and allow to come to room temperature (18-25°C).
2. Swirl gently to assure homogeneity of the contents. Avoid foaming.
3. Use controls as a patient sample in accordance with the requirements of the test method.
4. Immediately recap the controls and return to 2-8°C when not in use.

STABILITY AND STORAGE

1. Urine control remains stable until expiry date given on the product label when stored unopened at 2-8°C.
2. Store urine control tightly capped, avoid contamination.
3. Discard the controls if excessive turbidity or if there is any evidence of microbial contamination.
4. Product once opened and stored tightly capped at 2-8°C, will be stable for 30 days.

NOTES

- Mean creatinine values listed may decrease over the entire shelf life of the control.
- Use urine control taking into consideration analytical range of particular diagnostic kits for determination of parameters included in values sheet.
- The serum from each donor contributing urine for Urine Control was tested against human immunodeficiency virus (HIV), hepatitis B virus (HBV), hepatitis C virus (HCV) and found to be non-reactive. However Urine Control should be handled as though capable of transmitting infectious disease.
- This product is not intended for use as a standard.

CORMAY URINE CONTROL LEVEL 1 CORMAY URINE CONTROL LEVEL 2

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