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Data wydania: 01.2021

CORMAY URINE PROTEINS CALIBRATORS

Cat. No **5-181** (EN)
CORMAY URINE PROTEINS CALIBRATORS 1 x 1 ml

INTENDED USE

Calibrators for determination of total proteins in urine.
The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

CORMAY URINE PROTEINS CALIBRATORS kit is dedicated for calibration of total protein determinations in urine and cerebrospinal fluid on automatic analysers and manually.
The exact concentration of CORMAY URINE PROTEINS CALIBRATORS for each calibrator lot is shown on vial label.

RECONSTITUTION

Calibrator is ready for use. Before each use mix gently by inverting the vial.

STABILITY AND STORAGE

Unopened calibrator stored at 2-8°C is stable until expiry date given on the product label. After opening, tightly capped calibrator is stable for 30 days when stored at 2-8°C.

Protect from light! Do not freeze!

WARNINGS AND NOTES

▪ A raw material used to manufacture calibrator has been tested for the HBsAg, antibodies to HCV, HIV1, HIV2 and Syphilis with negative result. However the calibrator should be handled as potentially infectious material.

TRACEABILITY

Calibrator value is traceable to certified reference material ERM-DA470k/IFCC (CRM-470).

WASTE MANAGEMENT

Please refer to local legal requirements.

Date of issue: 01.2021

CORMAY URINE PROTEINS CALIBRATORS

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