CORMAY CK-MB CALIBRATOR

Nr kat. **5-182** (PL)

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Cat. No **5-182** (EN)

INTENDED USE

CORMAY CK-MB CALIBRATOR kit is intended for calibration of CK-MB activity determinations on automatic analysers and manual.

The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased CK-MB serum level is a strong marker of myocardial infarction.

PACKAGE

CORMAY CK-MB CALIBRATOR 1 x 1 ml

The exact activity of CK-MB is printed on the label of each bottle.

INSTRUCTION FOR USE

Open the bottle carefully and add exactly **1 ml** of distilled/deionised water, close bottle and let stand for **30 minutes**. Dissolve contents completely by rotating gently, avoiding the formation of foam.

STABILITY AND STORAGE

- Lyophilised serum remains stable when stored at 2-8°C until expiry date given on the product label.
- Stability of the component in the reconstituted calibrator:

Temp. 25°C	8 hours
Temp. 4°C	5 days
Temp. -20°C	4 weeks (when frozen once)
- Store calibrator tightly capped when not in use.

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- The calibrator has been tested for the HIV antibody, HBsAg and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.

TRACEABILITY

Values assigned from the reference material – IFCC-455.

WASTE MANAGEMENT

Please refer to local legal requirements.

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