CORMAY RF LATEX

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INTENDED USE

Agglutination test for detection of rheumatoid factor intended to use for the qualitative and semi-quantitative manual assays. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Rheumatoid factor (RF) is an autoantibody against the Fc portion of the IgG molecule, commonly seen in sera at a high concentration in some conditions, particularly in patients with rheumatoid arthritis. The measurement of RF value is useful in evaluating the diagnosis, effects of therapy and prognosis of RA, systemic lupus erythematosus, chronic hepatopathy, etc.

METHOD PRINCIPLE

Latex particles coated with human gamma-globulin are agglutinated when mixed with samples containing RF.

REAGENTS

Package

CORMAY RF Latex 100
1 x 5 ml
1 x 1 ml
1 x 1 ml
2 x 25 pcs.
2 x 9 pcs.

RF-Latex
RF-Control (+) (red cap)
RF-Control (-) (blue cap)
Stirrers
Slides (6 circles each)

Reagent preparation and stability

The reagents are ready to use.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Do not freeze.

Concentrations in the test

latex particles coated with human gamma-globulin	pH 8.2
human serum solution with a RF concentration	> 30 IU/ml
animal serum solution preservatives	

WARNINGS AND NOTES

- Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.
- Individuals suffering from infectious mononucleosis, hepatitis, syphilis as well as elderly people may give positive results.
- Results obtained with a RF-Latex method do not compare with those obtained with Waaler Rose test. Differences in the results between methods do not reflect differences in the ability to detect rheumatoid factors.
- Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

- The RF-latex sensitivity is calibrated against the RF International Standard from NIBSC 64/002.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- RF-Control (+) and RF-Control (-) meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

RF-Control (+) and RF-Control (-) contain 2-methyl-2H-isothiazol-3-one.

Warning

H317 May cause an allergic skin reaction.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

ADDITIONAL EQUIPMENT

- mechanical rotator with adjustable speed at 80-100 r.p.m.
- general laboratory equipment.

SPECIMEN

Serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

It is recommended to perform the assay with freshly collected samples.

PROCEDURE

Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 50 µl of the sample and one drop of each positive and negative controls into separate circles on the slide test.
- Swirl the RF-Latex reagent gently before using and add one drop (50 µl) next to the sample to be tested.
- Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
- Place the slide on a mechanical rotator at 80-100 r.p.m. **for 2 minutes.** False positive results could appear if the test is read later than two minutes.

Reading and interpretation

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates a RF concentration equal or greater than 8 IU/ml.

Semi-quantitative method

- Make serial two fold dilutions of the sample in 0.9% NaCl solution.
- Proceed for each dilution as in the qualitative method.

Reading and interpretation

The result (titer), in the semi-quantitative method, is defined as the highest dilution showing a positive result.

The approximate RF concentration in the patient sample is calculated as follows:

RF concentration (IU/ml) = 8 x the highest dilution

REFERENCE VALUES

adults	< 8 IU/ml
--------	-----------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

PERFORMANCE CHARACTERISTICS

- Analytical sensitivity:** 8 (6-16) IU/ml.
- Prozone effect:** no prozone effect up to 1500 IU/ml.
- Diagnostic sensitivity:** 100%.
- Diagnostic specificity:** 100 %
- Interferences**
Haemoglobin up to 10 g/l, bilirubin up to 20 mg/dl, triglycerides up to 10 g/l do not interfere with the test.

WASTE MANAGEMENT

Please refer to local legal requirements.

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