RF

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INTENDED USE

Diagnostic kit for determination of rheumatoid factor (RF), concentration concentration intended to use in automatic analysers according to their user manual.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Rheumatoid factor (RF) is an autoantibody against human IgG commonly seen in sera at a high concentration in some conditions, particularly in patients with rheumatoid arthritis. The measurement of RF value is useful in evaluating the diagnosis, effects of therapy and prognosis of RA, systemic lupus erythematosus, chronic hepatopathy, etc.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between RF in a sample and denatured human IgG which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of RF in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

	CORMAY RF	CORMAY RF 60
1-Reagent	1 x 43.5 ml	4 x 60 ml
2-Reagent	1 x 14 ml	2 x 40 ml
	OS-RF	B50-RF
1-REAGENT	2 x 41 ml	2 x 48.5 ml
2-REAGENT	2 x 16 ml	1 x 18.5 ml

The reagents are stable up to the kit expiry date printed on the package when stored at 2-10°C. The reagents are stable for 14 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of latex particles sensitized with denatured human IgG (pH 7.3) 0.17 w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.

- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

ADDITIONAL EQUIPMENT

- automated clinical chemistry analyzer capable of accommodating two-reagent assays;
- general laboratory equipment;

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

Applications for analysers are available on request.

The reagent is ready to use.

REFERENCE VALUES

serum, plasma	< 18 IU/ml
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY RF CALBRATORS kit (Cat. No 4-277) is recommended.

Calibrators and 0.9% NaCl should be used for calibration.

The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analysers Hitachi 912 and Hitachi 917. Results may vary if a different instrument is used.

- Sensitivity:** 8.39 IU/ml.
- Linearity:** up to 120 IU/ml.

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 10 g/l, bilirubin up to 0.66 g/l, intralipid up to 20 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	32,08	0,25	0,78
level 2	112,42	0,32	0,28
Reproducibility (day to day) n = 63	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	20,0	1,03	5,12
level 2	48,0	0,89	1,86

Method comparison

A comparison between RF values determined at **Hitachi 912** (y) and at **Cobas Integra 400 PLUS** (x) using 28 samples gave following results:

$$y = 0.8631x + 0.2352 \text{ IU/ml;}$$

$$R = 0.9999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Galvin J.P. et al.: Particle enhanced photometric immunoassay systems., Clin Lab. Assay (Pap.Annu.Clin.Lab.Assays Conf.), 4th,73 (1983).
- Singer J.M. et al.: The latex fixation test. I. Application to the serologic diagnosis of rheumatoid arthritis, Amer. J. Med., 888 (1956)
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 960, (2006).
- Internal reference range studies.

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