ACCENT-200 BILE ACIDS

Nr kat. 7-198 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of total bile acids concentration used in automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.


INTRODUCTION

Bile acids are the main product of degradation of endogenous cholesterol formed in the liver. Total bile acids are metabolized in the liver and are a valuable indicator of normal or abnormal liver function. Serum total bile acids are increased in patients with viral hepatitis, liver cirrhosis and liver cancer.

METHOD PRINCIPLE

Enzymatic method with 3- α -hydroxysteroid dehydrogenase (3- α -HSD).

Bile acids under the influence of 3-hydroxysteroid dehydrogenase (3- α -HSD) in the presence of thio-NAD are converted to 3-ketosteroids and thio-NADH. The reaction is reversible and 3- α -HSD can convert 3-ketosteroids and NADH to bile acids and NAD.

The rate of thio-NADH formation can be monitored at 405 nm and is proportional to the bile acids activity.

REAGENTS

Package

1-Reagent 2 x 14.5 ml
2-Reagent 2 x 5.4 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 6 weeks.

Concentrations in the test

1-Reagent
Thio-NAD > 0.1 mmol
Buffer
2-Reagent
3- α -HSD > 2 kU/l
NADH > 0.1 mmol
Buffer

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Avoid contact with skin and mucous membranes.
- Yellow or yellow-brown color of the reagent does not affect the reagents performance.
- Do not interchange reagents from different lot numbers.
- Samples from patients treated with ursodeoxycholic acid (UDCA) are not suitable for the determination of total bile acid concentrations.

SPECIMEN

Serum.

Total bile acids concentration is increased after meals, therefore samples should be collected under fasting conditions. Serum samples are stable for a 7 days at 4 °C or for 3 month at - 20 °C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

REFERENCE VALUES ³

serum	2.5 – 6.8 $\mu\text{mol/l}$ (1.25 – 3.4 $\mu\text{g/ml}$)
-------	--

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY BILE ACIDS CONTROLS (Cat. No 5-149) with each batch of samples.

For the calibration of automatic analyser BS-120 the CORMAY BILE ACIDS CALIBRATOR (Cat. No 3-125) is recommended.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY BILE ACIDS CALIBRATOR (Cat. No 3-125) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 6 weeks with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-220S and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity**
0.38 $\mu\text{mol/l}$ (0.19 $\mu\text{g/ml}$) – ACCENT-220S
0.70 $\mu\text{mol/l}$ (0.35 $\mu\text{g/ml}$) – ACCENT-MC240

- Linearity**
up to 150 $\mu\text{mol/l}$ (75 $\mu\text{g/ml}$) – ACCENT-220S
up to 130 $\mu\text{mol/l}$ (65 $\mu\text{g/ml}$) – ACCENT MC240

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 50 mg/dl, ascorbic acid up to 50 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [$\mu\text{mol/l}$]	SD [$\mu\text{mol/l}$]	CV [%]
ACCENT 220S n=10	level 1	31.46	0.30	0.95
	level 2	90.27	0.74	0.82
ACCENT MC240 n=20	level 1	30.44	0.18	0.59
	level 2	97.43	1.08	1.11
Reproducibility (day to day)		Mean [$\mu\text{mol/l}$]	SD [$\mu\text{mol/l}$]	CV [%]
ACCENT 220S n=10	level 1	28.00	0.59	2.10
	level 2	94.46	2.02	2.14
ACCENT MC240 n=80	level 1	30.2	1.34	4.4
	level 2	99.5	3.86	3.9

Method comparison

A comparison between bile acids values determined at ACCENT-220S (y) and at BS-800 (x) using 47 samples gave following results:

$$y = 1.0161x + 0.1825 \mu\text{mol/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between bile acids values determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 40 samples gave following results:

$$y = 0.9475x - 0.3129 \mu\text{mol/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule	
Rule	Two-point Linear
Sensitivity	0
Replicates	3
Interval (day)	42
Difference Limit	0
SD	0
Blank Response	0 500000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters	
Test Name	Bile Acids
Test No	77
Full Name	Bile Acids
Reference No	77
Analy. Type	Kinetic
Pri. Wave.	405 nm
Secon. Wave.	670 nm
Trend	Ascending
Reac. Time	4 16
Incuba. Time	10
Unit	µmol/l
Precision	0.01

R1	270
R2	90
Sample Volume	3
R1 Blank	
Mixed Reag. Blank	
Concentration	1.3 120
Linearity Limit	0.2
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

Calibration Rule	
Rule	Two-point Linear
Sensitivity	0
Replicates	3
Interval (day)	42
Difference Limit	0
SD	0
Blank Response	0 500000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters	
Test Name	Bile Acids
Test No	77
Full Name	Bile Acids
Standard No	77
Reac. Type	Kinetic
Pri. Wave.	405 nm
Sec. Wave.	670 nm
Direction	Increase
Reac. Time	2 16
Incuba. Time	10
Unit	µmol/l
Precision	0.01

R1	240
R2	80
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Linearity Range	0.38 150
Linearity Limit	0.2
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

Calibration Rule	
Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	42
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters	
Test Name	Bile Acids
Test No	77
Full Name	Bile Acids
Standard No	77
Reac. Type	Fixed-time
Pri. Wave.	405 nm
Sec. Wave.	630 nm
Direction	Increase
Reac. Time	10 18
Incuba. Time	16
Unit	µmol/l
Precision	0.01

R1	270
R2	90
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Linearity Range	3.7 130
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

Calibration Rule	
Rule	One-point Linear
Sensitivity	
Replicates	3
Interval (day)	
Difference Limit	
SD	
Blank Response	0 50000
Error Limit	
Coefficient	

• ACCENT S120

Chem	BILE ACIDS	No.	077	Sample Type	SERUM										
Chemistry	BILE ACIDS	Print name	BILE ACIDS	Reaction Direction	positive										
Reaction Type	Kinetic	Sec Wave	578 nm	Decimal	0.01										
Pri Wave	405 nm	Incubation Time	10	Reaction Time	3 10										
Unit	µmol/L	Blank Time		Standard	Sample Vol 2.5 µL Aspirated µL Diluent µL Reagent Vol R1 200 µL R2 40 µL										
Decreased	2.5 µL 20 µL 180 µL	Increased	µL µL µL	Sample Blank	<input type="checkbox"/>										
Auto Rerun	<input checked="" type="checkbox"/>	Linearity range (Standard)	0.7 144	Linearity Limit	0.2										
Linearity Range (Decreased)		Linearity Range (Increased)		Substrate Depletion	40000										
R1 Blank Abs	-40000 40000	Blank Response	-40000 40000	Mixed Blank Abs	-40000 40000										
Twin Chemistry		Control Chemistry		On-board Stability	Day(s)										
Prozone Check	<input type="checkbox"/>	O1		O2		V1		O3		O4		V2			
Q5		Q6		V3		PC1		PC2		Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
Pretreat Sample Vol	µL	Pretreat Sample Vol	µL	CALIBRATION SETTINGS		AUTO CALIBRATION									
Math model	Two-point linear	Factor		Replicates	2	Bottle Changed									
Cal Time	Hour	Slope Diff		SD		Lot Changed									
Sensitivity		Repeatability	40000	Deter Coeff		Cal Time									