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LITERATURA

1. Padhi, S., Nayak, A. K., Behera, A. Type II diabetes mellitus: a review on recent drug based therapeutics. Biomedicine & Pharmacotherapy, 131, 110708 (2020).
2. Pisarczyk-Wiza D., Zozulińska-Ziółkiewicz D. Glikokortykosteroidy a zaburzenia metabolizmu glukozy. Diabetologia Kliniczna 4,3, 110-116 (2015).
3. Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
4. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
5. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
6. Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
7. 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland

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
(EN)

INTENDED USE



GLUCOSE reagent is intended to determine quantitatively glucose level in plasma, serum and cerebrospinal fluid. It is intended to diagnosis, monitoring and as an aid in the diagnosis of clinical conditions connected with abnormal glucose level. GLUCOSE reagent is intended to use on automatic analysers ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180, ACCENT S120 / BS-230, ACCENT MC240 / BS-240Pro, ACCENT M320 / BS-360E, BS-120, ACCENT 400 and ACCENT Neo200. It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY^{1,2}

Determinations of glucose level in serum, plasma and cerebrospinal fluid are mainly used in diagnosis and treatment monitoring of diabetes mellitus. Measurement of glucose concentration is also used as an aid in diagnosis of clinical conditions associated with hypoglycemia and hyperglycemia. Increased glucose levels can be observed in transient hyperglycemia, acute stress reaction, Cushing's syndrome, glucocorticoid therapy, hyperthyroidism. Decreased concentrations of glucose can be associated with neonatal hypoglycemia or hypothyroidism.

METHOD PRINCIPLE³

Colorimetric, enzymatic method with glucose oxidase.

The colour intensity is proportional to the glucose concentration.

REAGENTS

Package

1-REAGENT 4 x 35 ml

Number of tests:

ACCENT-200 (II GEN)	500
ACCENT-220S	500
ACCENT S120	620
ACCENT MC240	620
ACCENT M320	620
BS-120	500

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT

phenol	5 mmol/l
glucose oxidase (GOD)	333.3 µkat/l
peroxidase (POD)	38.33 µkat/l
4- aminoantipyrene (4-AA)	0.75 mmol/l
phosphate buffer	
metal ion chelator	
stabilizer	
preservatives	

REAGENT STABILITY

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent stored on board of the analyzer at 2-10°C is stable for 7 weeks (ACCENT-200 (II GEN)) or 12 weeks (ACCENT MC240).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- EUH208 Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
- EUH210 Safety data sheet available on request.

SPECIMEN^{4,5,6}

EDTA or heparinized plasma in tubes containing sodium fluoride or sodium iodoacetate additive/ serum, free from hemolysis, cerebrospinal fluid.

Plasma / Serum. Serum and plasma specimens should be separated from cells within 30 minutes after collection.

Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level.

Serum and plasma can be stored up to 8 hours at 25°C or 72 hours at 4°C.

Plasma is the specimen recommended for the glucose determination in the blood.

Cerebrospinal fluid. Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

1-REAGENT is ready to use.

Deionised water is recommended as a reagent blank.

For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: ASAT - GLUCOSE, CK - GLUCOSE, AMYLASE EPS - GLUCOSE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analyzers CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 7 weeks (ACCENT-200 (II GEN)) or 12 weeks (ACCENT MC240).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ^{4, 7, 8, 9}

	mg/dl	mmol/l
plasma, serum ^{7,8,9}	70 – 99	3.9 – 5.5
cerebrospinal fluid ⁴	40 – 70	2.2 – 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 (II GEN) and/or ACCENT MC240 and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

LoB (Limit of Blank):

0.0 mg/dl (0.0 mmol/l) – BS-400

LoD (Limit of Detection):

0.3 mg/dl (0.017 mmol/l) – BS-400

LoQ (Limit of Quantitation):

5.5 mg/dl (0.31 mmol/l) - ACCENT-200 (II GEN)

6.0 mg/dl (0.33 mmol/l) - ACCENT MC240

Linearity:

up to 500 mg/dl (27.75 mmol/l) - ACCENT-200 (II GEN)

up to 600 mg/dl (33.30 mmol/l) - ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Measurement range:

ACCENT-200 (II GEN):

5.5 mg/dl (0.31 mmol/l) - 500 mg/dl (27.75 mmol/l)

ACCENT MC240:

6.0 mg/dl (0.33 mmol/l) - 600 mg/dl (33.30 mmol/l)

Specificity / Interferences

a) samples containing low glucose concentration

Haemoglobin up to 0.63 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

b) samples containing high glucose concentration

Haemoglobin up to 2.50 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 (II GEN) n = 20	level 1	84.6	0.45	0.5
	level 2	271.0	1.56	0.6
ACCENT MC240 n = 20	level 1	82.5	0.48	0.6
	level 2	282.1	1.37	0.5
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 (II GEN) n = 80	level 1	87.3	1.11	1.3
	level 2	281.7	4.38	1.6
ACCENT MC240 n = 80	level 1	85.5	1.61	1.9
	level 2	285.0	10.79	3.8

Method comparison

A comparison between glucose values determined at ACCENT-200 (II GEN) (y) and at BECKMAN COULTER AU680 (x) using 60 plasma samples gave following results:

$y = 0.9719x + 0.8776$ mg/dl;

R = 0.999 (R – correlation coefficient)

A comparison between glucose values determined at ACCENT-200 (II GEN) (y) and at BECKMAN COULTER AU680 (x) using 30 serum samples gave following results:

$y = 0.9757x + 1.6903$ mg/dl;

R = 0.999 (R – correlation coefficient)

A comparison between glucose values determined at ACCENT-200 (II GEN) (y) and at BECKMAN COULTER AU680 (x) using 30 CSF samples gave following results:

$y = 0.961x - 0.3386$ mg/dl;

R = 0.998 (R – correlation coefficient)

A comparison between glucose values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 60 plasma samples gave following results:

$y = 1.0018x + 3.2919$ mg/dl;

R = 0.998 (R – correlation coefficient)

A comparison between glucose values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 30 serum samples gave following results:

$y = 1.0273x + 0.0778$ mg/dl;

R = 1.000 (R – correlation coefficient)

A comparison between glucose values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 30 CSF samples gave following results:

$y = 1.0604x - 1.2171$ mg/dl;

R = 0.988 (R – correlation coefficient)

WASTE MANAGEMENT ¹⁰

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: b15 01 02
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ¹¹

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- Padhi, S., Nayak, A. K., Behera, A. Type II diabetes mellitus: a review on recent drug based therapeutics. Biomedicine & Pharmacotherapy, 131, 110708 (2020).
- Pisarczyk-Wiza D., Zozulińska-Ziólkiewicz D. Glikokortykosteroidy a zaburzenia metabolizmu glukozy. Diabetologia Kliniczna 4,3, 110-116 (2015).
- Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
- 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical Chemistry 48, 3, 436-472 (2002).
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clinical Chemistry 50, 1704-1706 (2004).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 06

Current version: 07

Sections updated: WARNING AND NOTES.

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- Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
- 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical Chemistry 48, 3, 436-472 (2002).
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clinical Chemistry 50, 1704-1706 (2004).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	49
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters		R1		R2	
Test	GLUC			250	
No	1				
Full Name	Glucose	Sample Volume		3	
Standard No	1	R1 Blank			
Reac. Type	Endpoint	Mixed Rtg. Blank			
Pri. Wave.	510 nm	Linearity Range		12 450	
Sec. Wave.	670 nm	Linearity Limit			
Direction	Increase	Substrate Limit			
Reac. Time	0 25	Factor			
Incuba. Time	16	<input type="checkbox"/> Prozone check			
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>			
Precision	Integer	PC <input type="checkbox"/>		Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters		R1		R2	
Test	GLUC			250	
No	1				
Full Name	Glucose	Sample Volume		3	
Standard No	1	R1 Blank			
Reac. Type	Endpoint	Mixed Rtg. Blank			
Pri. Wave.	510 nm	Linearity Range		12 480	
Sec. Wave.	670 nm	Linearity Limit			
Direction	Increase	Substrate Limit			
Reac. Time	-2 30	Factor			
Incuba. Time		<input type="checkbox"/> Prozone check			
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>			
Precision	Integer	PC <input type="checkbox"/>		Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

ACCENT-200 GLUCOSE

• ACCENT S120

Chem <input type="text" value="GLUC"/>	No. <input type="text" value="001"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="GLUCOSE"/>	Print name <input type="text" value="GLUCOSE"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="510 nm"/>	Sec Wave <input type="text" value="670 nm"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="0"/>	
	Reaction Time <input type="text" value="28"/> <input type="text" value="30"/>	
Standard <input type="text" value="2.4"/> <input type="text" value="2.4"/> <input type="text" value="20"/>	Aspirated <input type="text" value="20"/> <input type="text" value="180"/>	Diluent <input type="text" value="180"/>
Decreased <input type="text" value="2.4"/> <input type="text" value="20"/>	Reagent Vol R1 <input type="text" value="200"/> <input type="text" value="200"/>	R2 <input type="text" value="200"/>
Increased <input type="text" value="2.4"/> <input type="text" value="20"/>		
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="12"/> <input type="text" value="470"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-40000"/> <input type="text" value="40000"/>
R1 Blank Abs <input type="text" value="-40000"/> <input type="text" value="40000"/>	On-board Stability <input type="text" value="84"/> Day(s)
Blank Response <input type="text" value="-40000"/> <input type="text" value="40000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="40000"/>
Deter Coeff <input type="text"/>

• ACCENT MC240

Chem <input type="text" value="GLUC"/>	No. <input type="text" value="001"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="GLUCOSE"/>	Print name <input type="text" value="GLUCOSE"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="Positive"/>	
Pri Wave <input type="text" value="505 nm"/>	Sec Wave <input type="text" value="660 nm"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="0"/>	
	Reaction Time <input type="text" value="28"/> <input type="text" value="30"/>	
Standard <input type="text" value="2"/> <input type="text" value="2"/> <input type="text" value="20"/>	Aspirated <input type="text" value="20"/> <input type="text" value="180"/>	Diluent <input type="text" value="180"/>
Decreased <input type="text" value="2"/> <input type="text" value="20"/>	Reagent Vol R1 <input type="text" value="200"/> <input type="text" value="200"/>	R2 <input type="text" value="200"/>
Increased <input type="text" value="2"/> <input type="text" value="20"/>		
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="6"/> <input type="text" value="600"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text" value="84"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

• ACCENT M320

Chem	<input type="text" value="GLUC"/>	No.	<input type="text" value="001"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="GLUCOSE"/>	Print name	<input type="text" value="GLUCOSE"/>	Reaction Direction	<input type="text" value="positive"/>
Reaction Type	<input type="text" value="Endpoint"/>	Sec Wave	<input type="text" value="660 nm"/>	Decimal	<input type="text" value="0.1"/>
Pri Wave	<input type="text" value="505 nm"/>	Incubation Time	<input type="text" value="0"/>	Reaction Time	<input type="text" value="32"/> <input type="text" value="34"/>
Unit	<input type="text" value="mg/dl"/>	Reagent Vol	<input type="text" value="200"/> μ L	R1	<input type="text" value="200"/> μ L
Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	R2	<input type="text" value=""/>	R2	<input type="text" value=""/>
Standard	<input type="text" value="2.5"/> μ L	Aspirated	<input type="text" value="20"/> μ L	Diluent	<input type="text" value="180"/> μ L
Decreased	<input type="text" value="2.5"/> μ L	Increased	<input type="text" value=""/>	Sample Blank	<input type="checkbox"/>
Increased	<input type="text" value=""/>	Auto Rerun	<input type="checkbox"/>		

Linearity range (Standard)	<input type="text" value="6"/>	<input type="text" value="500"/>	Linearity Limit	<input type="text" value=""/>	
Linearity Range (Decreased)	<input type="text" value=""/>	<input type="text" value=""/>	Substrate Depletion	<input type="text" value=""/>	
Linearity Range (Increased)	<input type="text" value=""/>	<input type="text" value=""/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	
R1 Blank Abs	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	On-board Stability	<input type="text" value="84"/> Day(s)	
Blank Response	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text" value=""/>	
Twin Chemistry	<input type="text" value=""/>	Enzyme Linear Extension	<input type="checkbox"/>		
Prozone Check	<input type="checkbox"/>				
O1	<input type="text" value=""/>	O2	<input type="text" value=""/>	V1	<input type="text" value=""/>
Q5	<input type="text" value=""/>	Q6	<input type="text" value=""/>	V3	<input type="text" value=""/>
Sample Pretreatment	<input type="text" value=""/>	Control Pretreatment	<input type="text" value=""/>	Calibrator Pretreatment	<input type="text" value=""/>
		Pretreat Sample Vol	<input type="text" value=""/>	μ L	Pretreat Sample Vol
					μ L

CALIBRATION SETTINGS

Math model	<input type="text" value="Multi-point linear"/>	Bottle Changed	<input type="text" value=""/>
Factor	<input type="text" value=""/>	Lot Changed	<input type="text" value=""/>
Replicates	<input type="text" value="2"/>	Cal Time	<input type="text" value=""/>

AUTO CALIBRATION

Bottle Changed

Lot Changed

Cal Time

ACCEPTANCE LIMITS

Cal Time	<input type="text" value="2016"/>	Hour	<input type="text" value=""/>
Slope Diff	<input type="text" value=""/>	SD	<input type="text" value=""/>
Sensitivity	<input type="text" value=""/>	Repeatability	<input type="text" value="35000"/>
Deter Coeff	<input type="text" value=""/>		

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