ACCENT-200 ASAT

Nr kat. **7-214** (PL)

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LITERATURA

1. Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
2. Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7): 1131-1134 (2020).
3. Bergmeyer H.U., IFCC methods for the measurement of catalytic concentrations of enzymes: Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2), Clin. Chem. Acta 105, 147-172 (1980).
4. Bergmeyer HU, Hørder M, Rej R. International Federation of Clinical Chemistry (IFCC) Scientific Committee, Analytical Section: approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2). J Clin Chem Clin Biochem. 1986 Jul;24(7):481-95. PMID: 3734711.
5. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 911 (1989).

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
(EN)

INTENDED USE


ASAT reagent is intended to determine quantitatively aspartate aminotransferase (ASAT) activity in serum. It is used for monitoring and as an aid to diagnosis of clinical conditions associated with abnormal aspartate aminotransferase activity. ASAT reagent is intended to use on automatic analysers: ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180, ACCENT S120 / BS-230, ACCENT MC240 / BS-240Pro, ACCENT M320 / BS-360E, BS-120, ACCENT 400 and ACCENT Neo200. It is only for *in vitro* diagnostics, for healthcare professional users.

INTRODUCTION^{1,2}

Determinations of elevated aspartate aminotransferase activity (ASAT, AST, GOT) in serum are mainly used to aid in the diagnosis and monitor of liver diseases e.g. hepatitis, hepatic cirrhosis, drug-induced liver injury, hepatic metastasis, mononucleosis. Also skeletal muscle diseases such as muscle trauma, surgery, burns, muscular dystrophy, heat stroke and other diseases such as acute hemolytic anemia, acute pancreatitis and patients with higher risk of COVID-19 adverse outcome. Decreased levels of ASAT can be observed in acute renal disease, beriberi, diabetic ketoacidosis, pregnancy, chronic renal dialysis.

METHOD PRINCIPLE³

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to aspartate aminotransferase activity.

REAGENTS

Package

1-REAGENT	2 x 35 ml
2-REAGENT	1 x 17.5 ml

Number of tests

ACCENT-200 (II GEN)	310
ACCENT-220S	310
ACCENT S120	490
ACCENT MC240	490
ACCENT M320	370
BS-120	310

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT

L-aspartate	300 mmol/l
LDH	1.95 U/ml
MDH	1.15 U/ml
buffer Tris	
pH adjuster	
stabilizers	
preservative	

2-REAGENT

2-oxoglutarate	62 mmol/l
NADH	1.4 mmol/l
buffer	
pH adjuster	
preservatives	

REAGENT STABILITY

The reagents, stored at 2-8°C is stable up to the kit expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (ACCENT-200 (II GEN), ACCENT M320, ACCENT S120).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

- H315 Causes skin irritation.
- H319 Causes serious eye irritation.
- P280 Wear protective gloves, protective clothing, eye protection or face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN^{5, 6, 7, 8}

Serum free from hemolysis. Hemolysis should be avoided, since ASAT activity in erythrocytes is 15 times higher than in normal serum. ASAT activity remains stable in specimen up to 4 days at 15-25°C, up to 7 days at 2-8°C or up to 3 months at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples! Follow tube manufacturers' instructions carefully when using collection tubes .

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

PROCEDURE

Reagents are ready to use.

Deionised water is recommended as a reagent blank.

For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: ASAT - CHOL, ASAT - GLUCOSE, ASAT - LDH, ASAT - CK-MB. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200 (II GEN), ACCENT M320, ACCENT S120). Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ⁹

serum	37°C	
female	up to 31 U/l	up to 0.53 µkat/l
male	up to 35 U/l	up to 0.60 µkat/l

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-200 (II GEN), ACCENT MC240 i/lub BS-400. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
1.8 U/l (0.03 µkat/l) – BS-400
- LoD (Limit of Detection):**
3.0 U/l (0.05 µkat/l) – BS-400
- LoQ (Limit of Quantitation):**
7.0 U/l (0.12 µkat/l) – ACCENT-200 (II GEN), ACCENT MC240
- Linearity:**
up to 1000 U/l (16.67 µkat/l) – ACCENT-200 (II GEN), ACCENT MC240
- Measurement range:**
7.0 U/l (0.12 µkat/l) – 1000 U/l (16.67 µkat/l) – ACCENT-200 (II GEN), ACCENT MC240
- Specificity / Interferences**
Ascorbic acid up to 62 mg/l, triglycerides up to 1000 mg/dl, bilirubin up to 20 mg/dl, haemoglobin up to 0.31 g/dl (in samples containing low activity of ASAT) and up to 2.5 g/dl (in samples containing high activity of ASAT) do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 (II GEN) n=20	level 1	44.5	1.28	2.9
	level 2	218.0	1.80	0.8
ACCENT MC240 n=20	level 1	43.5	1.51	3.5
	level 2	204.9	2.01	1.0
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 (I GEN) n=80	level 1	42.2	1.92	4.5
	level 2	209.2	4.41	2.1
ACCENT MC240 n=80	level 1	43.6	1.65	3.8
	level 2	211.3	4.53	2.1

Method comparison

A comparison between ASAT values determined at **ACCENT-200 (II GEN)** (y) and at **BECKMAN COULTER AU68** (x) using 60 serum samples gave following results:

$$y = 0.9761x - 2.6497 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between ASAT values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU480** (x) using 60 serum samples gave following results:

$$y = 0.9415x - 0.4718 \text{ U/l};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT ¹⁰

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ¹¹

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cornay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- Pagana K, Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
- Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
- Bergmeyer H.U., IFCC methods for the measurement of catalytic concentrations of enzymes: Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2), Clin. Chem. Acta 105, 147-172 (1980).
- Bergmeyer HU, Hørdler M, Rej R. International Federation of Clinical Chemistry (IFCC) Scientific Committee, Analytical Section: approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2). J Clin Chem Clin Biochem. 1986 Jul;24(7):481-95. PMID: 3734711.
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 911 (1989). MAT BIOLOG
- NCCLS. Procedures for the Handling and Processing of Blood Specimens: Approved Guideline—Third Edition. NCCLS document H18-A3 (ISBN 1-56238-555-0). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA (2004).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 413, 1800 (2018).
- World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
- Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2257, (2006).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.

11. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 06	Current version: 07
Sections updated: WARNING AND NOTES.	

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- Bergmeyer H.U., IFCC methods for the measurement of catalytic concentrations of enzymes: Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2), Clin. Chem. Acta 105, 147-172 (1980).
- Bergmeyer HU, Hörder M, Rej R. International Federation of Clinical Chemistry (IFCC) Scientific Committee, Analytical Section: approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2). J Clin Chem Clin Biochem. 1986 Jul;24(7):481-95. PMID: 3734711.
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 911 (1989). MAT BIOLOG
- NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition. NCCLS document H18-A3 (ISBN 1-56238-555-0). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA (2004).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 413, 1800 (2018).
- World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
- Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2257, (2006).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters	ASAT II GEN	R1	200
Test	ASAT II GEN	R1	200
No	16	R2	50
Full Name	ASAT	Sample Volume	12
Standard No	16	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	340 nm	Linearity Range	7 1000
Sec. Wave.	450 nm	Linearity Limit	0.2
Direction	Decrease	Substrate Limit	5200
Reac. Time	3 15	Factor	
Incuba. Time	11	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters	ASAT II GEN	R1	200
Test	ASAT II GEN	R1	200
No	16	R2	50
Full Name	ASAT	Sample Volume	10
Standard No	16	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	340 nm	Linearity Range	7 1000
Sec. Wave.	450 nm	Linearity Limit	0.2
Direction	Decrease	Substrate Limit	5400
Reac. Time	3 12	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT M320

Chem	<input type="text" value="ASAT II GEN"/>	No.	<input type="text" value="016"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="ASAT"/>	Print name	<input type="text" value="ASAT"/>	Reaction Direction	<input type="text" value="negative"/>
Reaction Type	<input type="text" value="Kinetic"/>	Pri Wave	<input type="text" value="340 nm"/>	Sec Wave	<input type="text" value="450 nm"/>
Unit	<input type="text" value="U/l"/>	Decimal	<input type="text" value="0.1"/>	Incubation Time	<input type="text" value="20"/>
Blank Time	<input type="text"/>	Reaction Time	<input type="text" value="10"/>	<input type="text" value="25"/>	
Standard	<input type="text" value="8"/> μ L	Aspirated	<input type="text"/>	Diluent	<input type="text"/>
Decreased	<input type="text" value="8"/> μ L	<input type="text" value="20"/> μ L	<input type="text" value="180"/> μ L	Reagent Vol	<input type="text" value="160"/> μ L
Increased	<input type="text"/>	<input type="text"/>	<input type="text"/>	R2	<input type="text" value="40"/> μ L
		<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun		
Linearity range (Standard)	<input type="text" value="7"/>	<input type="text" value="1200"/>	Linearity Limit	<input type="text" value="0.8"/>	
Linearity Range (Decreased)	<input type="text"/>	<input type="text"/>	Substrate Depletion	<input type="text" value="4000"/>	
Linearity Range (Increased)	<input type="text"/>	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/>	<input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	On-board Stability	<input type="text" value="84"/> Day(s)	
Blank Response	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>	
Twin Chemistry	<input type="checkbox"/> Enzyme Linear Extension				
<input type="checkbox"/> Prozone Check					
Q1	<input type="text"/>	Q2	<input type="text"/>	V1	<input type="text"/>
Q3	<input type="text"/>	Q4	<input type="text"/>	V2	<input type="text"/>
O5	<input type="text"/>	O6	<input type="text"/>	V3	<input type="text"/>
PC1	<input type="text"/>	PC2	<input type="text"/>	PC3	<input type="text"/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
<input type="text"/>	Pretreat Sample Vol <input type="text"/>	μ L	Pretreat Sample Vol <input type="text"/>	μ L	
CALIBRATION SETTINGS			AUTO CALIBRATION		
Math model	<input type="text" value="Multi-point linear"/>				
Factor	<input type="text"/>	Replicates	<input type="text" value="2"/>	<input type="checkbox"/> Bottle Changed	
ACCEPTANCE LIMITS			<input type="checkbox"/> Lot Changed		
Cal Time	<input type="text" value="2016"/> Hour	<input type="checkbox"/> Cal Time			
Slope Diff	<input type="text"/>	SD	<input type="text"/>		
Sensitivity	<input type="text"/>	Repeatability	<input type="text" value="35000"/>		
Deter Coeff	<input type="text"/>				

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