ACCENT-200 HAPTOGLOBIN

Nr kat. 7-215 (PL)

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Cat. No 7-215 (EN)

INTENDED USE

Diagnostic kit for determination of haptoglobin concentration used in automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Haptoglobin is an acute phase protein whose primary function consists in binding free haemoglobin in serum. The complex is removed within minutes by the reticulo-endothelial system where its components are metabolised to free aminoacids and iron. Haptoglobin consequently plays a major role in preventing the loss of haemoglobin in urine and the consequent iron loss from the iron pool. The haptoglobin levels are increased during the acute phase and in such conditions as burns or nephrotic syndrome. Haptoglobin levels are abnormally high in intravascular hemolysis and when haemoglobin turnover is increased such as during haemolytic anaemia, transfusion reactions and malaria.

METHOD PRINCIPLE

The haptoglobin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to haptoglobin concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 35 ml
 2-Reagent 1 x 7 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-Reagent: anti-human haptoglobin antiserum, inorganic salt, buffer, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, fresh serum is recommended. Samples may be stored up to 2 weeks at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

adults	0.26 – 1.85 g/l
newborns	0.05 – 0.48 g/l
> 60 year M	0.35 – 1.64 g/l
> 60 year F	0.40 – 1.75 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems ACCENT-200, ACCENT-200 II GEN, ACCENT MC240 and ACCENT M320 the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers Hitachi or ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- LOQ:** 0.002 g/l
- Linearity:** up to 3.86 g/l.
- Interferences:**
 Hemoglobin up to 0.06 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 1669.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.49	0.007	1.44
level 2	1.31	0.008	0.63
level 3	2.14	0.017	0.77
Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.62	0.018	2.84
level 2	1.48	0.037	2.48
level 3	2.27	0.050	2.21

A comparison between haptoglobin values determined at ACCENT MC240 (y) and at BS-400 (x) using 58 serum samples gave following results:

$$y = 0.9792 x + 0.0232 \text{ g/l;}$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT MC240

Chem	HAPTO	No.	046	Sample Type	SERUM
Chemistry	HAPTOGLOBIN	Print name	HAPTO	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700nm	Decimal	0.01
Pri Wave	340nm	Incubation Time	21	Reaction Time	9 10
Unit	g/l	Sample Vol	2	Aspirated	
Blank Time	-3 -1	Diluent		Reagent Vol	R1 160 R2 32
Standard	2	Decreased	2	Increased	
Decreased	2	Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>
Increased		Linearity range (Standard)	0.05 3.20	Linearity Limit	
Linearity Range (Decreased)		Linearity Range (Increased)		Substrate Depletion	
Linearity Range (Increased)		R1 Blank Abs	-35000 35000	Mixed Blank Abs	-35000 35000
Blank Response	-35000 35000	Blank Response	-35000 35000	On-board Stability	
Twin Chemistry		Prozone Check	<input type="checkbox"/>	Enzyme Linear Extension	<input type="checkbox"/>
Q1		Q2		Q3	
Q5		Q6		PC1	
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
Pretreat Sample Vol		Pretreat Sample Vol		Pretreat Sample Vol	
Math model	Logit-Log 5P	Factor		Replicates	2
ACCEPTANCE LIMITS	Cal Time	Hour		Slope Diff	
Sensitivity	Repeatability	SD		35000	
Deter Coeff					
AUTO CALIBRATION	Bottle Changed	Lot Changed		Cal Time	

• ACCENT M320

ACCENT-200 HAPTOGLOBIN

Chem	<input type="text" value="HAPTO"/>	No.	<input type="text" value="046"/>	Sample Type	<input type="text" value="SERUM"/>	
Chemistry	<input type="text" value="HAPTOGLOBIN"/>			Print name	<input type="text" value="HAPTO"/>	
Reaction Type	<input type="text" value="Endpoint"/>			Reaction Direction	<input type="text" value="positive"/>	
Pri Wave	<input type="text" value="340nm"/>			Sec Wave	<input type="text" value="700 nm"/>	
Unit	<input type="text" value="g/l"/>			Decimal	<input type="text" value="0.001"/>	
Blank Time	<input type="text" value="-3"/>	<input type="text" value="-1"/>		Incubation Time	<input type="text" value="25"/>	
				Reaction Time	<input type="text" value="10"/>	<input type="text" value="11"/>
	<input type="text" value="2"/>	<input type="text" value="20"/>	<input type="text" value="180"/>	Reagent Vol	<input type="text" value="160"/>	<input type="text" value="32"/>
Standard	<input type="text" value="2"/>	<input type="text" value="20"/>	<input type="text" value="180"/>	R1	<input type="text" value="160"/>	<input type="text" value="32"/>
Decreased	<input type="text" value="2"/>	<input type="text" value="20"/>	<input type="text" value="180"/>	R2	<input type="text" value="32"/>	
Increased	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>			
	<input type="checkbox"/>	Sample Blank	<input checked="" type="checkbox"/>	Auto Rerun		

Linearity range (Standard)	<input type="text" value="0.05"/>	<input type="text" value="3.19"/>	Linearity Limit	<input type="text" value=""/>							
Linearity Range (Decreased)	<input type="text" value=""/>	<input type="text" value=""/>	Substrate Depletion	<input type="text" value=""/>							
Linearity Range (Increased)	<input type="text" value=""/>	<input type="text" value=""/>	Mixed Blank Abs	<input type="text" value="-35000"/>	<input type="text" value="35000"/>						
R1 Blank Abs	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	On-board Stability	<input type="text" value=""/>	Day(s)						
Blank Response	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text" value=""/>							
Twin Chemistry	<input type="text" value=""/>		<input type="checkbox"/>	Enzyme Linear Extension							
	<input type="checkbox"/>	Prozone Check									
Q1	<input type="text" value=""/>	Q2	<input type="text" value=""/>	V1	<input type="text" value=""/>	Q3	<input type="text" value=""/>	Q4	<input type="text" value=""/>	V2	<input type="text" value=""/>
Q5	<input type="text" value=""/>	Q6	<input type="text" value=""/>	V3	<input type="text" value=""/>	PC1	<input type="text" value=""/>	PC2	<input type="text" value=""/>		<input type="text" value=""/>
<input type="checkbox"/>	Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment						
	<input type="text" value=""/>		<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/>	<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>

CALIBRATION SETTINGS	AUTO CALIBRATION		
Math model	<input type="text" value="Logit-Log 5P"/>	<input type="checkbox"/>	Bottle Changed
Factor	<input type="text" value=""/>	<input type="checkbox"/>	Lot Changed
Replicates	<input type="text" value="2"/>	<input type="checkbox"/>	Cal Time

ACCEPTANCE LIMITS			
Cal Time	<input type="text" value=""/>	Hour	
Slope Diff	<input type="text" value=""/>		
Sensitivity	<input type="text" value=""/>	Repeatability	<input type="text" value="35000"/>
Deter Coeff	<input type="text" value=""/>		

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