ACCENT-200 CERULOPLASMIN

Nr kat. 7-218 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of ceruloplasmin concentration used in automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Ceruloplasmin is an α_2 -glycoprotein containing 6-7 atoms of copper per molecule. It has been considered for a long time as a copper transporter but since recently, it has been shown to be a serum ferroxidase playing a major role in oxidising iron (II) to iron (III) in serum and at the cell surface, thereby regulating its binding by transferrin.

Ceruloplasmin is a late acute phase reactant. Low levels are found in malnutrition, nephrosis, severe liver diseases such as primary biliary cirrhosis and in Wilson's disease, an autosomal recessive defect in the regulation of copper metabolism.

METHOD PRINCIPLE

The ceruloplasmin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to ceruloplasmin concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 35 ml
 2-Reagent 1 x 7 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: MES buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human ceruloplasmin antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store reagents closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁷

Serum.

Specimen without lipemia or hemolysis is recommended. Samples remain stable for 3 days at 2-8°C or 4 weeks at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁷

age	g/l	mg/dl
1 day – 3 months	0.05 – 0.18	5 – 18
6 – 12 months	0.33 – 0.43	33 – 43
1 – 7 years	0.24 – 0.56	24 – 56
> 7 years	0.18 – 0.45	18 – 45

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi or ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- LOQ:** 0.948 mg/dl (0.00948 g/l)
- Linearity:** up to 105.12 mg/dl (1.0512 g/l).

Specificity / Interferences

Hemoglobin up to 0.19 g/dl, bilirubin up to 22 mg/dl triglycerides up to 1669,5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.25	0.002	0.6
level 2	0.40	0.003	0.6
level 3	0.52	0.004	0.7

Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.19	0.006	3.0
level 2	0.36	0.006	1.6
level 3	0.49	0.006	1.3

Method comparison

A comparison between ceruloplasmin values determined at ACCENT MC240 (y) and at BS-400 (x) using 60 serum samples gave following results:

$$y = 0.9308x + 0.0231 \text{ g/l;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT MC240

Chem	CERULO	No.	045	Sample Type	SERUM
Chemistry	CERULOPLASMIN	Print name	CERULO	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700nm	Decimal	0.001
Pri Wave	340nm	Incubation Time	21	Reaction Time	9 10
Unit	g/l	Blank Time	-3 -1	Standard	4.5 μL
		Sample Vol	Aspirated	Diluent	Reagent Vol
		Decreased	4.5 μL	20 μL	180 μL
		Increased	μL	μL	μL
			μL	μL	μL
			Sample Blank	V	Auto Rerun

Linearity range (Standard)	0.02	1.26	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000	
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)	
Blank Response	-35000	35000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
<input type="checkbox"/> Prozone Check					
Q1	Q2	V1	Q3	Q4	V2
Q5	Q6	V3	PC1	PC2	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL	

CALIBRATION SETTINGS		AUTO CALIBRATION	
Math model	Logit-Log 5P	<input type="checkbox"/> Bottle Changed	
Factor	Replicates 2	<input type="checkbox"/> Lot Changed	
ACCEPTANCE LIMITS		<input type="checkbox"/> Cal Time	
Cal Time	Hour		
Slope Diff	SD		
Sensitivity	Repeatability 35000		
Deter Coeff			