ACCENT-200 CK

Nr kat. **7-220** (PL)

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Cat. No **7-220** (EN)

INTENDED USE

Diagnostic kit for determination of creatine kinase activity intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT-400 and ACCENT Neo200.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased level of CK is usually the result of muscle injury, myocardial or pulmonary infarction.

METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC).

The rate of absorbance changes at $\lambda=340$ nm is directly proportional to creatine kinase activity.

REAGENTS

Package
1-REAGENT 1 x 30 ml
2-REAGENT 1 x 7 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks (ACCENT-200) or 12 weeks (ACCENT M320, ACCENT S120).

Concentrations in the test

1-REAGENT
imidazole buffer 100 mmol/l
glucose 20 mmol/l
N-acetylcysteine 20 mmol/l
magnesium acetate 10 mmol/l
EDTA 2 mmol/l
NADP 2 mmol/l
ADP 2 mmol/l
AMP 5 mmol/l
HK > 2.5 U/ml

2-REAGENT
diadenosinepentaphosphate 10 µmol/l
glucose-6-phosphate-dehydrogenase (G6P-DH) > 1.5 U/ml
creatine phosphate 30 mmol/l
preservative

Warnings and notes

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- Do not freeze reagents.
- Protect from direct sunlight and avoid contamination!
- Do not use reagents past the expiry date.
- Do not interchange caps among reagents.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains imidazole.

Danger

H360D: May damage the unborn child.
P201: Obtain special instructions before use.
P308+P313: IF exposed or concerned: Get medical advice.
P405: Store locked up.

SPECIMEN

Serum, free from hemolysis.
CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light.
Specimens can be stored up to 4-8 hours at 15-25°C or 1-2 days at 2-8°C or 1 month at -20°C.
Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.
Deionized water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: CK - CHOL, CK - GLUCOSE, CK - CREATININE, CK - MG, CK-MB - CK. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ⁸

serum	37°C	
female	< 167 U/l	< 2.78 µkat/l
male	< 190 U/l	< 3.17 µkat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analysers: ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 4 weeks (ACCENT-200) or every 12 weeks (ACCENT M320, ACCENT S120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity**
45 U/l (0.75 µkat/l) - ACCENT 200
- LOQ (Limit of Quantitation)**
3.7 U/l (0.062 µkat/l) – ACCENT MC240
- Linearity**
up to 2000 U/l (33.3 µkat/l) - ACCENT 200
up to 2710 U/l (45.17 µkat/l) – ACCENT MC240

For higher activity, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 0,156 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=20	level 1	129.45	1.19	0.92
	level 2	530.00	7.16	1.35
ACCENT MC240 n=20	level 1	124.9	1.21	1.0
	level 2	424.7	3.43	0.8
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=20	level 1	131.28	4.30	3.28
	level 2	481.35	10.37	2.15
ACCENT MC240 n=80	level 1	120.7	2.54	2.1
	level 2	425.1	6.34	1.5

Method comparison

A comparison between CK values determined at ACCENT-200 (y) and at COBAS INTEGRA 400 (x) using 56 samples gave following results:

$$y = 0.9913x + 0.1166 \text{ U/l};$$

$$R = 0.9996 \quad (R - \text{correlation coefficient})$$

A comparison between CK values determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 63 serum samples gave following results:

$$y = 0.8972x + 12.506 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	<input type="text" value="One-Point Linear"/>
Sensitivity	<input type="text" value="1"/>
Replicates	<input type="text" value="3"/>
Interval (day)	<input type="text" value="28"/>
Difference Limit	<input type="text" value="0"/>
SD	<input type="text" value="0"/>
Blank Response	<input type="text" value="0"/> <input type="text" value="50000"/>
Error Limit	<input type="text" value="0"/>
Coefficient	<input type="text" value="0"/>

• **ACCENT-220S**

Parameters

Test	<input type="text" value="CK"/>	R1	<input type="text" value="200"/>
No	<input type="text" value="19"/>	R2	<input type="text" value="40"/>
Full Name	<input type="text" value="Creatine kinase"/>	Sample Volume	<input type="text" value="10"/>
Standard No	<input type="text" value="19"/>	R1 Blank	<input type="text"/>
Reac. Type	<input type="text" value="Kinetic"/>	Mixed Rtg. Blank	<input type="text"/>
Pri. Wave.	<input type="text" value="340 nm"/>	Linearity Range	<input type="text" value="7.4"/> <input type="text" value="2100"/>
Sec. Wave.	<input type="text" value="450 nm"/>	Linearity Limit	<input type="text" value="0.2"/>
Direction	<input type="text" value="Increase"/>	Substrate Limit	<input type="text"/>
Reac. Time	<input type="text" value="8"/> <input type="text" value="15"/>	Factor	<input type="text"/>
Incuba. Time	<input type="text" value="21"/>	<input type="checkbox"/> Prozone check	
Unit	<input type="text" value="U/l"/>	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	<input type="text" value="0.1"/>	PC <input type="text"/>	Abs <input type="text"/>

Calibration Rule

Rule	<input type="text" value="Multi-point Linear"/>
Sensitivity	<input type="text" value="1"/>
Replicates	<input type="text" value="3"/>
Interval (day)	<input type="text" value="28"/>
Difference Limit	<input type="text" value="0"/>
SD	<input type="text" value="0"/>
Blank Response	<input type="text" value="0"/> <input type="text" value="50000"/>
Error Limit	<input type="text" value="0"/>
Coefficient	<input type="text" value="0"/>

• **BS-120**

Parameters

Test	<input type="text" value="CK"/>	R1	<input type="text" value="200"/>
No	<input type="text" value="19"/>	R2	<input type="text" value="40"/>
Full Name	<input type="text" value="Creatine kinase"/>	Sample Volume	<input type="text" value="10"/>
Standard No	<input type="text" value="19"/>	R1 Blank	<input type="text"/>
Reac. Type	<input type="text" value="Kinetic"/>	Mixed Rtg. Blank	<input type="text"/>
Pri. Wave.	<input type="text" value="340 nm"/>	Linearity Range	<input type="text" value="10"/> <input type="text" value="1500"/>
Sec. Wave.	<input type="text" value="450 nm"/>	Linearity Limit	<input type="text" value="0.2"/>
Direction	<input type="text" value="Increase"/>	Substrate Limit	<input type="text"/>
Reac. Time	<input type="text" value="10"/> <input type="text" value="18"/>	Factor	<input type="text"/>
Incuba. Time	<input type="text" value="16"/>	<input type="checkbox"/> Prozone check	
Unit	<input type="text" value="U/l"/>	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	<input type="text" value="Integer"/>	PC <input type="text"/>	Abs <input type="text"/>

Calibration Rule

Rule	<input type="text" value="Multi-point Linear"/>
Sensitivity	<input type="text" value="1"/>
Replicates	<input type="text" value="3"/>
Interval (day)	<input type="text"/>
Difference Limit	<input type="text" value="0"/>
SD	<input type="text" value="0"/>
Blank Response	<input type="text" value="0"/> <input type="text" value="50000"/>
Error Limit	<input type="text" value="0"/>
Coefficient	<input type="text" value="0"/>

• **ACCENT S120**

Chem <input type="text" value="CK"/>	No. <input type="text" value="019"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="CK"/>		Print name <input type="text" value="CK"/>
Reaction Type <input type="text" value="Kinetic"/>		Reaction Direction <input type="text" value="positive"/>
Pri Wave <input type="text" value="340 nm"/>		Sec Wave <input type="text" value="450 nm"/>
Unit <input type="text" value="U/L"/>		Decimal <input type="text" value="0.1"/>
Blank Time <input type="text"/>		Incubation Time <input type="text" value="21"/>
		Reaction Time <input type="text" value="10"/> <input type="text" value="18"/>
Standard <input type="text" value="4.5"/> <input type="text" value="µL"/>	Aspirated <input type="text" value="20"/> <input type="text" value="µL"/>	Diluent <input type="text" value="180"/> <input type="text" value="µL"/>
Decreased <input type="text" value="4.5"/> <input type="text" value="µL"/>		Reagent Vol
Increased <input type="text" value="µL"/>		R1 <input type="text" value="150"/> <input type="text" value="µL"/>
		R2 <input type="text" value="30"/> <input type="text" value="µL"/>
	<input type="checkbox"/> Sample Blank	<input type="checkbox"/> V
		Auto Rerun

Linearity range (Standard)	<input type="text" value="3.8"/> <input type="text" value="2700"/>	Linearity Limit	<input type="text" value="0.2"/>
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text" value="40000"/>
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-40000"/> <input type="text" value="40000"/>
R1 Blank Abs	<input type="text" value="-40000"/> <input type="text" value="40000"/>	On-board Stability	<input type="text" value="84"/> Day(s)
Blank Response	<input type="text" value="-40000"/> <input type="text" value="40000"/>	Reagent Alarm Limit	<input type="text"/>
Twin Chemistry	<input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check			
Q1 <input type="text"/>	Q2 <input type="text"/>	V1 <input type="text"/>	Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>
Q5 <input type="text"/>	Q6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/> PC2 <input type="text"/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment	
	Pretreat Sample Vol <input type="text"/> <input type="text"/>	µL	Pretreat Sample Vol <input type="text"/> <input type="text"/>
		µL	

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text"/>
Sensitivity <input type="text"/>
Deter Coeff <input type="text"/>
SD <input type="text"/>
Repeatability <input type="text" value="40000"/>

ACCENT-200 CK

• ACCENT MC240

Chem <input type="text" value="CK"/>	No. <input type="text" value="019"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="CK"/>	Print name <input type="text" value="CK"/>	
Reaction Type <input type="text" value="Kinetic"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="450 nm"/>	
Unit <input type="text" value="U/L"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text"/>	Incubation Time <input type="text" value="21"/>	Reaction Time <input type="text" value="10"/> <input type="text" value="18"/>
Standard <input type="text" value="4.5"/> <input type="text" value="μL"/>	Aspirated <input type="text"/>	Diluent <input type="text"/>
Decreased <input type="text" value="4.5"/> <input type="text" value="μL"/>	<input type="text" value="20"/> <input type="text" value="μL"/>	<input type="text" value="180"/> <input type="text" value="μL"/>
Increased <input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> Sample Blank	<input type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="3.7"/> <input type="text" value="2710"/>	Linearity Limit <input type="text" value="0.2"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text" value="35000"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text" value="84"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
O1 <input type="text"/> O2 <input type="text"/> V1 <input type="text"/> O3 <input type="text"/> O4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> <input type="text" value="μL"/>	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> <input type="text" value="μL"/>	<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> <input type="text" value="μL"/>

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text"/>
Sensitivity <input type="text"/>
Deter Coeff <input type="text"/>
SD <input type="text"/>
Repeatability <input type="text" value="35000"/>

• ACCENT M320

Chem <input type="text" value="CK"/>	No. <input type="text" value="019"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="CK"/>	Print name <input type="text" value="CK"/>	
Reaction Type <input type="text" value="Kinetic"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="450 nm"/>	
Unit <input type="text" value="U/L"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text"/>	Incubation Time <input type="text" value="25"/>	Reaction Time <input type="text" value="20"/> <input type="text" value="30"/>
Standard <input type="text" value="4.5"/> <input type="text" value="μL"/>	Aspirated <input type="text"/>	Diluent <input type="text"/>
Decreased <input type="text" value="4.5"/> <input type="text" value="μL"/>	<input type="text" value="20"/> <input type="text" value="μL"/>	<input type="text" value="180"/> <input type="text" value="μL"/>
Increased <input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> Sample Blank	<input type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="8"/> <input type="text" value="2200"/>	Linearity Limit <input type="text" value="0.2"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text" value="35000"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text" value="84"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> <input type="text" value="μL"/>	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> <input type="text" value="μL"/>	<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> <input type="text" value="μL"/>

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text"/>
Sensitivity <input type="text"/>
Deter Coeff <input type="text"/>
SD <input type="text"/>
Repeatability <input type="text" value="35000"/>

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