ACCENT-200 LIPOPROTEIN (a)

Nr kat. 7-228 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of lipoprotein (a) concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lipoprotein(a) is a complex, cholesterol-carrying particle in the blood related to LDL which, when present at high levels, may be associated with the development of atherosclerosis and coronary heart disease, independent of LDL cholesterol and apoB. The structural component of Lp(a) distinguishing it from LDL is apolipoprotein(a), a large protein attached by disulfide bonding to the apoB-100 component of LDL. The similarity of the apo(a) sequence to those of plasminogen and hepatocyte growth factor suggests that the role of lipoprotein(a) in promoting the development of atherosclerosis may come from its capability to:

- 1) interfere in the breakdown of blood clots,
- 2) stimulate atherosclerotic cell proliferation.

Lp(a) levels is largely due to hereditary factors and considered to be useful in assessment of atherosclerosis risks.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between Lp(a) in a sample and anti-Lp(a) antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (700 nm), with the magnitude of the change being proportional to the quantity of Lp(a) in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-Reagent 1 x 25 ml
2-Reagent 1 x 12.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package.

Concentrations in the test

suspension of latex particles sensitized with anti-Lp(a) antibodies (rabbit) (pH 7.3) 0.4 w/v%
glycine buffer solution (pH 9.0)
preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

serum, plasma	< 30 mg/dl
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY Lp(a) CONTROL N (Cat. No 4-492) and CORMAY Lp(a) CONTROL P (Cat. No 4-493).

For the calibration of automatic analysers systems the CORMAY Lp(a) CALIBRATORS kit (Cat. No 4-281) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared each time the test is performed.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser ACCENT-200, ACCENT-200 II GEN and TBA-30R. Results may vary if a different instrument is used.

▪ Sensitivity

1.1 mg/dl – ACCENT-200
2.3 mg/dl – ACCENT-200 II GEN

▪ Linearity

up to 84 mg/dl – ACCENT-200
up to 89 mg/dl – ACCENT-200 II GEN

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

▪ Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 6000 mg/dl, intralipid up to 500 mg/dl, RF up to 500 IU/ml do not interfere with the test.

▪ Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT 200 n = 10	level 1	18.41	0.21	1.13
	level 2	45.37	0.40	0.88
ACCENT-200 II GEN n = 10	level 1	19,23	0,32	1,65
	level 2	54,77	0,22	0,40

▪ Method comparison

A comparison between Lp (a) values determined at ACCENT-200 (y) and at BS-400 (x) using 25 serum samples gave following results:

$$y = 0.9157x + 0.9956 \text{ mg/dl};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between Lp (a) values determined at ACCENT-200 II GEN (y) and at BS-400 (x) using 41 serum samples gave following results:

$$y = 1.0452x - 0.5364 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between CORMAY reagent (y) and commercially available assay (x) using 57 plasma samples gave following results:

$$y = 1.079x - 0.16 \text{ mg/dl};$$

$$R = 0.990 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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