## ACCENT-200 FERRITIN

Nr kat. 7-230

(PL)

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Cat. No **7-230** (EN)

### INTENDED USE

Diagnostic kit for determination of ferritin concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Ferritin is an iron-containing protein with a molecular weight of approximately 450 kD. It is found mainly in the human liver and spleen, where its function is to eliminate and store iron in the body, and is also found in small amounts in human serum. This amount varies according to the movement of iron in the body, and hepatitis and malignant tumors, may be seen to increase due to cell destruction or tumor cell production, independent of iron reserves. Consequently, the measurement of ferritin is considered to be useful in the diagnosis, treatment, assessment of disease progression, and postoperative prognosis for such disease conditions.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ferritin in a sample and anti-ferritin antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ferritin in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

**Package**  
 1-REAGENT 1 x 25 ml  
 2-REAGENT 1 x 11 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents are stable for 10 weeks on board the analyser at 2-10°C.

### Concentrations in the test

suspension of latex particles sensitized with anti-ferritin (rabbit) antibodies 0.07 w/v% (pH 7.3)  
 glycine buffer solution (pH 8.3)  
 preservative

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C.

- Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- EUH210 Safety data sheet available on request.

### SPECIMEN

Serum.

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.  
 0.9% NaCl is recommended as a reagent blank.

### Actions required:

When performing assays in analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S there is a probability of **cross-contamination** affecting the tests results: FERRUM II GEN - FERRITIN, LACTATE - FERRITIN, UREA - FERRITIN. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>6</sup>

serum	ng/ml
male	20 – 250
female	10 – 120

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL II (Cat. No 4-290) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY FERRITIN CALBRATORS kit (Cat. No 4-491) is recommended. 0.9% NaCl should be used as a calibrator 0. The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or HITACHI 917 and/or ACCENT MC240. Results may vary if a different instrument is used.

- LOQ:** 7.7 ng/ml – ACCENT-200
- Sensitivity:** 9.1 ng/ml – ACCENT MC240

- Linearity:**  
 up to 850 ng/ml – ACCENT-200  
 up to 1060 ng/ml – ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.98 g/dl, bilirubin up to 62 mg/dl, RF up to 520 IU/ml, intralipid up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [ng/ml]	SD [ng/ml]	CV [%]
ACCENT-200 n=20	level 1	99.78	1.50	1.50
	level 2	397.55	5.35	1.35
ACCENT MC240 n=20	level 1	99.57	1.98	1.99
	level 2	388.09	4.62	1.19
Reproducibility (day to day)		Mean [ng/ml]	SD [ng/ml]	CV [%]
ACCENT-200 n=80	level 1	101.3	4.03	4.0
	level 2	413.5	8.30	2.0
ACCENT MC240 n=80	level 1	97.8	3.49	3.6
	level 2	392.4	11.98	3.1

### Method comparison

A comparison between ferritin values determined at **ACCENT-200** (y) and at **ADVIA 1800** (x) using 60 serum samples gave following results:

$$y = 0.8883x + 2.7104 \text{ ng/ml};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between ferritin values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 49 samples gave following results:

$$y = 0.9296x - 7.1136 \text{ ng/ml};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements

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#### Calibration Rule

Rule	Logistic Log-5P
Sensitivity	1
Replicates	1
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

### • ACCENT-220S

#### Parameters

Test	FERRI	R1	180
No	52	R2	80
Full Name	Ferritin	Sample Volume	5
Standard No	52	R1 Blank	
Reac. Type	Fixed-time	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Concentration	9 570
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	2 20	Factor	
Incuba. Time	15	<input type="checkbox"/> Prozone check	
Unit	ng/ml	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

#### Calibration Rule

Rule	Logit-Log 5P
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

### • ACCENT-200 II GEN

#### Parameters

Test Name	FERRI	R1	180
Test No	52	R2	80
Full Name	Ferritin	Sample Volume	8
Reference No	52	R1 Blank	
Analy. Type	Fixed-time	Mixed Reag. Blank	
Pri. Wave.	578 nm	Concentration	9 830
Secon. Wave.		Linearity Limit	
Trend	Ascending	Substrate Limit	
Reac. Time	1 10	Factor	
Incuba. Time	15	<input type="checkbox"/> Prozone check	
Unit	ng/ml	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

#### Calibration Rule

Rule	Logistic Log-5P
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

### • ACCENT S120

Chem	FERRI	No.	052	Sample Type	SERUM				
Chemistry	FERRITIN	Print name	FERRITIN	Reaction Direction	positive				
Reaction Type	Fixed-time	Sec Wave		Decimal	0.1				
Pri Wave	578 nm	Incubation Time	16	Reaction Time	2 15				
Unit	ng/ml	Standard	5 5 5 5	Aspirated	20	Diluent	180	Reagent Vol	R1 180 R2 80
Blank Time		Decreased	5	Increased		Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>

Linearity range (Standard)	8.7	850	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000							
R1 Blank Abs	-40000	40000	On-board Stability	Day(s)							
Blank Response	-40000	40000	Reagent Alarm Limit								
Twin Chemistry		<input type="checkbox"/> Prozone Check	Enzyme Linear Extension	<input type="checkbox"/>							
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment		Pretreat Sample Vol		μL	Pretreat Sample Vol		μL

#### CALIBRATION SETTINGS

Math model	Logistic Log-5P		
Factor		Replicates	2

#### AUTO CALIBRATION

<input type="checkbox"/> Bottle Changed
<input type="checkbox"/> Lot Changed
<input type="checkbox"/> Cal Time

#### ACCEPTANCE LIMITS

Cal Time		Hour	
Slope Diff		SD	
Sensitivity		Repeatability	40000
Deter Coeff			

## ACCENT-200 FERRITIN

### • ACCENT MC240

Chem <input type="text" value="FERRI"/>	No. <input type="text" value="052"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="FERRITIN"/>	Print name <input type="text" value="FERRITIN"/>	
Reaction Type <input type="text" value="Fixed-time"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="570 nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="ng/ml"/>	Decimal <input type="text" value="0.1"/>	
Incubation Time <input type="text" value="21"/>	Reaction Time <input type="text" value="2"/> <input type="text" value="15"/>	
Blank Time <input type="text"/>		
Standard <input type="text" value="7"/> <input type="text" value="7"/> $\mu$ L	Aspirated <input type="text" value="20"/> $\mu$ L	Diluent <input type="text" value="180"/> $\mu$ L
Decreased <input type="text" value="7"/> $\mu$ L	<input type="text" value="20"/> $\mu$ L	<input type="text" value="180"/> $\mu$ L
Increased <input type="text"/> $\mu$ L	<input type="text"/> $\mu$ L	<input type="text"/> $\mu$ L
<input type="text"/> Sample Blank	<input checked="" type="checkbox"/> V	Auto Rerun

Linearity range (Standard) <input type="text" value="9,1"/> <input type="text" value="1000"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="text"/> Sample Pretreatment	<input type="text"/> Control Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> $\mu$ L	<input type="text"/> Calibrator Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> $\mu$ L	

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model <input type="text" value="Logistic Log-5P"/>	Factor <input type="text"/>	Replicates <input type="text" value="2"/>	
			<input type="checkbox"/> Bottle Changed
			<input type="checkbox"/> Lot Changed
			<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time <input type="text"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>	

### • ACCENT M320

Chem <input type="text" value="FERRI"/>	No. <input type="text" value="052"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="FERRITIN"/>	Print name <input type="text" value="FERRITIN"/>	
Reaction Type <input type="text" value="Fixed-time"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="570 nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="ng/ml"/>	Decimal <input type="text" value="0.1"/>	
Incubation Time <input type="text" value="16"/>	Reaction Time <input type="text" value="3"/> <input type="text" value="22"/>	
Blank Time <input type="text"/>		
Standard <input type="text" value="5"/> $\mu$ L	Aspirated <input type="text" value="20"/> $\mu$ L	Diluent <input type="text" value="180"/> $\mu$ L
Decreased <input type="text" value="5"/> $\mu$ L	<input type="text" value="20"/> $\mu$ L	<input type="text" value="180"/> $\mu$ L
Increased <input type="text"/> $\mu$ L	<input type="text"/> $\mu$ L	<input type="text"/> $\mu$ L
<input type="text"/> Sample Blank	<input checked="" type="checkbox"/> V	Auto Rerun

Linearity range (Standard) <input type="text" value="9,3"/> <input type="text" value="1040"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="text"/> Sample Pretreatment	<input type="text"/> Control Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> $\mu$ L	<input type="text"/> Calibrator Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> $\mu$ L	

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model <input type="text" value="Logistic Log-5P"/>	Factor <input type="text"/>	Replicates <input type="text" value="2"/>	
			<input type="checkbox"/> Bottle Changed
			<input type="checkbox"/> Lot Changed
			<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time <input type="text"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>	

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