ACCENT-200 RF

Nr kat. 7-237

(PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of rheumatoid factor levels used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT 400, ACCENT Neo200 and ACCENT M320.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Rheumatoid factor (RF) is an autoantibody against human IgG commonly seen in sera at a high concentration in some conditions, particularly in patients with rheumatoid arthritis. The measurement of RF value is useful in evaluating the diagnosis, effects of therapy and prognosis of RA, systemic lupus erythematosus, chronic hepatopathy, etc.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between RF in a sample and denatured human IgG which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of RF in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-Reagent 1 x 37 ml
 2-Reagent 1 x 10.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks (ACCENT 200), for 9 weeks (ACCENT S120) or for 6 weeks (ACCENT M320).

Concentrations in the test

suspension of latex particles sensitized with denatured human IgG (pH 7.3) 0.17 w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

0.9% NaCl is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S there is a probability of **cross-contamination** affecting the tests results: RF and CALCIUM, CALCIUM ARSENAZO. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ⁴

serum, plasma	< 18 IU/ml
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, the CORMAY RF CALBRATORS kit (Cat. No 4-277) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 4 weeks (ACCENT 200), every 3 weeks (ACCENT S120), every 2 weeks (ACCENT M320), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-220S and ACCENT MC240. Results may vary if a different instrument is used.

Sensitivity

11 IU/ml – ACCENT-220S
 11.8 IU/ml – ACCENT MC240

Linearity

up to 95 IU/ml – ACCENT-220S
 up to 115 IU/ml – ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 1.0 g/dl, bilirubin up to 66 mg/dl, intralipid up to 20 g/l do not interfere with the test.

Precision

Repeatability (run to run)		Mean [IU/ml]	SD [IU/ml]	CV [%]
ACCENT-220S n=10	level 1	15.36	0.9	5.89
	level 2	32.73	0.31	0.93
ACCENT MC240 n=20	level 1	14.72	0.73	4.94
	level 2	35.04	0.37	1.05
Reproducibility (day to day)		Mean [IU/ml]	SD [IU/ml]	CV [%]
ACCENT-220S n=20	level 1	17.16	0.33	1.94
	level 2	33.47	0.47	1.40
ACCENT MC240 n=80	level 1	17.7	1.25	7.0
	level 2	31.6	0.89	2.8

Method comparison

A comparison between RF values determined at ACCENT-220S (y) and at BS-400 (x) using 37 samples gave following results:

$$y = 0.9457x + 2.8273 \text{ U/l};$$

$$R = 0.982 \quad (R - \text{correlation coefficient})$$

A comparison between RF values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU 680 (x) using 60 serum samples gave following results:

$$y = 1.0006x - 0.6128 \text{ U/l};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Galvin J.P. et al.: Particle enhanced photometric immunoassay systems., Clin Lab. Assay (Pap.Annu.Clin.Lab.Assays Conf.), 4th, 73 (1983).
- Singer J.M. et al.: The latex fixation test. I. Application to the serologic diagnosis of rheumatoid arthritis, Amer. J. Med., 888 (1956)
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 960, (2006).
- Internal reference range studies.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	3
Interval (day)	28
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	RF	R1	210
No	31	R2	70
Full Name	RF	Sample Volume	8
Standard No	31	R1 Blank	
Reac. Type	Fixed-time	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	11 95
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	4 12	Factor	
Incuba. Time	20	<input type="checkbox"/> Prozone check	
Unit	IU/ml	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Logit-Log 5P
Sensitivity	1
Replicates	2
Interval (day)	28
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	RF	No.	031	Sample Type	SERUM
Chemistry	RF	Print name	RF	Reaction Direction	positive
Reaction Type	Fixed-time	Sec Wave		Decimal	0.1
Pri Wave	578 nm	Incubation Time	20	Reaction Time	3 8
Unit	IU/mL	Blank Time		Standard	6 μL
Sample Vol	6 μL	Aspirated	30 μL	Diluent	120 μL
Decreased	6 μL	Reagent Vol	180 μL	R1	180 μL
Increased	6 μL	R2	60 μL	Sample Blank	V
Auto Rerun					

Linearity range (Standard)	12.3	120	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000							
R1 Blank Abs	-40000	40000	On-board Stability	63 Day(s)							
Blank Response	-40000	40000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension								
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		Calibrator Pretreatment							
Pretreat Sample Vol		μL	Pretreat Sample Vol		μL						

CALIBRATION SETTINGS

Math model	Spline		
Factor		Replicates	2

AUTO CALIBRATION

<input type="checkbox"/> Bottle Changed	
<input type="checkbox"/> Lot Changed	
<input type="checkbox"/> Cal Time	

ACCEPTANCE LIMITS

Cal Time	504 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	40000
Deter Coeff			

ACCENT-200 RF

• ACCENT MC240

Chem <input type="text" value="RF"/>	No. <input type="text" value="031"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="RF"/>	Print name <input type="text" value="RF"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="605 nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="IU/mL"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="21"/>	
	Reaction Time <input type="text" value="19"/> <input type="text" value="28"/>	
Standard <input type="text" value="6"/> <input type="text" value="6"/> <input type="text" value="30"/> <input type="text" value="120"/>	Aspirated <input type="text"/>	Diluent <input type="text"/>
Decreased <input type="text"/>	Reagent Vol <input type="text" value="180"/>	R1 <input type="text" value="180"/> <input type="text" value="60"/>
Increased <input type="text"/>		R2 <input type="text"/>
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="11.8"/> <input type="text" value="115"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text" value="42"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Calibrator Pretreatment	
Pretreat Sample Vol <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> μ L

CALIBRATION SETTINGS	
Math model <input type="text" value="Spline"/>	
Factor <input type="text"/>	Replicates <input type="text" value="2"/>

ACCEPTANCE LIMITS	
Cal Time <input type="text" value="336"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>	

• ACCENT M320

Chem <input type="text" value="RF"/>	No. <input type="text" value="031"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="RF"/>	Print name <input type="text" value="RF"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="570 nm"/>	Sec Wave <input type="text" value="800 nm"/>	
Unit <input type="text" value="IU/mL"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="20"/>	
	Reaction Time <input type="text" value="28"/> <input type="text" value="30"/>	
Standard <input type="text" value="5"/> <input type="text" value="5"/> <input type="text" value="30"/> <input type="text" value="120"/>	Aspirated <input type="text"/>	Diluent <input type="text"/>
Decreased <input type="text"/>	Reagent Vol <input type="text" value="180"/>	R1 <input type="text" value="180"/> <input type="text" value="60"/>
Increased <input type="text"/>		R2 <input type="text"/>
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="6.1"/> <input type="text" value="117"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text" value="42"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Calibrator Pretreatment	
Pretreat Sample Vol <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> μ L

CALIBRATION SETTINGS	
Math model <input type="text" value="Spline"/>	
Factor <input type="text"/>	Replicates <input type="text" value="2"/>

AUTO CALIBRATION	
<input type="checkbox"/> Bottle Changed	
<input type="checkbox"/> Lot Changed	
<input type="checkbox"/> Cal Time	

ACCEPTANCE LIMITS	
Cal Time <input type="text" value="336"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>	

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